ADJUSTABLE CORTICAL FIXATION SYSTEM

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NOT251_170929

Version du 29/09/2017

This release is the last update of the instructions of use and replace the previous edition

ADJUSTABLE CORTICAL FIXATION SYSTEM
ANCHOR SYSTEM
SINGLE USE STERILE PRODUCT

DESCRIPTION

Fixation system for tendons and ligaments.

The adjustable fixation system comprises a loop for graft suspension and a titanium button for the cortex.

MATERIALS

- Materials implanted:

Cortical button: Non-ferromagnetic titanium alloy (Ti TA6V ELI)

Adjustable loop: 100% polyethylene

Extended plate for cortical fixation system (optional): Non-ferromagnetic titanium alloy (Ti TA6V ELI)

Neither animal nor human origin. Non-absorbable

- Materials in contact with the patient during the operation:

Traction wires: Siliconized polyethylene terephthalate + biocompatible colorant D&C green#6

INDICATIONS

Ligament or tendon anchor system for orthopaedic surgery, in particular for reconstruction of anterior cruciate ligaments

PERFORMANCE

This system facilitates ligament or tendon repair/reconstruction by graft fixation in the bone.

It has been approved for long-term implantation in the human body.

CONTRAINDICATIONS

- 1- Insufficient quantity or quality of bone, deficient bone stock or bone surface that could compromise appropriate fixation of the implant
- 2- Do not implant on an infected site or in case of allergies to the implant materials
- 3- Not to be used for the bone – patella tendon – bone technique.
- 4- Do not implant in children while they are growing.

POSSIBLE ADVERSE EFFECTS

Possible complications linked to the implant:

- 1- Breaking of the loop, breaking of the traction wires, breaking of the button
- 2- Failure when placed under tension
- 3- Migration or upturn of the button
- 4- Reaction to a foreign body

Possible complications linked to insertion of the implant:

- 1- Pain in the incision area or surgical site
- 2- Local and transient inflammatory reactions
- 3- Infection
- 4- Wound dehiscence

WARNING

1- Any decision to remove an implant must take into account the fact that a second surgical operation may involve a risk for the patient. Removal must be followed by suitable postoperative care.

2- The patient must be informed of the limitations inherent in the nature and function of the implant, in particular regarding the weight and constraints placed by the body on the system before the bone has sufficiently healed.

3- Ensure that the bone tunnel has the correct length and diameter.

PRECAUTIONS BEFORE USE

Before use, check the integrity of the implant and the packaging.

Do not use if the implant, labels and/or protective packaging is damaged.

Do not use the product after the expiry date indicated on the label. This could compromise the performance of the system and the patient's safety.

Unpacking and handling of the implant must be carried out in accordance with asepsis standards, to ensure sterility.

When opening the cardboard box, the surgeon must check that the batch number and model match to those indicated on the label on the inside cover. If they do not, use of the implant is strictly forbidden.

The harmlessness and compatibility of the implant have not been assessed in case of postoperative testing by magnetic resonance imaging (MRI). However, according to medical literature, the materials used in the implant (polyethylene and titanium) are recognised as non-ferromagnetic, and in view of the quantities used the system may be considered as compatible with MRI tests.

This implant can be used with standard instruments available on the market and specified for orthopaedic surgery for ligament or tendon anchorage.

Surgical instruments are subject to wear during normal use. Instruments that have been used many times or been subject to very high stress are likely to break. Instruments must only be used for their intended purpose. COUSIN BIOTECH recommends that instruments should be regularly inspected for wear and bending.

DO NOT REUSE

DO NOT RESTERILIZE

The implant is for single use, as specified on the label. It must not on any account be reused or re-sterilized (potential risks include, but are not limited to: loss of sterility of the product, risk of infection, loss of effectiveness of the product).

COUSIN BIOTECH does not provide any guarantee or recommendations regarding the use of a specific make of fixation device. The properties of these devices are subject to alteration by the manufacturer, over which COUSIN BIOTECH has no control.

This device must **only** be implanted by a qualified surgeon trained in the use of the product. Avoid using staples or sharp-edged screws.

SURGICAL TECHNIQUE

Preliminary information

The device is optimized for use of the 4-strand semitendinosus technique with retrodrilling and the all-inside and inside-out techniques.

In case of retrodrilling, COUSIN BIOTECH recommends a diameter of 4.5mm for the tunnel to insert the implant. In case of a full tunnel, it is preferable to use the additional extended plate.

Use of the device with a tunnel diameter of less than 4.5mm has not been validated. If the diameter is greater than 4.5mm, the surgeon can decide whether or not to use the extended plate.

Technique

1. Suture the graft and put it in the loop.

2. Create a bone tunnel to insert the implant with a diameter of 4.5mm (recommended diameter) and a graft cavity suited to its diameter.

3. Adjustable system: go to step 4

Fixed system: Choosing the size of the implant. Total length of the femoral tunnel (between the two cortices) - length of the graft tunnel +5cm. If the result is not equal to the size of the implant, round up to the larger implant size.

4. Adjustable system: Insert the suture threads and adjustment threads through the bone tunnel.

Fixed system: Pass the threads through the bone tunnel.

5. Pull on one of the sutures until the button comes out of the bone tunnel through the cortex.

6. Tip over the button against cortex surface by pulling on the graft.

7. Fixed system: go to step 10 and then step 12.

Adjustable system: While applying counter tension to the graft, pull the adjustment thread to draw the graft into the tunnel.

8. While still applying counter tension to the graft, slide the adjustment threads one after another to adjust the loop to the required length.

9. Use the coloured adjustment thread to lock the tension applied.

10. Fix the graft in the opposite bone tunnel using the required method.

11. If necessary, slide the adjustment threads alternately once again to obtain the required tension.

12. Remove the suture thread(s) from the button.

13. Select the adjustment threads.

In case of use of the extended plate

It may be necessary to use an extended plate to increase the contact surface of the button on the cortex (full tibial or femoral tunnel).

Repeat steps 1 to 7 above.

8. Remove one of the two sutures

9. Put the button on the extended plate. Ensure that the remaining suture is on the split side to avoid hindering the interlocking.

10. Fix the graft in the opposite bone tunnel using the required method.

11. If necessary, slide the adjustment threads alternately once again to obtain the required tension.

12. Remove the remaining suture thread from the button.

13. Select the adjustment threads.

STORAGE PRECAUTIONS

Store in a dry place away from sunlight, at room temperature and in the original packaging.

EXPLANTATION AND DISPOSAL OF THE DEVICE

COUSIN BIOTECH recommends that explantation and handling should be carried out in accordance with the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

All explanted devices must be returned for analysis according to the protocol in force. This is available on request from COUSIN BIOTECH. It should be noted that any implants that cannot be cleaned and disinfected before shipping must be placed in a sealed package.

Disposal of explanted medical devices must be carried out in accordance with the standards in force in the country, concerning the disposal of infectious waste.

There are no special recommendations for the disposal of non-implanted devices.

INFORMATION REQUESTS AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH undertakes to make every effort to produce and supply a high quality medical device. However if a health professional (customer, user, consultant etc.) has a complaint or grounds for dissatisfaction with a product, regarding its quality, safety or performance, he/she must inform COUSIN BIOTECH as soon as possible.

In the event of malfunctioning of an implant or if it contributes towards a serious adverse effect for the patient, the healthcare centre must follow the legal procedures in force in the country, and immediately inform COUSIN BIOTECH.

In all correspondence, please state the reference, batch number, name and address of a consultant, and provide a full description of the incident or complaint.

Brochures, documentation and the surgical technique are available on request from COUSIN BIOTECH and its distributors.

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www.cousin-biotech.com/IFU

IFU is available on our webpage for the corresponding product. You can view, save and print the IFU. Please choose the pages of the languages via the printing menu.

If requested, we will send you a free printed copy of the Instructions of Use in the language you require within 7 days.

[You can use the order form provided on our webpage for this purpose or by mail at ifurequest@cousin-biotech.com](#)

	en Batch number fr Numéro de lot de Chargennummer it Numero di lotto es Número de lote pt Número do lote ελ Αριθμός παρτίδας nl Batchnummer tu Parti numarası		en Caution (See instructions for use) fr Attention, voir notice d'instructions de Siehe Gebrauchsanweisung it Vedere manuale istruzioni es Veanse las instrucciones de uso pt Atenção, consultar nota de instruções ελ Δείτε τις οδηγίες χρήσης nl Zie gebruiksinstructies tu İçeriği inceleyin		en Manufacturer fr Fabricant de Hersteller it Produttore es Fabricante pt Fabricante ελ Κατασκευαστής nl Producent tu Üretici
	en Do not re-use fr Ne pas réutiliser de Nicht wiederverwenden it Non riutilizzare es No volver a utilizar pt Não reutilizar ελ Να μην επαναχρησιμοποιείται nl Niet hergebruiken tu Tekrar kullanmayın		en Reference on the brochure fr Référence du catalogue de Verweis auf der Broschüre it Referimento del catalogo es Referencia al folleto pt Referência do catálogo ελ Αριθμός καταλόγου nl Referentie in de brochure tu Broşür referansı		en Do not re-sterilize fr Ne pas restériliser de Nicht erneut sterilisieren it Non risterilizzare es No volver a esterilizar pt Não reesterilizar ελ Να μην επαναποστειρώνεται nl Niet opnieuw steriliseren tu Yeniden sterilize etmeyin

	en Keep in a dry place fr Conserver au sec de Trocken lagern it Conservare in un luogo asciutto es Manténgase en un lugar seco pt Manter em lugar seco ελ Φυλάξτε το προϊόν μακριά από την υγρασία nl Op een droge plaats bewaren tu Kuru bir yerde muhafaza edin		en Use before: year and month fr Utiliser jusqu'à : année et mois de Verbrauchen bis : Jahr und Monat it Utilizzare entro e non oltre: anno e mese es Úsese antes de: año y mes pt Utilizar até : ano e mês ελ Ημερομηνία λήξης: έτος και μήνας nl Te gebruiken vóór: jaar en maand tu Son kullanma tarihi: yıl ve ay
	en Keep away from sunlight fr A stocker à l'abri de la lumière du soleil de Vor Sonnenlicht geschützt lagern it Tenere lontano dalla luce solare es Mantener alejado de la luz solar pt Manter ao abrigo da luz solar ελ Προφυλάξτε το προϊόν από το φως nl Vrij van licht houden tu Işıktan uzak tutun		en Do not use if the packaging is damaged fr Ne pas utiliser si l'emballage est endommagé de Nicht verwenden, wenn Packung beschädigt ist it Non utilizzare se la confezione è danneggiata es No usar si el embalaje está dañado pt Não usar se a embalagem estiver danificada ελ Να μην χρησιμοποιείται, εάν η συσκευασία έχει υποστεί κάποια βλάβη nl Niet gebruiken indien de verpakking beschadigd is tu Paketi hasar görmüş ürünleri kullanmayın

	en Sterilized by ethylene oxide fr Produit stérile. Méthode de stérilisation : oxyde d'éthylène de Mit Ethylenoxid sterilisiert it Sterilizzato all'ossido di etilene es Esterilizado con óxido de etileno pt Produto estéril. Método de esterilização: óxido de etileno ελ Στείρο προϊόν. Μέθοδος αποστείρωσης: αιθυλενοξείδιο nl Gesteriliseerd met ethyleenoxide tu Etilen oksit ile sterilidir
---	---

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