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Made in France

Caution : Federal law (USA) restricts this device to sale,
distribution and use by or on the order of a physician.

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Sterile spinal dynamic posterior stabilization device

DESCRIPTION

The sterile spinal dynamic posterior stabilization device BDyn is intended to restore the stabilization of the non-cervical segment of the spine by preserving the anatomical lordosis and the deadening of the intervertebral joint. It is composed of the BDyn spinal shock absorber fixed on the vertebrae thanks to the dedicated polyaxial pedicle screws BDyn with plug and can be replaced according to the decision of the surgeon by a dedicated replacement rigid rod BDyn.

MATERIALS

BDyn implant: Titanium alloy Ti6Al4V ELI (ISO 5832-3; ASTM F 136), Polydimethyl siloxane (PDMS) MED 4770 unrestricted Polycarbonate urethane (PCU) Bionate® II 80A

BDyn ancillaries: Stainless steel, Titanium, Polypropylene
No human or animal origin – Non resorbable.

INDICATIONS

The sterile spinal dynamic posterior stabilization device BDyn is intended for posterior stabilization from thoracic vertebrae T10 to sacrum S1 with or without bone graft for the following indications:

- Degenerative intervertebral disc disease and/or articular facets confirmed by further examinations
- Spinal canal stenosis
- Degenerative spondylolisthesis grade 1
- Segmental instability

PERFORMANCE

The BDyn spinal shock absorber of the sterile spinal dynamic posterior stabilization device BDyn consists of a metallic cylindrical hollow part containing elastomer components made of silicone elastomers (polydimethyl siloxane (PDMS)) and long-term implantable polycarbonate urethane (PCU) which are bending out under the effect of a metallic piston rod connected with the vertebra of the treated segment by the polyaxial pedicle screws BDyn. The elastomer components assure the absorption of the mechanical loads applied on the intervertebral joint in compression, pulling, flexion-traction and lateral flexion. Some configurations enable fixation across several spinal segments.

CONTRAINDICATIONS

- Active infectious process or significant risk of infection (immunocompromise)
- Signs of local inflammation
- Fever or leukocytosis
- Morbid obesity
- Pregnancy
- Mental illness
- Grossly distorted anatomy caused by genital abnormalities
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of congenital abnormalities, elevation of the sedimentation rate unexplained by other diseases, elevation of the white blood count
- Suspected or documented metal allergy or intolerance
- Any case where the implant components selected for use would be too large or too small to achieve a successful result
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance
- Any patient unwilling to follow postoperative instructions
- Any case not describe in the indications
- Traumas (i.e. fracture or dislocation)
- Abnormal curvatures (i.e. scoliosis and/or hyper lordosis)
- Tumors.
- Spondylolisthesis grade 2 and more
- Pseudarthrosis and/or failed previous fusion
- Severe bone resorption, osteomalacia, severe osteoporosis

UNDESIRABLE SIDE EFFECTS

All of the possible adverse events associated with spinal surgery and without instrumentation are possible:

- Infection
- Pseudomeningocele, fistula, breach dura, persistent CSF leakage, meningitis
- Loss of neurological function, sensorial and/or motor, including complete or incomplete paralysis, dysesthesias, hyperesthesia, anesthesia, paresthesia, appearance of radiculopathy, and/or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, and/or visual deficits
- *Cauda equina* syndrome, neuropathy, transient or permanent neurological deficits, paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, and/or muscle loss
- Urinary retention or loss of bladder control or other types of urological system compromise
- Scar formation possibly causing by a neurological compromise or compression around nerves and/or pain
- Fracture, microfracture, resorption, damage or penetration of any spinal bone (including the sacrum, pedicles, and/or vertebral body) and/or bone graft or bone graft harvest site at, above and/or below the level of surgery
- Herniated *nucleus pulposus*, disc disruption or degeneration at, above, or below the level of surgery, canal adjacent stenosis
- Non-union or pseudarthrosis, delayed union. Mal-union
- Cessation of any potential growth of the operated portion of the spine
- Loss of or increase in spinal mobility or function
- Inability to perform the activities of daily living
- Bone loss or decrease in bone density
- Graft donor site complications including pain, fracture, or wound healing problems
- Ileus, gastritis, bowel obstruction or loss of bowel control or other types of gastrointestinal system compromise
- Hemorrhage, hematoma, occlusion, seroma, edema, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to blood vessels, or other types of cardiovascular system compromise
- Reproduction system compromise, sterility, sexual dysfunction
- Development of respiratory problems, e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc
- Change in mental status
- Death

All of the possible adverse events associated with spinal surgery and instrumentation are possible. A listing of potential adverse events linked to the medical device includes:

- Early or late loosening of any or all of the components
- Disassembly, bending and/or breakage of any or all of the components (screw breakage)
- Foreign body (allergic) reaction to implants, debris, corrosion products (from crevice, fretting, and/or general corrosion), including metallosis, tumor formation and/or autoimmune disease
- Pressure on the skin from components parts with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, fibrosis, neurosis, and/or pain
- Tissue or nerve damage caused by improper positioning and placement of implants or instruments
- Post-operative change in spinal curvature, loss of correction, height, and/or reduction

PRECAUTIONS FOR USE

The BDyn devices are delivered sterile, sterilized by ethylene oxide for the BDyn spinal shock absorber and by Gamma irradiation for the pedicle screws with plug and the replacement rigid rods.

Before any use, inspect the integrity of the packaging and device (including peelable pouches).

Do not use in the event of deterioration of the labels and/or the device and/or the packaging.

Do not use if the device is out of date.

The sterile spinal dynamic posterior stabilization device BDyn is only intended to adult human female or male more than 40 Kg.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

The implantation of the BDyn spinal shock absorber has to be carried out only with the dedicated polyaxial pedicle screws BDyn manufactured by the company COUSIN BIOTECH.

The eventual substitution of the BDyn spinal shock absorber has to be carried out only with the dedicated replacement rigid rods BDyn manufactured by the company COUSIN BIOTECH.

You should never use implants made of stainless steel and titanium alloy in the same construct.

A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where many extenuating circumstances may compromise the results. This device system is not intended to be the sole means of spinal support. No spinal implant can withstand body loads without the support of bone. In this event, bending, loosening, disassembly and/or breakage of the device(s) will eventually occur.

Preoperative and operating procedures, including precise knowledge of suitable surgical techniques, and proper selection of the good reference of the device adapted to the patient and its narrow setting up are important considerations in the successful use of the device by the surgeon.

The BDyn device must be implanted only by a qualified surgeon, having knowledge in the use of the product and who has the knowledge of the anatomy, spinal surgery, pedicle screws fixation technique and specific BDyn device surgical technique

Postoperative precautions

The BDyn device has not been tested for heating or migration in the MR environment.

The BDyn device has been evaluated for safety and compatibility in the MR (Magnetic Resonance) environment. According to medical literature and in view of the materials used in the implant (titanium alloy TiA614V, PDMS and PCU), the device can be considered as MR safe.

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

In accordance with the labelling of this product, BDyn implants are for single use; they must not in any case be reused or re-sterilised (potential risks include, but are not limited to, loss of sterility of the products, risk of infection, loss of product efficacy, recurrence, etc.)

SURGICAL TECHNIQUE

The placement of the spinal dynamic posterior stabilization device BDyn begins by the preparation of the pedicles of the vertebrae and by the pedicle aim and the fixation of the pedicle screws.

When tapping the vertebrae to prepare the hole before the insertion of the pedicle screw, flexion efforts on the tap are to be avoided. It must stay in the lined-up in order to make as straightest a pre-hole as possible to ensure as secure a screw placement as possible, with a straight guidance way. Thus, when tapping the vertebrae in order to stay in the line-up and to create a straight pre-hole.

Viewed in the sagittal plane, the pedicle screws must be as parallel as possible in order to let the biggest available space to place the cylindrical body of BDyn spinal shock absorber and facilitate its insertion between the head of the polyaxial pedicle screws. It is preferable to be sure that the adequate implant, with an appropriate size is used with the instrumentation which suit. A 2 mm space between the top of the cylindrical body of the BDyn spinal shock absorber and the lower part of the screw head must be respected to insure the shock absorption necessary for the good working of the BDyn device. For that, the surgeon can use the trial device M50 to check if there is enough space to place the BDyn spinal shock absorber.

The mobile piston rod of the BDyn spinal shock absorber must be fixed in the head of the screw concentrically to the cylindrical body.

The surgery ends by the final tightening by using the dedicated torque limiting of the plugs on the head of the polyaxial pedicle screws by using systematically the suitable antirotator.

It is important to manipulate correctly the implants. The components of the implants to be arched do not have to be against-arched in the opposite direction. In any case they have to be chipped. These operations could effectuate the concentration of internal stresses which could become the place of an eventual failure of the implant.

PRESENTATION AND WARNING ABOUT ANCILLARIES

- The BDyn ancillaries are Class I medical devices, intended for temporary use and are re-usable.

Unlike the BDyn implants, the ancillaries included in the instrument set are furnished non sterile.

Before use, the ancillaries must be:

- cleaned by the appropriate mode of cleaning
- sterilized in an autoclave by moist heat

After use, the ancillaries have to be cleaned and decontaminated according to a suitable protocol especially by taking into account the reduction of the risk of transmission of unconventional transmissible agents - UTA.

In case of a patient suspected of having Creutzfeld-Jacob Disease (CJD), perform suitable decontamination.

In the event of a confirmed case of CJD, incinerate the ancillaries.

CLEANING OF BDYN ANCILLARIES

Decontamination and cleaning are performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

1. Manual Pre – cleaning process

- The detachable instruments have to be dismantled.
- The instruments have to be brushed with a soft nylon brush until all visible residuals are removed.
- Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position.
- Then, the instruments have to be treated in an ultrasonic bath with 0.5% solution of Neodisher Mediclean at 40°C for 10 min.
- After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position.

2. Cleaning in a washer / disinfecter

For automated cleaning, the instruments have to be placed in the instrument rack of the washer / disinfecter (type Miele G 7735 CD or equivalent). The following steps have to be performed:

- 1 min pre-cleaning with cold tap water and drain,
- 3 min pre-cleaning with cold tap water and drain,
- 5 min cleaning with deionized water at 45°C with a 0.5 % solution of Neodisher Mediclean and drain,
- 3 min neutralization/rinsing with deionized water and drain,
- 2 min final rinse with deionized water and drain.

After the cleaning operation, the instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization.

In case of non respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

STERILIZATION OF BDYN ANCILLARIES

Sterilization is performed under the responsibility of the healthcare centers.

The following prescribed method and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

Check if the instruments are placed in the right position inside the trays as described on the position map fixed on the tray before the sterilization.

COUSIN BIOTECH recommends sterilization by moist heat (steam). Ancillaries can be sterilized in a:

- moist heat autoclave used at 134°C during 18 minutes,
- pre-vacuum moist heat autoclave used after 3 pre-vacuum phases, at 134°C during 5 minutes followed by a drying time of 10 min.

Always verify the integrity and the correct operation of the ancillaries after each cleaning/sterilization cycle and before each use.

In case of non respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent sterilization process for which it will be fully responsible.

STORAGE PRECAUTIONS

- The BDyn devices must be stored in a clean, dry and tempered place under atmospheric pressure, away from sunlight and rays in its original packaging. Ensure that the box is not dropped, hit or crushed.
- The ancillaries must be stored in their dedicated container or equivalent packaging to prevent any deterioration.

EXPLANATION AND ELIMINATION OF DEVICES

Explanation and handling should be done following recommendations of ISO 12891-1:2011 « Implants for surgery – Retrieval and analysis of surgical implants » Part 1: « Retrieval and Handling ».

Any explanted device must be send back, for analysis, following the current protocol. This protocol is available on demand to COUSIN BIOTECH. It is important to note that any implant that should not have been cleaned and disinfected before expedition must be contained in a sealed package.

The elimination of explanted medical device must be conducted in accordance with standards in the country for the disposal of infectious hazards waste.

The elimination of a non-implanted device is not the subject of specific recommendations.

The surgeon will decide the necessity of the resection of the BDyn, its replacement or the rachis fusion according to the known and recognized techniques.

INFORMATION REQUEST AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the healthcare center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time. For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim. Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

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e-IFU


IFU is available on our webpage for the corresponding product. You can view, save and print the IFU. Please choose the pages of the languages via the printing menu.

www.cousin-biotech.com/IFU

If requested, we will send you a free printed copy of the Instructions of Use in the language you require within 7 days.

You can use the order form provided on our webpage for this purpose or by mail at ifurequest@cousin-biotech.com

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