INTRASPINE®

DEVICE FOR THE INTER SPINOUS SPACE WITH LAMINAR SUPPORT

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	4
de	<u>Gebrauchsanweisung</u>	Seite	6
it	<u>Istruzioni per l'uso</u>	Pagina	8
es	<u>Instrucciones de uso</u>	Página	10
pt	<u>Nota de instruções</u>	Página	12
ελ	<u>Οδηγίες χρήσης</u>	σελίδα	14
pl	<u>Instrukcja obsługi</u>	Strona	16
cs	<u>Pokyny na použitie</u>	Strana	18
ru	<u>Инструкция по применению</u>	Страница	20

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This release is the last update of the instructions of use and replace the previous edition

IntraSPINE®
DEVICE FOR THE INTER-SPINOUS SPACE WITH LAMINAR SUPPORT
SINGLE USE STERILE PRODUCT

DESCRIPTION

The IntraSPINE® prosthesis is an interlaminar medical device that comprises a dimethyl siloxane (silicon) wedge (available in several sizes) adapted to the intervertebral space and a ligament in polyethylene terephthalate.

The wedge is covered with woven polyethylene terephthalate with an area covered by sheeting on the front face in the medullary zone. Tension is applied to the ligament by a titanium loop (TA6V).

MATERIALS

IntraSPINE® prosthesis: Dimethyl siloxane, Polyethylene terephthalate, Titanium.

IntraSPINE® ancillaries and trial prostheses: Stainless steel, Acetal copolymer, Dimethyl siloxane.

Non-human and non-animal origin – Non-resorbable.

INDICATIONS

The IntraSPINE® prosthesis is recommended in lumbar surgery in the following cases:

- Arthropathic facet-syndrome
- Foraminal stenosis
- Degenerative discopathy
- Intervertebral ligament insufficiency

PERFORMANCE

The IntraSPINE® prosthesis is an elastic support device with a cushioning effect. Its interlaminar implantation enables it to be close to the centre of rotation, thereby producing the cushioning effect and providing tension on the interspinous ligamentary structures due to the elasticity of the supra-spinous wedge.

The principal function of the ligaments surrounding the upper and lower vertebrae in relation to the area receiving the prosthesis is to assist the posterior ligament during flexion.

IntraSPINE® has been validated for a long term implantation in the human body.

CONTRAINDICATIONS

1) Do not use the prosthesis in the following cases:

- Allergy to one of the components
- Growing child
- Infected site
- Pregnant woman

2) Factors likely to compromise successful implantation:

- Severe osteoporosis
- Major deformations of the spinal column
- Local bone tumors
- Systemic or metabolic disorders
- Infectious diseases
- Obesity
- Drug addictions
- Intense physical activity, for example competitive sports or hard labor

UNDESIRABLE SECONDARY EFFECTS

1) Complications likely related to the device:

- Infection of the prosthesis, rupture of the ligament, loosening, fixed wing breakage, permanent elongation of the ligament, displacement of the wedge, spondylolisthesis to the dynamically stabilized segment, persistent pain attributed to the dynamically stabilized level and required revision surgery with interbody fusion at that level;

Spinous process fractures: intraoperative, post-operative or during device removal (removal due to infection) and can be associated with dislocation of the implant;

Late reoperation due to absence of recovery;

In certain cases the medical instrument had to be removed and pedicle screw fusion had to be performed.

2) Complications unlikely related to the device but related to surgery or to the patient's conditions:

Dural tear, cerebral spinal fluid leak, hematoma, cardiovascular problems, deep venous thrombosis and embolism, acute myocardial infection, superficial wound infection, short term unexplained visual disturbance, self-limiting pseudo-radicular pain in the other leg, late-onset lumbar spinal canal hematoma formation.

POSTOPERATIVE PRECAUTIONS

The patient must be informed of the postoperative cares and of precautions.

The use of a brace could be necessary but is not mandatory. This decision is the surgeon one.

PRECAUTIONS FOR USE

Before usage, verify the integrity of the prosthesis, trial prostheses, ancillaries and packaging.

Do not use in the event of deterioration of the prosthesis and/or the labels, and/or the trial prostheses and/or the ancillaries and/or the packaging.

Do not use if the prosthesis has exceeded its use-by date.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the prosthesis are equal to the data on the label of the inner pouch. Otherwise the use of the prosthesis is strictly prohibited.

Prevent the prosthesis, the trial prostheses and ancillaries from coming in to contact with objects that may alter their surface.

Operation conditions: Temperature from +32°C up to 42°C. Implants should be resistant to biologic fluids and the discharge of body tissues, which they contact during use. Therefore, the IntraSPINE® prosthesis has been validated for a long term implantation and, currently, Cousin Biotech has 7 years of feedback concerning these products.

The IntraSPINE® prosthesis, trial prostheses and ancillaries must only be used by a qualified surgeon who is trained in spinal surgery and to the surgical technique of IntraSPINE®. The information contained in this notice is necessary but not sufficient to control the surgical technique.

The prescription of the device is decided on by the surgeon, the only person qualified to do so.

Each IntraSPINE® prosthesis, trial prostheses and ancillaries is part of a complete system and must only be used with the corresponding elements from the original COUSIN BIOTECH system.

IntraSPINE® is made with soft materials, any contact with harder one may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials.

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). However, the implant is made with Polyethylene terephthalate, Dimethyl siloxane and Titanium which are non-ferromagnetic in these quantities according to literature. Consequently, the IntraSPINE® is compatible and safe with MRI.

Use of the system of COUSIN BIOTECH ancillaries and specific trial prostheses is required for implantation of the IntraSPINE® prosthesis. Implants must be handled and stored with great care, in a dry place, protected from light and at room temperature.

Transportation: Transportation of the devices should be done in the closed vehicles in accordance with the rules for those transport vehicles. There are no special transportation requirements for the device.

IMPORTANT : DO NOT REUSE AND DO NOT RESTERILISE THE INTRASPINE® PROSTHESIS

As indicated on the product labelling, the prosthesis is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be done following recommendations of ISO 12891-1:2011 « Implants for surgery – Retrieval and analysis of surgical implants » Part 1: « Retrieval and Handling ».

Any explanted device must be sent back, for analysis, following the current protocol. This protocol is available on demand to COUSIN BIOTECH. It is important to note that any implant that should not have been cleaned and disinfected before expedition must be contained in a sealed package.

The elimination of explanted medical device must be conducted in accordance with standards in the country for the disposal of infectious hazards waste.


The elimination of a non-implanted device is not the subject of specific recommendations.

WARNING ABOUT ANCILLARIES AND TRIAL PROSTHESES

The IntraSPINE® ancillaries are Class I medical devices, intended for temporary use and are re-usable.

The IntraSPINE® trial prostheses are Class IIa medical devices, intended for temporary use and are re-usable.

Unlike the IntraSPINE® prosthesis, ancillaries and trial prostheses (included with the instrument kit/set) are furnished unsterile.

KITACINTRA

Before first use

Ancillaries and trial prostheses must be:

- Cleaned in the appropriate way
- Sterilised in an autoclave

After use, cleanse and decontaminate the IntraSPINE® ancillaries and the trial prostheses according to a suitable protocol (reduction of the risk of transmission of unconventional transmissible agents - UTA).

In case of a patient suspected of having Creutzfeld-Jacob Disease (CJD), perform suitable decontamination.

In the event of a confirmed case of CJD, incinerate the ancillaries and the trial prostheses.

STERILITY OF INTRASPINE® ANCILLARY ELEMENTS

Sterilization and the prior steps such as cleansing and decontamination are performed under the responsibility of the health care centers; prescribed methods and materials must be used (Reduction of risk of transmission of UTA, French health ministry (DGS/R13/2011/449 dated 01/12/2011)

COUSIN BIOTECH recommends sterilization by steam autoclave. Ancillaries and the trial prostheses can be sterilized in a steam autoclave at 134°C for 18 minutes. Always verify the integrity and the correct operation of the ancillaries and the trial prostheses after each cleaning/sterilization cycle and before each usage.

STORAGE PRECAUTIONS

Ancillaries and trial prostheses must be stored in their dedicated container or equivalent packaging to prevent any deterioration.

Keep in a dry place away from light, at room temperature, in its original packaging.

INFORMATION REQUEST AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible.

In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time.

For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor, or COUSIN BIOTECH at contact@cousin-biotech.com.

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