PREMIUM®

NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT

en	Instructions for use	Page	2
fr	Notice d'instructions	Page	4
de	Gebrauchsanweisung	Seite	6
it	Istruzioni per l'uso	Pagina	8
es	Instrucciones de uso	Pagina	10
pt	Nota de instruções	Pagina	12
ελ	Οδηγίες χρήσης	Σελίδα	14
fi	Käyttöohjeet	Sivu	16
nl	Gebruiksaanwijzing	Pagina	18
ro	Instrucțiuni de Utilizare	Pagină	20

COUSIN BIOTECH

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NOT221 / 170519

Version du 19/05/2017

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

PREMIUM®

NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT

SINGLE USE STERILE PRODUCT

DESCRIPTION

PREMIUM® meshes are non resorbable parietal reinforcement implants.

IMPLANTED MATERIALS

• Mesh : polypropylene (knitted monofilament). Blue thread (depending on references): monofilament polypropylene, non absorbable, colorant [phthalocyaninato(2-)] copper

Origin neither human nor animal – Non resorbable.

INDICATIONS

Repair of hernia or other fascial defects that require the addition of a reinforcing or bridging material to obtain the desired surgical result.

PERFORMANCE

Meshes are particularly adapted either to laparoscopic surgery, or to conventional surgery.

The PREMIUM® is designed for extraperitoneal implantation only.

Meshes are biocompatible and non resorbable, They have the advantage of having a shape memory, being cut-able, a very high resistance to suture, having a great flexibility and allow an optimum integration and colonisation.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects:

- | | | | |
|---------------------|----------------|----------------------|-----------|
| - Discomfort / pain | - Recurrence | - Adhesion formation | - Erosion |
| - Infection | - Inflammation | - Fistsula formation | - Seroma |

USAGE PRECAUTIONS

PREMIUM® meshes are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. Mesh fixation points should be at least 1cm from the edge of the mesh with 1cm spacing between fixation points.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labeling, the PREMIUM® mesh is for single use only. It can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

STORAGE

To be stored in a dry place away from light and at room temperature in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

DESCRIPTION OF THE SURGICAL PROCEDURES:

Open Inguinal hernia repair : Lichtenstein

- The patient is placed in the supine position
- Local or general anaesthesia
- Iliac incision of approximately 5 cm
- Dissection is performed accordingly to the Lichtenstein technique
- The hernia sac is identified and reduced
- The transversalis fascia can be tightened with a non-absorbable suture. With direct hernias it helps to avoid any immediate post-operative recurrence
- PREMIUM® mesh can be cut to fit the patient's anatomy
- The prosthesis is progressively put in place
 - begin by placing PREMIUM® on the pubic spine
 - the upper part of the prosthesis is placed on the conjoint tendon
 - carry placing the prosthesis towards the deep inguinal ring so that the first arm is placed behind
 - then the second arm is placed so that it encircles the spermatic cord
 - Fix the PREMIUM® mesh with non-resorbable sutures
 - apex is sutured to the pubic tubercle
 - the lower border of the mesh is sutured to the free edge of the inguinal ligament
 - continuous suture extends up just medial to the anterior superior iliac spine
 - suture the two tails of the mesh together around the spermatic cord
 - the infero-medial corner of the mesh is attached overlapping the pubic tubercle.
 - the mesh is anchored to the conjoined tendon by interrupted sutures

- the cord then rests on the prosthesis
- suction drain can be placed beneath the external oblique aponeurosis, especially in large inguinal hernias, where an extensive dissection was performed
- the aponeurosis of external oblique is then closed with absorbable sutures
- before closure of the surgical incision, its edges are infiltrated with a long-acting local anesthetic, such as Naropein®
- During peri-operative care of the patient, prophylactic antibiotics is usually given for 48 – 72 hours postoperatively

Laparoscopic procedure (TAPP or TEP)

1/ Example of a TAPP repair

- Local or general anaesthesia
- A pneumoperitoneum is created in the usual fashion
- intra-abdominal findings are reported [intra-abdominal pathology and inguinal hernia defects and sacs].
- The two additional Trocars are inserted under direct vision.
- Dissection is initiated and carried out following classical TAPP technique
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected.
- The PREMIUM® mesh is rolled like a cigarette and inserted through the Trocar into the intra-abdominal cavity and deployed over the inguinal region. The blue thread on the oval mesh is the anatomic landmark for the pubic bone. For the meshes with flaps, the cross shows the upper flap. The PREMIUM® mesh is attached to Cooper's Ligament, around and lateral to the Inferior Epigastric Vessels using tacks, staples or other mesh fixing device, according to the Surgeon's usual technique. The operator should be meticulous to avoid the triangle of doom, triangle of pain, iliac vessels and to place the fixations lateral to the inguinal ring.
- The operator should check the PREMIUM® mesh is well anchored to the surrounding structures.
- The peritoneum is closed meticulously and no defect between the peritoneum and the abdominal wall should be left open. In addition, it should cover the entire Mesh.
- The trocars are removed under direct vision. The fascia of the sub-umbilical trocar site is closed as needed.

2/ Example of a TEP repair

- General anaesthesia
- Create the pro-peritoneal space.
- The 10 mm trocar is inserted.
- Create the pneumo-pro-peritoneum.
- Two 5 mm ports are inserted under direct vision.
- If a peritoneal defect is created at this time it should be immediately closed.
- The anatomical landmarks are identical to when performing a TAPP repair.
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected. This is done by gentle traction with atraumatic graspers.
- The PREMIUM® mesh is rolled like a cigarette and inserted into the pro-peritoneal space and deployed over the inguinal region. The blue thread on the oval mesh is the anatomic landmark for the pubic bone. For the meshes with flaps, the cross shows the upper flap.
- The PREMIUM® mesh is first attached or secured to Cooper's Ligament, and the superior aspect of the pubic ramus. The PREMIUM® mesh is then fixed on the posterior aspect of the Linea Alba. The anchoring is continued around and lateral to the Inferior Epigastric Vessels. A few fixations are used to fix the PREMIUM® mesh lateral to the internal ring.
- The PREMIUM® mesh should gently and generously cover the Iliac Vessels without major gaps.
- The operator should check the Mesh is well anchored to the surrounding structures. The Mesh should not migrate and remain in place.
- The pro-peritoneal space will be checked for any peritoneal defects. If any, they should be closed. The 5 mm trocars are removed under direct vision. The pro-peritoneal space will collapse.
- If there are any questions about a missed peritoneal defect, a completion laparoscopy should be performed.
- Trocars are removed and the fascial defect closed with the appropriate suture. The skin edges are approximated in the usual manner.

Open Extra-Peritoneal Ventral Hernia repair

Example: Retrorectus mesh repair: Rives-Stoppa technique.

This technique utilizes the hernia sac to separate the mesh from the intra-abdominal contents. Superior to the umbilicus, dissection is performed above the posterior rectus fascia and under the rectus muscle.

Below the umbilicus, dissection occurs in the preperitoneal space due to the lack of a

posterior rectus sheath. The PREMIUM® mesh is placed in the newly formed pre-peritoneal space, and fixated to the muscle layer above. The peritoneum must be meticulously checked and no defect in the peritoneum should be left open to prevent any contact between the mesh and intra-abdominal contents. Redon's drains under aspiration are usually placed.

The anterior sheath is closed. If tension is present, relaxing incisions may be used.

Skin closure. Using the Redon's drains, fill the pre-peritoneal space with an antiseptic solution and clamp the drains for 1 hour. A belt can be used to wrap the patient's abdomen for a few days to reinforce the abdominal wall.

Description of the prosthesis:

PREMIUM® meshes are constructed from knitted polypropylene monofilament with a diameter of approximately 0.1 mm.

PREMIUM® meshes are specifically indicated for laparoscopic and open surgery.

PREMIUM® exists in different sizes and shape, it is adapted to surgeon's technique and to patient's anatomy.

Size of the prostheses:

The PREMIUM® product range is composed of meshes of various sizes and shapes :

- Pre-cut meshes with keyhole for LICHTENSTEIN inguinal hernia repair
- Oval, pre-cut and rectangular meshes for TAPP/TEP inguinal hernia repair
- Square and rectangular meshes for extra-peritoneal open ventral hernia repairs

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