BIOLIG[®]

DEVICE FOR THE INTER-SPINOUS SPACE WITH SPINOUS SUPPORT

| | | | |
|----|---------------------------------------|--------|---|
| en | Instructions for use | Page | 2 |
| fr | Notice d'instructions | Page | 3 |
| de | Gebrauchsanweisung | Seite | 4 |
| it | Istruzioni per l'uso | Pagina | 5 |
| es | Instrucciones de uso | Pagina | 6 |
| pt | Nota de instruções | Pagina | 7 |
| cs | Návod k použití | Strana | 8 |

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NOT075/160531

Version du 31/05/2016

Instructions for use**BIOLIG®****DEVICE FOR THE INTER-SPINOUS SPACE WITH SPINOUS SUPPORT
SINGLE USE STERILE PRODUCT****DESCRIPTION**

BIOLIG® devices are ligaments for the inter-spinous space with spinous support

BIOLIG® range offers two possibilities for fastening system:

- Ligament with needle and titanium loop
- Ligament with needle, hardened loop and titanium crimp

IMPLANTED MATERIALS

Polyethylene terephthalate + Titanium + Barium sulfat + Polyvinyl chloride

(See on package)

MATERIALS IN CONTACT WITH THE PATIENT DURING THE SURGERY

Stainless steel,

Origin neither human nor animal – Not resorbable.

INDICATIONS

BIOLIG® devices for the interspinous space with spinous support are recommended for:

- Dorso lumbar instability
- Vertebral osteosynthesis ligament
- Spine surgery

PERFORMANCES

The main function of **BIOLIG®** ligament is to assist the posterior joint ligament during flexion.

Products may contain a radio-opaque wire (Barium sulfate, polyethylene terephthalate, and polyvinyl chloride) (see packaging).

CONTRAINDICATIONS

1) Do not use in the following cases:

- | | |
|------------------------------------|--------------------|
| - Allergy to any of the components | - Growing children |
| - Infected site | - Pregnancy |

2) Factors likely to adversely affect the success of the implantation:

- | | |
|---|---------------------------------|
| - Advanced osteoporosis | - Major spine deformity |
| - Local bone tumours | - Systemic or metabolic disease |
| - Infectious disease | - Obesity |
| - Addiction to drugs of abuse or medicinal products. | |
| - Intense physical activity, for example participation in competitive sports or heavy manual work | |

UNDESIRABLE SIDE EFFECTS

The **BIOLIG®** devices for the interspinous space with spinous support are susceptible to involve undesirable side effects such as:

- | | |
|--|--|
| - Infection of the prosthesis | - Rupture of ligament |
| - Cardiovascular problems | - Permanent elongation of the ligament |
| - Venous thrombosis and pulmonary embolism | - Hypersensitivity reaction |
| - Haematoma | |
| - Surgical revision with intersomatic arthrodesis in the operated area | |

PRECAUTIONS FOR USE

BIOLIG® devices for the interspinous space with spinous support are delivered sterile.

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the prosthesis has exceeded its use-by date.

Suture clips or sharp mounting screws must not be used.

This device can be used with other devices providing that the customer or manufacturer has tested the two products to ensure they are fully compatible.

This device must be implanted **only** by a qualified surgeon and formed with the use of the product (knowledge of anatomy and spinal surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As indicated on the product labelling, the prosthesis is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, additional surgery).

STORAGE OF THE PROSTHESIS

To be stored in a dry place away from light and at room temperature in its original packaging.

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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