BIOMESH SOFT PROLAPS'[®]

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NOT217/170904

Version du 04/09/2017

This release is the last update of the instructions of use and replace the previous edition.

Instructions for use

BIOMESH® SOFT PROLAPS' Implant for treatment of prolapse

DESCRIPTION

The system comprises a single-use, sterile surgical instrument, known as an ancillary, and a fabric implant comprising a central body with 4 extension arms. They are single use and delivered sterile.

It is designed for the treatment of female genital prolapse by vaginal approach.

- ◆ The arms of these prosthesis are extended by self-tightening loops for their placing.
- ◆ The distal/perineal arms (white threads) and the two proximal/sacro-spinal/sacro-sciatic arms (yellow threads) are differentiated by means of the 2 colours of threads.

The end of the ancillary needle has a slot for fixing self-tightening loops.

When placing the self-tightening loop in the slot, always make sure it is as tight as possible and that the thread is in the slot designed for this purpose so that the prosthesis does not come apart from the ancillary when passing in the patient's tissues. (See figure 1 below)


Fig 1

MATERIALS

- **Prosthesis:** polypropylene (knitted monofilament) with dimethylsiloxane-based biocompatible ink
 - **Thread:** Polyethylene terephthalate + biocompatible dye + dimethyl siloxane
 - **Ancillary device:** Stainless steel needle, polyoxymethylene handle
- Origin neither human, nor animal _ Non resorbable.

INDICATIONS

BIOMESH® SOFT PROLAPS' implants are designed to be used for support purposes in the surgical treatment of posterior genital prolapse via a vaginal approach after failure of a previous surgical treatment or if a particular clinical element concerns a high risk of recurrence.

PERFORMANCE

The biocompatible and non biodegradable product, produces a reaction resulting in fibrosis which assumes a reinforcing role.

These prosthesis are easy to insert due to their specially designed shape, excellent flexibility and optimum, rapid colonization and integration.

CONTRAINDICATIONS

Do not use the implant in the following circumstances:

- Allergy to any of the components
- In a place with progressive or latent infection or signs of tissue necrosis
- For women receiving anticoagulant therapy
- Urinary or genital tract infection or non-aseptic conditions
- Pregnancy
- Growing children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

SIDE EFFECTS

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection.

The known risks of surgery for the treatment of prolapse include acute or chronic pain, scarification of the vagina, retraction of the mesh, infection, erosion, urgency of novo, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

POSSIBLE UNDESIRABLE EFFECTS

- Acute or chronic pain
- Urination difficulties
- Constipation
- Hematoma
- Infection
- Erosion
- Dyspareunia
- Disorders of micturition
- Overactive bladder, urgency
- De novo SUI
- Prolapse in other compartment
- Rectal and vesical perforation
- Perforation or damage to blood vessels, bladder, rectum, intraperitoneal organs or nerves
- As with all implants, local irritation and/or a reaction to the foreign body may occur.

OPERATING TECHNIQUE

(Detailed operating technique available on request)

- Surgery under local, regional or general anesthesia.
- Incision: colpotomy with appropriate technique.
- Dissection: dissection of prolapse with appropriate technique.
- Positioning of the device and possible fixing.
- Closure of vagina with absorbable suture using appropriate technique.

PRECAUTIONS FOR USE

BIOMESH® SOFT PROLAPS' devices must **always** be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in gynaecological and urological surgery) who has a good knowledge of the device, its intended use, the instruments and the surgical technique.

BIOMESH® SOFT PROLAPS' prosthesis and ancillary are delivered sterile (sterilised using ethylene oxide). Inspect the packaging to be sure it is intact (do not use if protective packaging is damaged). Do not use if the device is damaged or out of date.

Cousin Biotech advises surgeons to follow the best practices for using reinforcement implants for the treatment of vaginal prolapse. In France, the French Health Authority specifies that the use of a surgical implant to treat vaginal prolapse is favourable in certain situations, notably:

- as a second-intention treatment after the failure of previous surgery
- if a specific clinical element suggests a high risk of relapse

IMPORTANT: DO NOT REUSE – DO NOT RESTERILIZE

As indicated on the product labelling, the prosthesis and ancillary are designed for single use. They should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse)

STORAGE OF THE PROSTHESIS AND ANCILLARIES

To be stored in a dry place away from light, at room temperature, in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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