BIOMESH SOFT PROLAPS[®]

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NOT120 / 170315

Version du 15/03/2017

BIOMESH® SOFT PROLAPS' Implant for treatment of prolapse

DESCRIPTION

BIOMESH® SOFT PROLAPS' prostheses consist of knitted polypropylene monofilament.

BIOMESH® SOFT PROLAPS' prostheses are available as rectangular meshes or precut meshes. There are intended for the treatment by a vaginal approach of the pelvic prolapse of the women.

- ◆ This prosthesis is equipped with traction threads on the tips of 4 arms to facilitate insertion and positioning.
- ◆ The main body of the prosthesis is marked with a black line to view for centering.
- ◆ The left and right arms are marked with 2 different shapes (circle, triangle) for differentiation. Top and bottom arms are differentiated by the number of markings (1 or 2).

COUSIN BIOTECH *INTRAMESH LIFT*® (FANCOFIL, FANCOFEN) instruments can be used to implant *BIOMESH® SOFT PROLAPS'* devices. Refer to the *INTRAMESH LIFT*® instructions for further information.

MATERIALS

100% polypropylene: single thread knitted, weight 38 g/m².

Traction threads: Polyethylene terephthalate + biocompatible dye

Dimethyl siloxane-based biocompatible ink

Origin neither human, nor animal _ Non resorbable.

INDICATIONS

BIOMESH® SOFT PROLAPS' implants are designed to be used for support purposes in the surgical treatment of anterior genital prolapse via a vaginal approach after failure of a previous surgical treatment or if a particular clinical element concerns a high risk of recurrence.

PERFORMANCE

The biocompatible and non biodegradable product, produces a reaction resulting in fibrosis which assumes a reinforcing role. These prostheses are easy to insert due to their specially designed shape, excellent flexibility and optimum, rapid colonization and integration.

CONTRAINDICATIONS

Do not use the implant in the following circumstances:

- Allergy to any of the components
- In a place with progressive or latent infection or signs of tissue necrosis
- For women receiving anticoagulant therapy
- Urinary or genital tract infection or non-aseptic conditions
- Pregnancy
- Growing children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

SIDE EFFECTS

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection.

The known risks of surgery for the treatment of prolapse include acute or chronic pain, scarification of the vagina, retraction of the mesh, infection, erosion, urgency of novo, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

POSSIBLE UNDESIRABLE EFFECTS

- Acute or chronic pain
- Urination difficulties
- Constipation
- Hematoma
- Infection
- Erosion
- Dyspareunia
- Disorders of micturition
- Overactive bladder, urgency
- De novo SUI
- Prolapse in other compartment
- Rectal and vesical perforation
- Perforation or damage to blood vessels, bladder, rectum, intraperitoneal organs or nerves
- As with all implants, local irritation and/or a reaction to the foreign body may occur.

RETRIEVAL AND DISPOSAL OF DEVICES

Explantation and handling should be done according to the recommendations of ISO 12891-1:2011 « Implants for surgery – Retrieval and analysis of surgical implants » Part 1: « Retrieval and Handling ».

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The disposal of the removed medical device must be done in accordance with the standards of the country concerned for the disposal of infectious waste.

The disposal of a non-implanted device is not the subject of specific recommendations.

OPERATING TECHNIQUE

(Detailed operating technique available on request)

- Surgery under local, regional or general anesthesia.
- Incision: colpotomy with appropriate technique.
- Dissection: dissection of prolapse with appropriate technique.
- Positioning of the device and possible fixing.
- Closure of vagina with absorbable suture using appropriate technique.

PRECAUTIONS FOR USE

BIOMESH® SOFT PROLAPS' devices must **always** be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in gynaecological and urological surgery) who has a good knowledge of the device, its intended use, the instruments and the surgical technique.

BIOMESH® SOFT PROLAPS' prosthesis and ancillary are delivered sterile (sterilised using ethylene oxide). Inspect the packaging to be sure it is intact (do not use if protective packaging is damaged). Do not use if the device is damaged or out of date.

Cousin Biotech advises surgeons to follow the best practices for using reinforcement implants for the treatment of vaginal prolapse. In France, the French Health Authority specifies that the use of a surgical implant to treat vaginal prolapse is favourable in certain situations, notably:

- as a second-intention treatment after the failure of previous surgery
- if a specific clinical element suggests a high risk of relapse

IMPORTANT: DO NOT REUSE – DO NOT RESTERILIZE

As indicated on the product labelling, the prosthesis and ancillary are designed for single use. They should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse)

STORAGE

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

PACKAGING, STORAGE AND STORING OF ANCILLARIES

The ancillaries must be stored conditioned in order to avoid any deterioration.

Sterilization and the prior steps such as cleansing and decontamination are performed under the responsibility of the health care centers; prescribed methods and materials must be used (Reduction of risk of transmission of UTA, French health ministry (DGS/R13/2011/449 dated 01/12/2011))

COUSIN BIOTECH recommends a steam sterilization. The instruments are compatible with a steam treatment with 134°C during 18min.

Always make sure that ancillaries are complete and in good working order after each cleaning/ sterilization cycle and before use.

If there is a possibility that the patient has contracted Creutzfeldt-Jakob Disease (CJD), proceed with an appropriate cleaning and decontamination method.

If the patient has contracted CJD, the instruments must be incinerated.

INFORMATION REQUEST AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time .

For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor, or COUSIN BIOTECH at contact@cousin-biotech.com.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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