BIORING®

REPAIR KIT WITHOUT PORT

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	4
de	<u>Gebrauchsanweisung</u>	Seite	6
nl	<u>Gebruikinstructies</u>	Page	8

**COUSIN BIOTECH**

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Date de marquage CE: 02/2009
NOT249/161125
Version du 25/11/2016

Instructions for use

ADJUSTABLE GASTRIC BANDING BIORING®

REPAIR KIT WITHOUT PORT

SINGLE USE STERILE PRODUCT

DESCRIPTION

Kit including a catheter, a connector and a connection system.

The various parts of the kit were designed for fitting perfectly to the BIORING® adjustable gastric banding.

MATERIALS

The catheter is made of silicone (Dimethyl Siloxane) and baryum sulphate labelled with a biocompatible ink. The connection system is made of silicone (Dimethyl Siloxane) and baryum sulphate. The connector is made of stainless steel.

No human nor animal origin. Not resorbable

INDICATIONS

Repair Kit to be used when a BIORING® gastric banding has been implanted and in the event of peroperative and postoperative complications, repeated surgery on patients following complications (leakage, catheter damaged).

PERFORMANCES

The connection system secures the connection between the port and the catheter by an on-clipping system

Connector: When the catheter on the BIORING® adjustable gastric banding is damaged (cut, leakage), the connector and the catheter from the kit can be used to replace the damaged system.

Always refer to the instructions for the BIORING® adjustable gastric banding marketed by COUSIN BIOTECH.

The product must be implanted only by a qualified surgeon and well formed for the use of it.

CONTRA-INDICATIONS

Do not use the implant in the following circumstances:

- Allergy to any of the components
- Growing children
- Do not implant to fit with a gastric band that is not a BIORING® adjustable gastric banding
- Inflammatory disease or congenital or acquired anomalies of the gastrointestinal tract.
- Severe organic disease.
- Potential bleeding disorder of the upper gastrointestinal tract or patients who regularly or systematically take aspirin or non-steroidal anti-inflammatory drugs.
- Diagnosis or preexisting symptoms or family history of connective tissue autoimmune disease.
- Pregnancy.
- Patients with a drug or alcohol addiction.
- Patients who refuse to accept the dietary restrictions required by the procedure.
- Serious psychiatric disorders.

Possible Complications

Complications that may arise from the use of an adjustable gastric band include risks associated with the anaesthesia, the surgical technique, and the degree of intolerance of the patient to any foreign body implanted in the body.

Perioperative complications (due to the surgical technique):

- Lesions of the spleen
- Perforation of the stomach or oesophagus
- other bleeding

In the immediate postoperative period:

- Nausea, vomiting (especially in the first few days after the intervention, and when the patient eats more than recommended).
- Bloating due to gas
- Dysphagia, dehydration
- Thromboembolic complications
- Infectious complications
- Postoperative bleeding

Delayed postoperative complications:

- Dilation of the upper pouch
- Slippage of the band
- Gastroesophageal reflux
- Pyrosis
- Gastritis of the upper pouch
- Migration of the band
- Displacement of the injection port
- Constrictive obstruction (due for example to oedema, food, an inadequate calibration, slippage of the band, or to torsion or dilation of the upper pouch)
- Deflation of the band (due to a leak in the band, injection port, or connection tube)
- Rapid weight loss may lead to symptoms of malnutrition, anemia, and associated complications
- Infection of the injection port
- Constipation

PRECAUTION FOR USE AND FOR STORAGE

- The components of the implantable port and of the connecting system have to be used only with the BIORING® adjustable gastric banding.
- The band can be inflated immediately using physiological saline (NaCl 0.9%) The volume, as well as the interval before filling, is at the discretion of the practitioner/surgeon.
- Inflating of the band: maximum volume of physiological saline according to sizes of the bands

Band sizes (diameter : mm)	Maximum volumen of physiological saline (cm3)
Ø 21/24	10
Ø 27/29	10
Ø 36/38	15

- To keep in a dry place without any light, and in an ambient atmosphere, in its packaging of origin
- Please check the perfect integrity of the packaging (do not use in case of damaged packaging protection and/or damage of the label)
- Do not use in case of damaged product or if the product has ran out date
- **DO NOT RE-USE** **DO NOT RESTERILIZE**

As specified on the product's labeling, the BIOMESH® mesh is for single use only. It can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

INFORMATION REQUEST AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible.

In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time .

For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor, or COUSIN BIOTECH at spine@cousin-biotech.com.

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