INTRAMESH® T1

NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT (INTRAPERITONEAL)

en	<u>Instructions for use</u>	page	2
fr	<u>Notices d'instructions</u>	page	3
de	<u>Gebrauchsanweisung</u>	Seite	4
it	<u>Istruzioni per l'uso</u>	Pagina	5
es	<u>Instrucciones de uso</u>	Pagina	6
pt	<u>Nota de instruções</u>	Pagina	7
ελ	<u>Οδηγίες χρήσης</u>	σελίδα	8
bg	<u>Инструкции за употреба</u>	страни	9
da	<u>Brugsvejledning</u>	Side	10
fi	<u>Käyttöohjeet</u>	Sivu	11
hu	<u>Felhasználási útmutató</u>	Pagina	12
nl	<u>Gebruiksaanwijzing</u>	Pagina	13
pl	<u>Instrukcja użycia</u>	Strona	14
ro	<u>Instrucţiuni de Utilizare</u>	Pagină	15
sk	<u>Návod k použiti</u>	Strana	16
cs	<u>Pokyny k použití</u>	Strana	17
tu	<u>Kullanım Kilavuzu</u>	Sayfa	18

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Made in France

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NOT009_171130
Version du 30/11/2017

This release is the last update of the instructions of use and replace the previous edition

**IMPLANT INTRAMESH® T1
NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT
SINGLE USE STERILE PRODUCT**

DESCRIPTION

INTRAMESH® T1 implants are parietal reinforcement implants.

IMPLANTED MATERIALS

Polypropylene and ePTFE for anti-adherent side (See on the package).

Biocompatible ink with dimethyl siloxane (See on the package).

Origin neither human nor animal – Non resorbable.

INDICATIONS

Repair and parietal reinforcement of umbilical and ventral hernias.

PERFORMANCE

These products can be used in laparoscopic surgery or open surgery. INTRAMESH® T1 products have a no-adhesive side, which is designed to reduce visceral adherences during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect. Marking in black ink, composed of a cross and dots is on the non-adhesive side. Meshes are biocompatible and non resorbable. They have the advantage of having a shape memory, a great flexibility and allow an optimum and rapid integration and colonisation.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects:

- Discomfort
- Infection
- Recurrence
- Adhesions

USAGE PRECAUTIONS

INTRAMESH® T1 implants are delivered sterile. Before any use, inspect the integrity of packaging and device (of which blister / pelables sachets). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a qualified surgeon trained to the use of the product (knowledge of anatomy and visceral surgery).

As indicated on the product labelling, the implant, is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of products sterility, risk of infection, loss of effectiveness, recurrence).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

Do not cut out the mesh and suture it at least 1cm of the edge. The trocar to be used must be chosen according to the size of the implant to be used. This is to allow passage without excessive force. Always roll the implant with the ePTFE side inside for the passage in the trocar. For use in keyhole surgery, we recommend the use of non-notched atraumatic forceps.

STORAGE

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor

IMPLANT INTRAMESH® T1
IMPLANT DE RENFORCEMENT PARIETAL NON RESORBABLE
PRODUIT STERILE A USAGE UNIQUE

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