# BIOMESH<sup>®</sup> SOFT PROLAPS'

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**NOT157/180123**

Version du 23/01/2018

This release is the last update of the instructions of use and replace the previous edition

## BIOMESH® SOFT PROLAPS'

### DESCRIPTION

*BIOMESH® SOFT PROLAPS'* prostheses consist of knitted polypropylene monofilament.

*BIOMESH® SOFT PROLAPS'* prostheses are available as rectangular meshes or precut meshes. There are intended for the treatment by a vaginal approach of the pelvic prolapse of the women.

- ♦ *BIOMESH® SOFT PROLAPS'* prostheses delivered as rectangular meshes, are supplied to be trimmed by the surgeon to the right shape and size for the patient's problem and anatomy, as well as the surgeon's personal surgical technique. In case, a pattern can be provided (according to the references), in order to facilitate cutting with the desired form. This drawing is only the indicative shape of the mesh to be cut out, the mesh carried out remains entirely under the responsibility of the surgeon.
- ♦ Certain precut prostheses of the range just as the rectangular meshes, can present traction threads assembled at the ends of certain arms allowing their installation and the adjustment of their position. These prostheses are also delivered with free traction threads, which can be assembled by the surgeon with the arms of the prosthesis at the operational time which is appropriate to him. In this case, the threads of traction comprise a loop tightened by a slip knot in which it is enough to make pass the arms of the prosthesis and to tighten the node in order to block the thread of traction on the arm of the prosthesis.

COUSIN BIOTECH *INTRAMESH LIFT®* (FANCOFIL, FANCOFEN) instruments can be used to implant *BIOMESH® SOFT PROLAPS'* devices. Refer to the *INTRAMESH LIFT®* instructions for further information.

### MATERIALS

100% polypropylene: single thread knitted, weight 38 g/m<sup>2</sup>.

Traction threads (if applicable): PET

Origin neither human, nor animal \_ Non resorbable.

### INDICATIONS AND PERFORMANCES

#### INDICATIONS

*BIOMESH® SOFT PROLAPS'* implants are designed to be used for support purposes in the surgical treatment of genital prolapse via a vaginal approach.

#### PERFORMANCES

The biocompatible and non biodegradable products, produce a reaction resulting in fibrosis which assumes a reinforcing role six months after implantation. The features of these implants are their insertability with easy cutting, high resistance to suture, excellent flexibility, and optimum and rapid colonization and integration.

#### CONTRAINDICATIONS

Not to be implanted in the following circumstances:

- Allergy to any of the component materials
- Urinary tract infection
- Growing children
- For woman receiving anticoagulant therapy
- Pregnancy

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

#### SIDE EFFECTS

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection

Over-correction may lead to temporary or permanent obstruction of the lower urinary tract and to retention of urine.

The known risks of surgery for the treatment of incontinence include acute or chronic pain, scarification of the vagina, retraction of the mesh, infection, erosion, urgency of novo, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

#### OPERATING TECHNIQUE

(Detailed operating technique available on request)

- Surgery under local, regional or general anesthesia.
- Incision: colpotomy with appropriate technique.
- Dissection: dissection of prolapse with appropriate technique.
- Choice of prosthesis *BIOMESH® SOFT PROLAPS'*: once the dissected prolapse, choose a prosthesis according to the width of the prolapse to be treated and according to the practices of the surgeon + possible cutting of the prosthesis for a better adaptation
- Laying of the plate and possible fixing.
- Closure of vagina with absorbable suture using appropriate technique.

#### PRECAUTIONS FOR USE

*BIOMESH® SOFT PROLAPS'* devices must **always** be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in gynaecological and urological surgery) who has a good knowledge of the device, its intended use, the instruments and the surgical technique.

*BIOMESH® SOFT PROLAPS'* implants are delivered sterile (sterilised using ethylene oxide), inspect the packaging to be sure it is intact (do not use if the labels and/or protective packaging are damaged).

Do not use if the device is damaged or out of date.

Cousin Biotech advises surgeons to follow the best practices for using reinforcement implants for the treatment of vaginal prolapse. In France, the French Health Authority specifies that the use of a surgical implant to treat vaginal prolapse is favourable in certain situations, notably:

- as a second-intention treatment after the failure of previous surgery
- if a specific clinical element suggests a high risk of relapse

Do not reuse and do not resterilize the prosthesis.

#### STORAGE OF THE PROSTHESIS

To be stored in a dry place away from light, at room temperature, in its original packaging.

#### WARNINGS OF THE ANCILLARIES

Ancillaries are medical devices of Class I, intended for a temporary use and reusable. Clean and decontaminate ancillaries after use, employing an appropriate method (to reduce the risk of non-conventional transmissible agent transmission).

Sterilization as well as the preceding stages such as cleaning and decontamination, are under the responsibility of the centers of care which must use a material and methods validated (Reduction of the risk of transmission of ATNC, circular DGS/R13/2011/449 of the 01/12/2011).

COUSIN BIOTECH recommends a steam sterilization. The instruments are compatible with a steam treatment with 134°C during 18min. Always make sure that ancillaries are complete and in good working order after each cleaning/sterilization cycle and before use.

If there is a possibility that the patient has contracted Creutzfeldt-Jacob Disease (CJD), proceed with an appropriate cleaning and decontamination method.

If the patient has contracted CJD, the instruments must be incinerated.

#### IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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