SUBURETHRAL SUPPORT TAPE

en	<u>Instructions for use</u>	Page	2
	<u>Surgical Procedure</u>		13
fr	<u>Notice d'instructions</u>	Page	3
	<u>Technique Opératoire</u>		13
de	<u>Gebrauchsanweisung</u>	Seite	4
	<u>Operationstechnik</u>		13
it	<u>Istruzioni per l'uso</u>	Pagina	5
	<u>Tecnica Operatoria</u>		13
es	<u>Instrucciones de uso</u>	Pagina	6
	<u>Técnica Operatoria</u>		13
pt	<u>Nota de instruções</u>	Pagina	7
	<u>Técnica Operatoria</u>		13
ελ	<u>Οδηγίες χρήσης</u>	σελίδα	8
	<u>Χειρουργική Διαδικασία</u>		13
nl	<u>Gebruiksaanwijzing</u>	Page	8
	<u>Chirurgische procedure</u>		13
pl	<u>Instrukcja obsługi</u>	Strona	9
	<u>Technika operacyjna</u>		13
sk	<u>Návod k použiti</u>	Strana	10
	<u>Chirurgický postup</u>		13
cs	<u>Pokyny k použiti</u>	Strana	11
	<u>Chirurgický postup</u>		13

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NOT219/180123

Version du 23/01/2018

This release is the last update of the instructions of use and replace the previous edition

LIFT® Suburethral support tape

DESCRIPTION

The LIFT® kit comprises a sling and 3 ancillary devices (1 pair for the transobturator technique and 1 for the retropubic and suprapubic approach).

They are single use and delivered sterile.

The LIFT® prosthesis consists of a strip of knitted monofilament polypropylene, packaged in a plastic sheath.

It is prolonged by self-tightening loops for placing.

The ancillary devices have different shapes depending on the chosen approach: 1 curved needle for the retropubic and suprapubic approach, a pair of helical needles for the in-out and out-in transobturator approaches.

The end of the ancillary needles has a slot for fixing self-tightening loops.

A single-use, sterile guide is supplied with the implant to secure the passage of the ancillary device during the in-out transobturator approach.

When placing the self-tightening loop in the slot, always make sure it is as tight as possible and that the thread is in the slot designed for this purpose so that the prosthesis does not come apart from the ancillary when passing in the patient's tissues. (See figure 1 below)


Fig 1

MATERIALS

- **Prosthesis:** polypropylene (knitted monofilament)
 - **Thread:** Polyethylene terephthalate + biocompatible dye + dimethyl siloxane
 - **Ancillary device:** Stainless steel needle, polyoxymethylene handle.
 - **Guide:** Liquid crystal polymer + biocompatible dye.
- Not human nor animal origin – Non absorbable.**

INDICATIONS

The LIFT® prosthesis is designed for use as a suburethral support sling for the surgical treatment by vaginal approach in women.

PERFORMANCE

These biocompatible, not biodegradable prostheses cause reactive fibrosis that produces additional urethral support.

These prostheses are easy to insert due to their specially designed shape, excellent flexibility and optimum, rapid colonization and integration.

CONTRAINDICATIONS

Do not implant in the following cases:

- Allergy to any one of the components
- For women receiving anticoagulant therapy.
- Urinary tract infection
- Pregnant women
- Skeletally-immature children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

SIDE EFFECTS

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection or dyspareunia.

Over-correction may lead to temporary or permanent obstruction of the lower urinary tract and to retention of urine.

The known risks of surgery for the treatment of incontinence include acute or chronic pain, scarification of the vagina, retraction of the mesh, infection, erosion, urgency of novo, dyspareunia, migration of the device, the total failure of the intervention leading to relapse.

One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

USAGE PRECAUTIONS

The LIFT® device should be implanted by qualified surgeons only (knowledge of gynaecological and urological anatomy and surgery) who are extremely familiar with the device.

The LIFT® implant, guide and ancillary devices are delivered sterile (ethylene oxide sterilization). Check that the packaging is intact before using (do not use if the labels and/or protective packaging are damaged). Do not use if the device is damaged.

The plastic sheath that covers the sling must be removed and must not be implanted.

Do not reuse or resterilize the implant, the guide and ancillary devices.

IMPORTANT: DO NOT REUSE – DO NOT RESTERILIZE

As indicated on the product labelling, the implant, guide and ancillary devices are designed for single use. They should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

STORAGE

To be stored in a dry place away from light, at room temperature, in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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
Transobturator Methode

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Via transobturatoriz

Installation of the patient, anaesthesia and infiltration performed as specified on the page named «Preliminary Information 1-2-3»

0. Trans-obturator incisions of approximately 1 cm each, positioned at least 1 cm above the horizontal line passing by the urethral meatus at 1 cm from the ischiopubic branch.
1. Vaginal incision and urethro-vaginal dissection.
2. Passage of helioid ancillary by way of the trans-obturator (blue handle). (the tunnelling device is inserted into one of the trans-obturator incisions).
3.
 - Hooking of LIFT® bandage onto needle.
 - Passage of the tape; and heterolateral passage: performance of the 3 previous operations on the heterolateral side.
 - The tape is deployed by pulling on the thread. There should remain a space of about one centimetre between the tape and the urethra.
4.
 - Cutting of the tape, flush with the skin at the level of the incisions. Hemostasis is checked after tensioning of the tape.
 - Closure of the vaginal wall by a crossed suture with absorbable thread and closure of the cutaneous incision.

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