ResCUBE™

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NOT263 / 180112

Version du 12/01/2018

Caution: Federal law (USA) restricts this device to sale,
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

ResCUBE™

Ligament System for Spine Stabilisation
Single use sterile product

DESCRIPTION

The ResCUBE device is an implant with the purpose of correction and stabilisation of the spine during fusion of the instrumented levels.

It consists of:

- A braid for correction and stabilisation of the levels, through tension;
- A malleable lamella for guidance of the braid during surgery and attached to the braid through a thermal-retractable sleeve. The malleable metal lamella and the sleeve linking it to the braid are not intended for long term implantation and must be removed after placement of the implant.
- A connector holding the tension on the braid

MATERIALS

Braid	Polyethylene Terephthalate	For long term implantation
Connector	Titanium TA6V ELI	
Malleable lamella	Stainless Steel	Not for long term implantation
Shrink Sleeve	Polyolefin	

All used materials are non resorbable and not of human or animal origin.

INDICATIONS

ResCUBE is indicated for the following surgical treatments:

- Spinal trauma surgery that can be treated with sublaminar, interspinous or facet ligamentoplasty.
- Spinal degenerative pathologies, as an adjunct to intervertebral fusion
- Spinal reconstructive surgery incorporated into constructions with the aim of correcting spinal deformities such as scoliosis, kyphosis and spondylolisthesis.

CONTRA-INDICATIONS

The ResCUBE device should not be used in the following cases:

- General or local infection in the intervention area
- Local inflammation in the intervention area
- Pregnancy
- Immune deficiency
- Any pathology affecting bone metabolism, such as osteopenia, osteoporosis, metastasis, tumour... can compromise the expected mechanical support of the implant, and therefore, the success of the surgery
- Insufficient tissue or muscle coverage
- Suspected or confirmed sensitivity to at least one of the components
- Excess weight or obesity can lead to mechanical overload on the implant, compromising the success of the surgery
- Any mental state or activity level incompatible with the postoperative follow-up can compromise the success of the surgery
- Smoking patients are more likely to present a non-satisfactory bone fusion or non-fusion
- Any other contra-indication to instrumented surgery of the spine

The described contra-indications can be absolute or relative and must be evaluated by the surgeon prior to any implant decision.

POTENTIAL UNDESIRABLE SECONDARY EFFECTS

In addition to the risks linked to instrumented spine surgery, the following effects can be observed in rare cases:

- Unusual sensations, discomfort or pain due to the presence of the implant
- Migration, deformation or breakage of one or several components of the implant that can cause neurological, mobility complications or pain
- Infection, allergic or inflammatory reactions due to the components of the implant, or to the presence of micro particles (metallosis)
- Halting the growth of instrumented levels
- Partial or total loss of correction
- Intra-operative or post-operative lesions, fissures or fractures, in case of insufficient bone density
- Pseudo arthritis
- Pulmonary complications
- Death

In case of undesirable secondary effects, a second surgery may be necessary, in order to remove or replace one or several implants.

PRECAUTIONS FOR USE

General precautions

ResCUBE devices must be implanted by specialist spinal surgeons, trained in the surgical technique. The information contained in this leaflet is necessary but not sufficient to master the surgical technique.

IMPORTANT: DO NOT RE-USE - DO NOT RESTERILISE

ResCUBE implants are supplied sterile and for single use. They must not be resterilised or re-used (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

Implants must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature.

Pre-operative precautions

The implantation decision must be taken by the surgeon, after informed consent of the patient, considering the indications, contra-indications and potential undesirable side effects.

Before the surgery, it is recommended to check the compatibility of devices which will be associated (size, materials...).

Peroperative precautions

Before opening, it is recommended to check the packaging is untampered with. Any damage on packaging may jeopardise the implant sterility. If the packaging is damaged, the implant must not be used, or resterilised.

Any implant whose expiration date has passed must not be used.

After opening, it is recommended to check the implant integrity.

Any implant with a damaged surface must not be implanted.

The placement of the implant must be done with the appropriate ancillaries, following the surgical technique recommended by COUSIN BIOTECH.

The malleable metal lamella and the sleeve linking it to the braid are not intended for long term implantation and must be removed after placement of the implant.

Post-operative precautions

The patient must be informed of the post-operative care and of precautions.

RETRIEVAL AND DISPOSAL OF DEVICES

Retrieval and handling should be done in compliance with ISO 12891-1:2011 « Implants for surgery – Retrieval and analysis of surgical implants » Part 1: « Retrieval and Handling ».

Any removed device must be sent back for analysis, following the current protocol. The latter is available on request to COUSIN BIOTECH. It is important to note that any implant that cannot be cleaned and disinfected before despatch must be contained in a sealed package.

The disposal of removed medical devices must be conducted according to national standards for the disposal of infectious hazards waste.

The disposal of a non-implanted device is not subject to specific recommendations.

INFORMATION REQUESTS AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply high quality medical devices. However, if a health professional (client, user, prescriber ...) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible.

In case of dysfunction of an implant or if it contributed to causing serious adverse effects for the patient, the health centre must follow the legal procedures in its country, and inform COUSIN BIOTECH as soon as possible.

For any correspondence, please specify the reference, batch number, details of a reference and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

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e-ifu




www.cousin-biotech.com/ifu

IFU is available on our webpage for the corresponding product. You can view, save and print the IFU. Please choose the pages of the languages via the printing menu.

If requested, we will send you a free printed copy of the Instructions of Use in the language you require within 7 days.

You can use the order form provided on our webpage for this purpose or by mail at ifurequest@cousin-biotech.com

	en fr de it es pt el ru	Batch number Numéro de lot Chargenbezeichnung Numero di lotto Número de lote Número do lote Αριθμός παρτίδας Номер партии		en fr de it es pt el ru	Caution (See instructions for use) Attention, (voir notice d'instructions) Vorsicht (siehe Verwendungsanleitung) Attenzione (Vedi istruzioni per l'uso) Precaución (véanse las instrucciones de uso) Atenção (Ver instruções de utilização) Προσοχή (βλέπε οδηγίες χρήσης) Внимание (см. инструкцию по применению)		en fr de it es pt el ru	Manufacturer Fabricant Fabrikant Produttore Fabricante Fabricante Κατασκευαστής Изготовитель
	en fr de it es pt el ru	Reference on the brochure Référence du catalogue Referenz des Kataloges Codice del catalogo Referencia del catalogo Referência do catálogo Αναφορά καταλόγου Обозначение в каталоге		en fr de it es pt el ru	Use before: year, month and day Utiliser jusque : année, mois et jour Zu verwenden bis: Jahr, Monat und Tag Utilizzare entro il: anno, mese e giorno Utilizar hasta: año, mes y día Utilizar até: ano, mês e dia Χρήση έως: έτος, μήνας και ημέρα Использовать до: год, месяц, число		en fr de it es pt el ru	Do not re-use Ne pas réutiliser Nicht wiederverwenden Non riutilizzare No reutilizar Não reutilizar Μην επαναχρησιμοποιείτε Не использовать повторно
	en fr de it es pt el ru	Do not re-sterilize Ne pas restériliser Nicht wieder sterilisieren Non risterilizzare No esterilizar de nuevo Não reesterilizar Μην επαναποστειρώνετε Не стерилизовать повторно		en fr de it es pt el ru	Keep in a dry place Conserver au sec An einem trockenen Ort lagern Conservare in un luogo asciutto Conservar en un lugar seco Manter em local seco Διατηρείστε σε στεγνό μέρος Хранить в сухом месте		en fr de it es pt el ru	Non pyrogenic Non pyrogène Nichtpyrogene Apirogeno Απυρόγαινα Não pirogênica μη πυρετογόνο Непирогенный
	en fr de it es pt el ru	Keep away from sunlight A stocker à l'abri de la lumière du soleil Von direkter Sonneneinstrahlung schützen Tenere lontano dalla luce del sole Mantener alejado de la luz solar Manter afastado da luz solar Φυλάξτε μακριά από τον ήλιο Хранить вдали от солнечного света		en fr de it es pt el ru	Do not use if the packaging is damaged Ne pas utiliser si l'emballage est endommagé Nicht benutzen, wenn die Verpackung beschädigt ist Non utilizzare se l'imballaggio è danneggiato No utilizar si el envase está dañado Não usar se a embalagem estiver danificada Μην το χρησιμοποιείτε εάν η συσκευασία έχει υποστεί ζημιά Не использовать в случае повреждения упаковки			
	en fr de it es pt el ru	Connector only compatible with Xmm rod L = Left connector / R = Right connector Connecteur utilisable sur une tige de Xmm uniquement L = Connecteur gauche / R = Connecteur droit Dieses Verbindungsstück darf nur mit einem Xmm-Stift verwendet werden L = Linkes Verbindungsstück / R = Rechtes Verbindungsstück Connettore utilizzabile solo su un'asta da X mm L = Connettore sinistro / R = Connettore destro Conector utilizable en una varilla de Xmm solamente L = Conector izquierdo / R = Conector derecho Conector apenas utilizável numa haste de X mm L = Conector Esquerdo / R = Conector Direito Συνδετήρας που λειτουργεί μόνο με ράβδο Xmm L = Αριστερός συνδετήρας / R = Δεξιός συνδετήρας Соединительный элемент используется только на стержне X мм L = Соединительный элемент левый / R = Соединительный элемент правый						
	en fr de it es pt el ru	Sterilized by gamma radiation. Produit stérile. Méthode de stérilisation : irradiation Sterilisierungsmethode : Gammastrahlung. Il prodotto è sterile. Sterilizzazione : raggi Gamma. Esterilizado por rayos gamma. Produto estéril. Método de esterilização: irradiação Στείρο προϊόν. Μέθοδος αποστείρωσης: ακτινοβολία Стерильно. Технология стерилизации: иррадиация.						
	en fr de it es pt el ru	CE mark and identification number of Notified Body. Product conforms to essential requirements of the Medical device directive 93/42/EEC Marque CE et numéro d'identification de l'organisme notifié. Produit conforme aux exigences essentielles de la directive 93/42/CEE CE-Kennzeichnung und Identifikationsnummer der benannten Stelle. Das Produkt entspricht den Anforderungen der Richtlinie 93/42/EWG über Medizinprodukte Marchio CE e numero di identificazione dell'Organismo Notificato. Il prodotto è conforme ai requisiti essenziali della Direttiva CEE/93/42 Marca CE y número de identificación de Organismo Notificado. Este producto cumple los requisitos esenciales de la directiva 93/42/EEC Marca CE e número de identificação de Organismo Notificado. O produto obedece aos principais requisitos da Directiva 93/42/CEE Προϊόν που συμμορφώνεται με τις απαιτήσεις της οδηγίας 93/42/ΕΟΚ Знак CE и идентификационный номер уполномоченного органа. Изделие соответствует основным требованиям директивы 93/42/CEE						

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Caution: Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

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Single use NAJA™ Ancillary **NAJA™ ancillaries**

DESCRIPTION

This ancillary is intended for the placement of the NAJA™ ligament system for spinal stabilisation. This ancillary is a Class Is medical device, intended for single use and is sterile supplied.

MATERIALS

The different elements of the ancillary are made from:

- IXEF®

- Stainless steel (for metallic components)

All the materials used are non-absorbable and are neither of human nor animal origin.

PRECAUTIONS FOR USE

The ancillary must stay stored in its original packaging, which ensures safety.

Carry out a visual examination of each instrument before each use in order to detect any possible deterioration

The instruments must not be used except with elements corresponding to the original COUSIN BIOTECH system.

Use of specific ancillary COUSIN BIOTECH system is required for implantation of the implant.

In addition, the instruments must not be used except by a surgeon with appropriate training in the placement of NAJA™ implants.

Sterile device (Sterilized using Gamma irradiation)

Store the ancillary in its original protective packaging, and do not remove them from the packaging until immediately before use.

Prior to use, check the product expiration date and verify the integrity of the sterile packaging. Do not use, if the package is damaged.

Single-use device

According to the labelling of this product, the ancillary is a single use device.

Products intended for single-use must not be re-used. Re-use or reprocessing (e.g. cleaning and resterilization) may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness or death. Furthermore, reuse or reprocessing of single-use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user.

Contaminated ancillary must not be reprocessed. Any ancillary that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol. Even though it may appear undamaged, the ancillary may have small defects and internal stress patterns that may cause material fatigue.

Storage precautions of the ancillary

The device must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature in its original packaging.

INFORMATION REQUEST AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible.

In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time.

For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor, or COUSIN BIOTECH at spine@cousin-biotech.com.

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