# SOFT LIFT®

## SUBURETHRAL SUPPORT TAPE

---

en	Instructions for use	Page	2	
	<b>Surgical Procedure</b>		12	13 -14
fr	Notice d'instructions	Page	3	
	<b>Technique Opératoire</b>		12	13 -14
de	Gebrauchsanweisung	Seite	4	
	<b>Operationstechnik</b>		12	13 -14
it	Istruzioni per l'uso	Pagina	5	
	<b>Tecnica Operatoria</b>		12	13 -14
es	Instrucciones de uso	Pagina	6	
	<b>Técnica Operatoria</b>		12	13 -14
pt	Nota de instruções	Pagina	7	
	<b>Técnica operatória</b>		12	13 -15
da	Brugsdvejledning	Side	8	
	<b>Operativ teknik</b>		12	13 -15
lt	Naudojimo instrukcija	Puslapis	9	
	<b>Operatyvinė technika</b>		12	13 -15
nl	Gebruiksaanwijzing	Pagina	10	
	<b>Chirurgische procedure</b>		12	13 -15

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Date de marquage CE : Aout 2002  
**NOT131/180531**

Version du 31/05/2018

## Suburethral support tape SOFT LIFT® And tunnelling device associated (FANCOFIL, FANCOFEN)

### **DESCRIPTION**

The SOFT LIFT® device consists of a woven, single-thread polypropylene tape, extended at either end with traction treads for installation and peroperative tension adjustment.

The special COUSIN BIOTECH tunnelling device (FANCOFIL, FANCOFEN) has to be used to implant the SOFT LIFT® device.

A guide with implant is provided to guide the passage of the ancillary through the transobturator approach.

**Origin neither human, nor animal - Non Resorbable.**

### **MATERIALS**

• **Prosthesis:** Tape: 100% polypropylene (woven, single-thread)

Traction threads: Polyethylene terephthalate

• **Tunnelling device:** Stainless steel, acetal copolymer

• **Guide:** Liquid crystal polymer + biocompatible dye.

### **INDICATIONS**

The SOFT LIFT® prosthesis is designed for use as a urethra support tape for the surgical treatment of feminine stress incontinence.

### **CONTRA INDICATIONS**

Not to be implanted in the following circumstances:

- Allergy to any of the component materials
- For woman receiving anticoagulant therapy
- Urinary tract infection
- Pregnancy
- Growing children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

### **SIDE EFFECTS**

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection

Over-correction may lead to temporary or permanent obstruction of the lower urinary tract and to retention of urine.

The known risks of surgery for the treatment of incontinence include acute or chronic pain, scarification of the vagina, retraction of the mesh, infection, erosion, urgency of novo, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

### **WARNINGS**

The SOFT LIFT® device must **always** be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in gynaecological and urological surgery) who has a good knowledge of the device, its intended use, the instruments and the surgical technique.

The SOFT LIFT® device is delivered sterile (sterilised using ethylene oxide), inspect the packaging to be sure it is intact (do not use if the labels and/or protective packaging are damaged). Do not use if the device is damaged or out-of-date.

Do not reuse and do not resterilize the prosthesis.

Guides can only be used once and should not be implanted.

To extract the tape of the packaging :

Roll the exterior white rubber band to free the green threads. Then, pull gently one of the green threads so that it slides and comes out of the tube.

### **STORAGE**

To be stored in a dry place away from light, at room temperature, in its original packaging.

### **PACKAGING, STORAGE AND STORING OF TUNNELLING DEVICE SOFT LIFT®**

The tunnelling device must be stored conditioned in order to avoid any deterioration.

Sterilization as well as the preceding stages such as cleaning and decontamination, are under the responsibility of the centers of care which must use a material and methods validated (Reduction of the risk of transmission of ATNC, circular DGS/R13/2011/449 of the 01/12/2011).

COUSIN BIOTECH recommends a steam sterilization. The instruments are compatible with a steam treatment with 134°C during 18min.

Always make sure that ancillaries are complete and in good working order after each cleaning/ sterilization cycle and before use.

### **HANDLING OF TUNNELLING DEVICE SOFT LIFT®**

Before use, check that the packaging is not damaged in any way that might suggest compromise of the contents.

The tunnelling device must be stored in a container specifically designed for this use in equivalent packaging.

To prevent that the instruments are not in contact with objects which can deteriorate their surface.

Proceed to a visual examination of each instrument, in order to detect any possible deterioration.

The tunnelling device SOFT LIFT® are only suitable for use with original COUSIN BIOTECH instruments.

Moreover, the instruments must be used by a surgeon correctly trained to the implantation of prostheses SOFT LIFT®.

### **Use of tunnelling device SOFT LIFT® during the intervention**

Refer to the following operative technique instructions for further information.

SOFT LIFT® tunnelling device are designed for use in implantation of the SOFT LIFT® prosthesis, using the traction threads at the ends of the tape. During the procedure, these threads are attached to the ends of needles which are used to position the tape according to the pelvic anatomy prior to implantation of the device.

SOFT LIFT® tunnelling device come in different shapes according to the route of implantation: curved needles for the vaginal or supra-pubic route, and spiral needles for the in-out and out-in trans-obturator routes.

### **Warnings on the tunnelling device**

-Tunnelling device SOFT LIFT® are medical devices of Class I, intended for a temporary use or reusable.

-Clean and decontaminate tunnelling device SOFT LIFT® after use, employing an appropriate method (to reduce the risk of non-conventional transmissible agent transmission).

If there is a possibility that the patient has contracted Creutzfeldt-Jakob Disease (CJD), proceed with an appropriate cleaning and decontamination method.

If the patient has contracted CJD, the instruments must be incinerated.

### **IMPORTANT**

**For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.**

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