4D Mesh®

SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT

en	Instructions for use	Page	2
fr	Notice d'instructions	Page	4
de	Gebrauchsanweisung	Seite	6
it	Istruzioni per l'uso	Pagina	8
es	Instrucciones de uso	Pagina	10
pt	Instruções de utilização	Página	12
ελ	Οδηγίες χρήσης	Σελίδα	14
cs	Návod k použití	Strana	16
da	Brugsvejledning	side	18
fi	Käyttöohjeet	Sivu	20
hu	Használati Útmutató	Pagina	22
nl	Gebruiksaanwijzing	Pagina	24
ro	Instrucțiuni de utilizare	Pagină	26
sk	Návod na použitie	Strana	28
sv	Bruksanvisning	Sida	30
tu	Kullanım Talimatları	Sayfa	32

COUSIN BIOTECH

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NOT233_180627

Version du 27/06/2018

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

SEMI-RESORBABLE PARIETAL REINFORCEMENT IMPLANT STERILE PRODUCT FOR SINGLE USE

DESCRIPTION

4DMESH[®] meshes are semi-resorbable parietal reinforcement implants.

MATERIALS

Polypropylene (non resorbable) – Poly L Lactic Acid (resorbable). Blue thread (depending on references): monofilament polypropylene, non absorbable, colorant [phthalocyaninato(2-)] copper.

No human nor animal origin. Semi-resorbable.

INDICATIONS

4DMESH[®] meshes are designed for the repair and reinforcement of inguinal and femoral hernias.

PERFORMANCES

The device 4DMESH[®] is designed to be extra-peritoneally implanted. These products are biocompatible and semi-resorbable. Their main benefits are: can be cut to fit patient's anatomy, high resistance to suture and abdominal pressure, macroporous, transparent, lightweight after PLLA resorption.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
- Recurrence
- Adhesion formation
- Erosion
- Infection
- Inflammation
- Fistsula formation
- Seroma

PRECAUTIONS FOR USE

4DMESH[®] must always be separated from abdominal cavity by peritoneum

4DMESH[®] prosthesis are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labelling, the 4DMESH[®] mesh is for single use only. It can not be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

STORAGE OF THE PROSTHESIS

To be stored in a dry place away from light and at room temperature in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

SURGICAL PROCEDURE

DESCRIPTION OF THE SURGICAL PROCEDURES:

Open Inguinal hernia repair : Lichtenstein

- The patient is placed in the supine position
- Local or general anaesthesia
- Iliac incision of approximately 5 cm
- Dissection is performed accordingly to the Lichtenstein technique
- The hernia sac is identified and reduced
- The transversalis fascia can be tightened with a non-absorbable suture. With direct hernias it helps to avoid any immediate post-operative recurrence
- 4DMESH[®] mesh can be cut to fit the patient's anatomy
- The prosthesis is progressively put in place
 - begin by placing 4DMESH[®] on the pubic spine
 - the upper part of the prosthesis is placed on the conjoint tendon
 - carry placing the prosthesis towards the deep inguinal ring so that the first arm is placed behind
 - then the second arm is placed so that it encircles the spermatic cord
 - Fix the 4DMESH[®] mesh with non-resorbable sutures
 - apex is sutured to the pubic tubercle
 - the lower border of the mesh is sutured to the free edge of the inguinal ligament
 - continuous suture extends up just medial to the anterior superior iliac spine
 - suture the two tails of the mesh together around the spermatic cord
 - the infero-medial corner of the mesh is attached overlapping the pubic tubercle.
 - the mesh is anchored to the conjoined tendon by interrupted sutures
 - the cord then rests on the prosthesis
 - suction drain can be placed beneath the external oblique aponeurosis, especially in large inguinal hernias, where an extensive dissection was performed
 - the aponeurosis of external oblique is then closed with absorbable sutures
 - before closure of the surgical incision, its edges are infiltrated with a long-acting local anaesthetic, such as Naropin[®]
 - During peri-operative care of the patient, prophylactic antibiotics is usually given for 48 – 72 hours postoperatively

Laparoscopic procedure (TAPP or TEP)

1/ Example of a TAPP repair

- General anaesthesia
- A pneumoperitoneum is created in the usual fashion
- intra-abdominal findings are reported [intra-abdominal pathology and inguinal hernia defects and sacs].
- The two additional trocars are inserted under direct vision.
- Dissection is initiated and carried out following classical TAPP technique
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected.
- The 4DMESH[®] mesh is inserted through the Trocar into the intra-abdominal cavity and deployed over the inguinal region. The blue thread on the oval and preshaped mesh is the anatomic landmark for the pubic bone. For the meshes with flaps, the cross shows the upper flap. The 4DMESH[®] mesh is attached to Cooper's Ligament, around and lateral to the Inferior Epigastric Vessels using tacks, staples or other mesh fixating device, according to the Surgeon's usual technique. The operator should be meticulous to avoid the triangle of doom, triangle of pain, iliac vessels and to place the fixations lateral to the inguinal ring. Mesh fixation points should be at least 1cm from the edge of the mesh with 1 cm spacing between fixation points.
- The operator should check the 4DMESH[®] mesh is well anchored to the surrounding structures.
- The peritoneum is closed meticulously and no defect between the peritoneum and the abdominal wall should be left open. In addition, it should cover the entire Mesh.
- The trocars are removed under direct vision. The fascia of the sub-umbilical trocar site is closed as needed.

2/ Example of a TEP repair

- General anaesthesia
- Create the pre-peritoneal space.
- The 10 mm trocar is inserted.
- Create the pneumo-pre-peritoneum.
- Two 5 mm ports are inserted under direct vision.
- If a peritoneal defect is created at this time it should be immediately closed.
- The anatomical landmarks are identical to when performing a TAPP repair.
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected. This is done by gentle traction with atraumatic graspers.
- The 4DMESH[®] mesh is inserted into the pre-peritoneal space and deployed over the inguinal region. The blue thread on the oval and preshaped mesh is the anatomic landmark for the pubic bone. For the meshes with flaps, the cross shows the upper flap.
- The 4DMESH[®] mesh is first attached or secured to Cooper's Ligament, and the superior aspect of the pubic ramus. The 4DMESH mesh is then fixed on the posterior aspect of the Linea Alba. The anchoring is continued around and lateral to the Inferior Epigastric Vessels. A few fixations are used to fix the 4DMESH mesh lateral to the internal ring. Mesh fixation points should be at least 1cm from the edge of the mesh with 1 cm spacing between fixation points.
- The 4DMESH[®] mesh should gently and generously cover the Iliac Vessels without major gaps.
- The operator should check the Mesh is well anchored to the surrounding structures. The Mesh should not migrate and remain in place.
- The pre-peritoneal space will be checked for any peritoneal defects. If any, they should be closed. The 5 mm trocars are removed under direct vision. The pre-peritoneal space will collapse.
- If there are any questions about a missed peritoneal defect, a completion laparoscopy should be performed.
- Trocars are removed and the fascial defect closed with the appropriate suture. The skin edges are approximated in the usual manner.

Size of the prosthesis:

The 4DMESH[®] product range is composed of meshes of various sizes and shapes:

- Pre-cut meshes with keyhole for LICHTENSTEIN inguinal hernia repair
- Oval, pre-shaped and rectangular meshes for TAPP/TEP inguinal hernia repair

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