4D *Ventral*®

SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT EXTRAPERITONEAL

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	3
de	<u>Gebrauchsanweisung</u>	Seite	4
it	<u>Istruzioni per l'uso</u>	Pagina	5
es	<u>Instrucciones de uso</u>	Pagina	6
pt	<u>Instruções de utilização</u>	Página	7
ελ	<u>Οδηγίες χρήσης</u>	Σελίδα	8
bg	<u>Инструкции за употреба</u>	страница	9
cs	<u>Návod k použití</u>	Strana	10
da	<u>Brugsvejledning</u>	side	11
fi	<u>Käyttöohjeet</u>	Sivu	12
nl	<u>Gebruiksaanwijzing</u>	Pagina	13
pl	<u>Instrukcja obsługi</u>	Strona	14
ro	<u>Instrucțiuni de utilizare</u>	Pagină	15
sk	<u>Návod na použitie</u>	Strana	16

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Date de marquage CE: 04 / 2016

NOT245/180627
Version du 27/06/2018

Caution: Federal law (USA) restricts this device to sale,
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT
EXTRAPERITONEAL
SINGLE USE STERILE PRODUCT

DESCRIPTION

4D Ventral[®] meshes are extraperitoneal semi resorbable parietal reinforcement implants.

IMPLANTED MATERIALS

Polypropylene (non-resorbable knitted monofilament) – Poly L Lactic Acid (resorbable knit)

Origin neither human nor animal – Semi resorbable.

INDICATIONS

Repair of ventral hernia or other fascial defects that require the addition of an extraperitoneal reinforcing or bridging material to obtain the desired surgical result.

PERFORMANCE

The 4D Ventral[®] implant is designed for extraperitoneal implantation only.

Meshes are biocompatible and semi resorbable. They have the advantage of being cut-able, of having a very high resistance to suture, of having a great flexibility.

CONTRAINDICATIONS

Do not implant intraperitoneally.

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- | | | | |
|-------------------|----------------|----------------------|-----------|
| - Discomfort/Pain | - Recurrence | - Adhesion formation | - Erosion |
| - Infection | - Inflammation | - Fistula formation | - Seroma |

USAGE PRECAUTIONS

4D Ventral[®] implants are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. Mesh fixation points should be at least 1cm from the edge of the mesh with 1cm spacing between fixation points.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labeling, the 4D Ventral[®] implants are for single use only. They can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

STORAGE

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

SURGICAL PROCEDURE**Extra-Peritoneal Hernia repair**

Example: Retrorectus mesh repair of a ventral hernia (Rives-Stoppa technique).

This technique utilizes the hernia sac to separate the mesh from the intra-abdominal contents. Superior to the umbilicus, dissection is performed above the posterior rectus fascia and under the rectus muscle. Below the umbilicus, dissection occurs in the preperitoneal space due to the lack of a posterior rectus sheath. The 4D Ventral[®] implant is placed in the newly formed pre-peritoneal space, and fixated to the muscle layer above. The peritoneum must be meticulously checked and no defect in the peritoneum should be left open to prevent any contact between the mesh and intra-abdominal contents. Redon's drains under aspiration are usually placed.

The anterior sheath is closed. If tension is present, relaxing incisions may be used.

Skin closure. Using the Redon's drains, fill the pre-peritoneal space with an antiseptic solution and clamp the drains for 1 hour. A belt can be used to wrap the patient's abdomen for a few days to reinforce the abdominal wall.

Size of the meshes:

The 4D Ventral[®] product range is composed of meshes of various sizes and shapes:

- Round (12cm diameter) for extra-peritoneal hernia repairs
- Square and rectangular meshes for extra-peritoneal hernia repairs (various sizes available to fit with surgeon's needs)

IMPLANT DE RENFORCEMENT PARIETAL SEMI RESORBABLE
EXTRAPERITONEAL
PRODUIT STERILE A USAGE UNIQUE

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