4D DOME[®]

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NOT035 / 180627

Version du 27/06/2018

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

Instructions for use

4DDOME®

SEMI-RESORBABLE REINFORCEMENT IMPLANT SINGLE USE STERILE PRODUCT

DESCRIPTION

The device 4DDOME® consists of a dome and semi-resorbable reinforcement mesh for the anterior and posterior reinforcement.

MATERIAL

Polypropylene (non resorbable) Poly-L-Lactic Acid (resorbable)
Origin neither human nor animal origin. Semi-resorbable

INDICATIONS

Repair and surgical reinforcement of inguinal and crural hernias.

PERFORMANCES

These implants, biocompatible and semi-resorbable make a reactional fibrosis which replaces the reinforcement after an implantation of 6 months. They have the advantages to have a very strong resistance to the suture, a quick and optimum integration and colonization and to have a strong resistance to the abdominal pressure.

The 4DDOME® device is particularly designed to be implanted extra-peritoneally.

CONTRAINDICATIONS

Do not implant in the following cases:

- allergy to any of the ingredients
- septic environment
- pregnancy
- growing child

POSSIBLE UNDESIRABLE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable effects:

- Discomfort
- Infection
- Recurrence

PRECAUTIONS FOR USE

The dome should always be installed with an anterior reinforcement.

The dome should not be cut.

4DDOME® devices are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which blister / peel-away pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As indicated on the product labelling, the 4DDOME® implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

STORAGE

To be stored in a dry place away from sunlight, and at room temperature, in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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