BIOMESH[®] A2 & A4

NON RESORBABLE EXTRAPERITONEAL PARIETAL REINFORCEMENT IMPLANT

en	Instructions for use	Page 2
fr	Notice d'instructions	Page 3
de	Gebrauchsanweisung	Seite 4
it	Istruzioni per l'uso	Pagine 5
es	Instrucciones de uso	Pagina 6
ελ	Οδηγίες χρήσης	σελίδα 7
cs	Navod k použití	Strana 8
sk	Pokyny k použiti	Strana 9

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Made in France

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NOT069/180725
Version du 25/07/2018

This release is the last update of the instructions of use and replace the previous edition

BIOMESH A2 & A4

Extraperitoneal reinforcement implant

STERILE PRODUCT FOR SINGLE USE

DESCRIPTION

BIOMESH[®] meshes are parietal reinforcement implants.

MATERIALS

Polyethylene terephthalate (see on the package)

Polyethylene terephthalate deep-dyed with dimethyl siloxane (see on the package)

Marked thread : Polyethylene terephthalate tinted with a biocompatible dye (See on the package).

Origin neither human nor animal – Non resorbable.

INDICATIONS

Repair and parietal reinforcement for inguinal, crural and ventral hernias.

PERFORMANCES

Meshes are adapted to both laparoscopic and conventional surgery.

The BIOMESH[®] device is designed for extraperitoneal implantation.

Meshes, biocompatible and non resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a shape memory, being cut-able, have a very high resistance to suture, have a great flexibility and allow an optimum and rapid integration and colonisation.

CONTRAINDICATIONS

Do not use in the following cases :

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy
- Surgery via a vaginal approach

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to involve possible undesirable effects:

- Discomfort
- Infection
- Recurrence

PRECAUTIONS FOR USE

BIOMESH[®] meshes are delivered sterile.

Before any use, inspect the integrity of packaging and device (of which blister / pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a **qualified** surgeon and trained to the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT : DO NOT REUSE - DO NOT RESTERILIZE

As indicated on the product labelling, the implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

STORAGE OF THE PROSTHESIS

To be stored in a dry place away from light and at room temperature in its original packaging.

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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