# SACROMESH SOFT PROLAPS'<sup>®</sup>

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Made in France

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**NOT123 / 180904**  
Version du 04/09/2018

**SACROMESH SOFT PROLAPS®**  
**UROGENITAL IMPLANT**  
Treatment of genital prolapsus by promontofixation

**DESCRIPTION**

Sacromesh Soft Prolaps' consist of knitted polypropylene monofilament. There are intended for sacrocolpopexy of the pelvic prolapse of the women, by abdominal approach.

**MATERIALS**

100% polypropylene: single thread knitted ,weight 38 g/m<sup>2</sup> .  
Marked thread : Polyethylene terephthalate tinted with a biocompatible dye (See on the package).  
Origin neither human, nor animal \_ No resorbable.

**INDICATIONS**

Sacromesh Soft Prolaps' prosthesis are designed for use as reinforcement and support of tissues when treating genital prolapse surgically by sacrocolpopexy.

They can be used in either conventional or laparoscopic surgery.

**PERFORMANCES**

The biocompatible and non-resorbable prosthesis have the advantage of being easy to install, show good resistance to sutures and traction, with excellent flexibility, optimal and rapid integration and colonization.

**CONTRAINDICATIONS**

Do not use the implant in the following circumstances:

- Allergy to any of the components
- In a place with progressive or latent infection or signs of tissue necrosis
- For women receiving anticoagulant therapy
- Urinary or genital tract infection or non-aseptic conditions
- Pregnancy
- Growing children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

**SIDE EFFECTS**

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection.

The known risks of surgery for the treatment of prolapse include acute or chronic pain, scarification of the vagina, retraction of the mesh, infection, erosion, urgency of novo, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

**POSSIBLE UNDESIRABLE EFFECTS**

- Acute or chronic pain
- Urination difficulties
- Constipation
- Hematoma
- Infection
- Erosion
- Dyspareunia
- Disorders of micturition
- Overactive bladder, urgency
- De novo SUI
- Prolapse in other compartment
- Rectal and vesical perforation
- Perforation or damage to blood vessels, bladder, rectum, intraperitoneal organs or nerves

As with all implants, local irritation and/or a reaction to the foreign body may occur.

**SURGICAL PROCEDURE**

(Detailed surgical procedure available on request)

Surgery under general anaesthetic.

Abdominal incision or installation of trocars.

Exposure of the lumbosacral angle.

Dissection of the prerectal fascia after incision of the anterior wall of Douglas' cul-de-sac. The dissection continues laterally as far as the posterior elevator muscles, using a vaginal valve if necessary.

The using mesh is fixed as distally as possible on the pelvic floor.

Dissection of the anterior face of vagina.

The anterior mesh is fixed by a running stitch, avoiding transfixation. Depending on the strip used and the surgeon's preference, the meshes is passed laterally relative to the uterine isthmus. The uterus is held in place, or two branches are branches are passed on either side of it (for precut mesh or mesh cut longitudinally by the surgeon).

Fastening the anterior and posterior mesh to the lumbosacral angle without excessive tension (presacral ligament)

Peritonization must be as complete as possible to avoid contact between the prosthesis and the intraperitoneal constituents.

**PRECAUTIONS FOR USE**

Sacromesh Soft Prolaps' devices must **always** be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in gynaecological and urological surgery) who has a good knowledge of the device, its intended use, the instruments and the surgical technique.

Sacromesh Soft Prolaps' implants are delivered sterile (sterilised using ethylene oxide), inspect the packaging to be sure it is intact (do not use if the labels and/or protective packaging are damaged).

Do not use if the device is damaged or out of date.

Do not reuse and do not resterilize the prosthesis.

**STORAGE OF THE PROSTHESIS**

To be stored in a dry place away from light, at room temperature, in its original packaging.

**IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

[Table of contents](#)

**SACROMESH SOFT PROLAPS®  
IMPLANT UROGENITAL**

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