BIOMESH[®] CA.B.S'Air[®]

NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT

EXTRAPERITONEAL

en	Instructions for use	Page	2
fr	Notice d'instructions	Page	4
de	Gebrauchsanweisung	Seite	6
it	Istruzioni per l'uso	Pagina	8
es	Instrucciones de uso	Pagina	10
ελ	Οδηγίες χρήσης	Σελίδα	12

COUSIN BIOTECH

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NOT013 / 181010

Version du 10/10/2018

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

BIOMESH® CA.B.S.'Air®
NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT
Expandable parietal prosthesis with inflatable balloon for visceral surgery
STERILE PRODUCT FOR SINGLE USE
EXTRA-PERITONEAL

Description

The BIOMESH® CA.B.S.'Air® prosthesis consists of (according to the references on the box):

- A parietal prosthesis made of a plug and a polypropylene mesh and two polyethylene terephthalate and ePTFE fixation threads crimped of stainless steel needles.
- A balloon which allows the deployment of the parietal prosthesis

Implanted materials :

Polypropylene, Polyethylene terephthalate, Biocompatible dye, ePTFE.

Materials in contact with the patient during the surgery :

Stainless steel, Poly vinyl chloride.

No human nor animal origin. Non resorbable.

Indications

The BIOMESH® CA.B.S.'Air® prosthesis is intended for the use in all forms of hernia repair requiring reinforcement with non-absorbable support material.

Performances

The BIOMESH® CA.B.S.'Air® prosthesis is particularly designed to be **extra-peritoneally** implanted for the mini invasive surgical treatment (small incision) of all form of hernia.

Biocompatible non resorbable products, cause fibrosis reaction taking over the reinforcement after six months of implantation. They have the advantage of having a very high resistance to suture, have a great flexibility and allow an optimum and rapid integration and colonisation.

Contraindications

Do not use in the following cases :

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

Undesirable side effects

Like any implantable medical device, this implant is susceptible to involve possible undesirable effects:

- Recurrence
- Infection
- Discomfort

Precautions for use

The BIOMESH® CA.B.S.'Air® prosthesis are delivered sterile (ethylene oxide sterilization).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**THE BALLOON MUST BE IMPERATIVELY WITHDRAWN**

As specified on the product's labeling, the BIOMESH® CA.B.S.'Air® mesh is for single use only. It can not be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

Storage of the prosthesis



To be stored in a dry place away from light and at room temperature in its original packaging.

Important

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

SURGICAL PROCEDURE

BIOMESH® CA.B.S.'Air® is not intended to be used in the laparoscopic and endoscopic surgery. The prosthesis should always be installed with an anterior reinforcement for inguinal and crural hernia.

SURGICAL PROCEDURE					
<p>1-Installation:</p> <ul style="list-style-type: none"> - The prosthesis equipped with its balloon crosses completely the hernial defect, -its positioning becomes then entirely retro muscular 		<p>2-Deployment:</p> <ul style="list-style-type: none"> - The expansion and the spreading out of the prosthesis are obtained by inflating the balloon, (30ml), volume of a syringe. - Parietal fixing is carried out by two points thanks to the pre-positionned switched threads. 		<p>3- Prosthesis in place:</p> <ul style="list-style-type: none"> -The prosthesis is fixed, -This balloon must imperatively be withdrawn 	

Surgical technique for inguinal hernia repair

Horizontal inguinal incision of about 5 cm, beginning immediately above the pubic spine. This incision extends down to the level of the external abdominal oblique muscle. Placement of retractors. Short incision of the external abdominal oblique muscle at the top of the external orifice of the inguinal canal. Exploration of the inguinal canal enables 1) direct and indirect inguinal hernias to be differentiated, and 2) the dissection and removal of prehernia lipomas.

Direct inguinal hernia:

Simple exposure of the patulous transversalis fascia. Opening of the transversalis fascia. Retraction of the preperitoneal fatty tissue in order to create the space for a correct centred placement of the prosthesis. The prosthesis is introduced folded into the preperitoneal space that was created. Inflation of the balloon in order to expand the prosthesis while ensuring its proper centring. Fixation of the prosthesis with two prepositioned sutures that are passed through the posterior wall of the inguinal canal. Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned. The sutures are tied without tension. Tension-free closure of the transversalis fascia in front of the prosthesis by a continuous, absorbable suture.

External abdominal oblique hernia:

Opening of the spermatic cord. Detection of the peritoneal sac of the hernia. The parts of the spermatic cord are located and isolated. Deep dissection of the peritoneal sac in order to facilitate its reintegration into the abdominal cavity and to expose sufficient preperitoneal space for the correct placement of the prosthesis. After retraction of the peritoneal sac of the hernia, the prosthesis is placed folded in the preperitoneal space. Inflation of the balloon enabling expansion of the prosthesis while ensuring that it remains well-centred relevant to the deep orifice of the inguinal canal. Fixation of the prosthesis with two prepositioned sutures that are passed through the posterior wall of the inguinal canal without tension. Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned.

The wall must be reinforced by placing a second prosthesis around the spermatic cord using the Lichtenstein technique. Closure of the incision of the external abdominal oblique muscle with absorbable suture material. Subcutaneous approximation and cutaneous closure with a continuous intradermal, absorbable suture.

Surgical technique for umbilical hernia repair

Hemicircumferential periumbilical incision centred on the hernia. Dissection of the peritoneal sac, enabling its abdominal reintegration by clearing sufficient space behind the aponeurosis for the prosthesis to be implanted. The folded prosthesis is slipped through the aponeurotic orifice of the hernia in order to position it in the preperitoneal space.. Inflation of the expansion balloon of the prosthesis while ensuring that the prosthesis remains centred relative to the hernia orifice. Fixation of the prosthesis using prepositioned sutures passing through the aponeurosis. Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned. If possible, aponeurotic approximation in front of the prosthesis without tension. Subcutaneous and cutaneous closure with absorbable suture material.

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