Adhesix[®] BIORING[®]

ADJUSTABLE GASTRIC BANDING & REPLACEMENT KIT POM PORT

en	Instructions for use	Page	2
fr	Notice d'instructions	Page	4
de	Gebrauchsanweisung	Seite	6
it	Istruzioni per l'uso	Pagina	8
es	Instrucciones de uso	Pagina	10
ελ	Οδηγίες χρήσης	σελίδα	12
hu	Használati utasítás	Oldal	20

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NOT178/180911

Version du 11/09/2018

**ADJUSTABLE GASTRIC BANDING ADHESIX® BIORING®
REPLACEMENT KIT POM PORT
STERILE PRODUCT FOR ONE USE**

DESCRIPTION

-Kit including a rinse catheter, a catheter, a connector, a plug, a huber needle, a port and a connection system.
-The various parts of the kit were designed for fitting perfectly to the **ADHESIX® BIORING®** adjustable gastric banding.

MATERIALS

The injection port is composed of polyoxymethylene impregnated with barium sulfate. Its septum is composed of dimethylsiloxane. A textile part is fixed to the port by means of a thread in polyethylene terephthalate. This knitting is in polypropylene soaked with an adhesive resorbable hydrogel (Polyvinylpyrrolidone – Polyethylene glycol). The septum is made of silicone. The catheter is made of silicone (dimethylsiloxane) and barium sulfate labelled with a biocompatible ink. The connection system is made of silicone (dimethylsiloxane) and barium sulfate. The plug is made of stainless steel and silicone (dimethylsiloxan). The connector is made of stainless steel.
No human or animal origin. Semi-resorbable

INDICATIONS

Kit to be used when a **ADHESIX® BIORING®** gastric banding has been implanted and in the event of preoperative and postoperative complications, repeated surgery on patients following complications (infection of the port, leakage, catheter damaged).

PERFORMANCES

- The port is an X-Ray- opaque injection site and it allows per-cutaneous adjustment of the diameter of the gastric banding
- The connection system secures the connection between the port and the catheter by an on-clipping system

The injection port is self-adhesive. It is not necessary to use sutures; simply position it on the muscle:

- inserted in the aponeurosis of the rectus abdominis muscle under the left costal edge,
- inserted subcutaneously under the left costal edge,
- inserted pre-sternally, on the xiphoid apophysis;
- in any other position of the surgeon's appreciation.

Plug: In the event of an infection at the junction between the port and the catheter, the port is replaced with the plug until a new port is fitted. The plug is connected to the free end of the catheter of the **ADHESIX® BIORING®** adjustable gastric banding.

Connector: When the catheter on the **ADHESIX® BIORING®** adjustable gastric banding is damaged (cut, leakage), the connector and the catheter from the kit can be used to replace the damaged system.

Always refer to the instructions for the **ADHESIX® BIORING®** adjustable gastric banding marketed by COUSIN BIOTECH.

The product must be implanted **only** by a qualified surgeon and well formed for the use of it.

CONTRA-INDICATIONS

Do not use the implant in the following circumstances:

- Allergy to any of the components
- Growing children
- Do not implant to fit with a gastric band that is not a **BIORING®** adjustable gastric banding
- Inflammatory disease or congenital or acquired anomalies of the gastrointestinal tract.
- Severe organic disease.
- Potential bleeding disorder of the upper gastrointestinal tract or patients who regularly or systematically take aspirin or non-steroidal anti-inflammatory drugs.
- Diagnosis or preexisting symptoms or family history of connective tissue autoimmune disease.
- Pregnancy.
- Patients with a drug or alcohol addiction.
- Patients who refuse to accept the dietary restrictions required by the procedure.
- Serious psychiatric disorders.

POSSIBLE COMPLICATIONS

Complications that may arise from the use of an adjustable gastric band include risks associated with the anaesthesia, the surgical technique, and the degree of intolerance of the patient to any foreign body implanted in the body.

Perioperative complications (due to the surgical technique):

- Lesions of the spleen Perforation of the stomach or oesophagus - other bleeding

In the immediate postoperative period:

- Nausea, vomiting (especially in the first few days after the intervention, and when the patient eats more than recommended)
- Bloating due to gas
- Dysphagia, dehydration
- Thromboembolic complications
- Infectious complications
- Postoperative bleeding

Delayed postoperative complications:

- Dilatation of the upper pouch,
- Slippage of the band
- Gastroesophageal reflux,
- Pyrosis
- Gastritis of the upper pouch,
- Migration of the band,
- Displacement of the injection port,
- Constrictive obstruction (due for example to oedema, food, an inadequate calibration, slippage of the band, or to torsion or dilatation of the upper pouch),
- Deflation of the band (due to a leak in the band, injection port, or connection tube),
- Rapid weight loss may lead to symptoms of malnutrition, anemia, and associated complications,
- Infection of the injection port,
- Constipation

PRECAUTION FOR USE AND FOR STORAGE

- The components of the implantable port and of the connecting system have to be used only with the **BIORING®** adjustable gastric banding
 - The band can be inflated immediately using physiological saline (NaCl 0.9%) The volume, as well as the interval before filling, is at the discretion of the practitioner/surgeon (maximum volume, see the next table).
- Inflating of the band: maximum volume of physiological saline according to sizes of the bands

Band sizes (diameter : mm)	Maximum volume of physiological saline (cm ³)
Ø 21/24	10
Ø 27/29	10
Ø 36/38	15

- To keep in a dry place without any light, and in an ambient atmosphere, in its packaging of origin
- Please check the perfect integrity of the packaging (do not use in case of damaged packaging protection and/or damage of the label)
- Do not use in case of damaged product or if the product has ran out date

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