CORTICAL FIXATION SYSTEM

ANCHOR SYSTEM

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	3
de	<u>Gebrauchsanweisung</u>	Seite	4
it	<u>Istruzioni per l'uso</u>	Pagina	5
es	<u>Instrucciones de uso</u>	Pagina	6
pt	<u>Nota de instruções</u>	Pagina	7
ελ	<u>Οδηγίες χρήσης</u>	σελίδα	8
bg	<u>Инструкции за употреба</u>	Страница	9
tu	<u>Kullanım Kilavuzu</u>	Sayfa	10

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Date de marquage CE: Nov 2005

NOT189/181217
Version du 17/12/2018

Instructions for use**CORTICAL FIXATION SYSTEM****ANCHOR SYSTEM****STERILE PRODUCT FOR SINGLE USE****Description**

Ligaments and attachment unit of the tendons.

Implanted materials

Titanium, Polyethylene terephthalate.

Materials in contact with the patient during the surgery

Polyethylene terephthalate, DC Green 6

Not of human or animal origin.

Indications

Tendinous and ligamentary anchoring device designed for orthopaedic reconstruction, especially the anterior cruciate ligament of the knee.

Contraindications

Insufficient or poor-quality bone tissue.

Sensitivity to foreign bodies.

Do not implant in the case of allergy to materials.

The device should not usually be implanted in growing children.

Possible adverse side effects

Possible complications linked to the implant:

1- Breaking of the loop, breaking of the traction wires, breaking of the button

2- Failure when placed under tension

3- Migration or upturn of the button

4- Reaction to a foreign body

Possible complications linked to insertion of the implant:

1- Pain in the incision area or surgical site

2- Local and transient inflammatory reactions

3- Infection

4- Wound dehiscence.

Precautions for use

Inspect the packaging to be sure it is intact; Do not use if the labels and/or protective packaging are damaged.

Do not use in case of damaged product.

DO NOT REUSE**DO NOT RESTERILIZE**

As specified on the product's labeling, the implant is for single use only. It can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency).

COUSIN BIOTECH does not make any guarantees or recommendations as regards the use of a particular brand of fixation devices. The properties of these fastenings can be the subject to modifications from the manufacturer, which cannot be controlled by COUSIN BIOTECH.

This device should be implanted **only** by a qualified surgeon who has been trained in the use of this product.

Never use staples or sharp-edged screws.

Surgical technique: semitendinosus technique in knee ligamentoplasty.

This system promotes normal healing but it is designed neither to substitute for normal bone tissue nor to bear the weight of the body if bone consolidation is incomplete.

The patient should be informed of the limitations inherent to the implant's nature and function, especially with respect to stress and the strain imposed on the system by the body before complete consolidation of the bone tissue.

Precautions for storage

To be stored in a dry place away from light, at room temperature, in its original packaging.

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