QUICK ANCHOR BIOLIG®

**NON RESORBABLE
REINFORCING JOINT LIGAMENT WITH AN ANCHORING SYSTEM**

en	Instructions for use	Page	2
fr	Notice d'instructions	Page	3
de	Gebrauchsanweisung	Seite	4
it	Istruzioni per l'uso	Pagina	5
es	Instrucciones per l'uso	Pagina	6
pt	Nota de instruções	Pagina	7

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NOT077/190412
Version du 12/04/2019

Instructions for use

QUICK ANCHOR *BIOLIG*[®] REINFORCING JOINT LIGAMENT WITH AN ANCHORING SYSTEM STERILE PRODUCT FOR ONE USE

DESCRIPTION

The purpose of ligaments is to surgically reinforce instabilities in the various articulations. They may be fixed in an intra or extra-articular basis. Ligament reinforcement should always be used to complement suture or grafting of the natural ligament..

MATERIALS

Anchors: TA6V ELI ISO 5832-3

Ligament: Polyethylene terephthalate

Needles: Stainless steel (5832-1) (contact<24 hours – Limited contact)

Handle: Polyoxymethylene and Stainless steel (5832-1) (not implanted)

Nor human neither animal origin

INDICATIONS

Reinforcing ligaments for orthopedic surgery

POSSIBLE ADVERSE SIDE EFFECTS

- Routine surgical complications
 - superficial or deep infection
 - anesthesia complications
 - local hematoma, blood loss
- Anchor breakage, anchoring failure, suture breakage
- Permanent elongation of the ligament / Rupture of the ligament
- Postoperative loss of motion, frozen joint
- Inflammatory / allergic reactions / arthro fibrosis (around 1%) / Bone resorption (osteolysis)
- Ongoing instability episodes
- Neurovascular injury, peripheral neuropraxia
- Wound dehiscence, serous leak
- Post-surgical pain, persistent pain
- Ligament injury

WARNINGS

Inspect the packaging to be sure it is intact; Do not use if the labels and/or protective packaging are damaged
Do not use if the device is damaged.

DO NOT REUSE

DO NOT RESTERILIZE

Avoid using the device in a septic environment.

Do not implant in an infected site or in the case of allergy to the materials.

COUSIN BIOTECH advises the use of the provided ancillary. This one is of single use and is not resterilization.

COUSIN BIOTECH does not make any guarantees or recommendations as regards the use of a particular brand of fixation devices. The properties of these fastenings can be the subject to modifications from the manufacturer, which cannot be controlled by COUSIN BIOTECH. It is necessary to avoid fasten or the screws on sharp board.

It is not advised to implant the device in a child while he is still growing.

This device should be implanted **exclusively** by a qualified surgeon who has been trained in the use of this product.

Attention with the sterile handle of one use.

STORAGE

To be stored in a dry place away from light, at room temperature, in its original packaging

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