# BIOMESH<sup>®</sup> CA.B.S'Air<sup>®</sup> Light

## NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT EXTRA-PERITONEAL USE

---

en	<a href="#">Instructions for use</a>	Page	2
fr	<a href="#">Notice d'instructions</a>	Page	4
de	<a href="#">Gebrauchsanweisung</a>	Seite	6
it	<a href="#">Istruzioni per l'uso</a>	Pagina	8
es	<a href="#">Instrucciones de uso</a>	Pagina	10

### COUSIN BIOTECH

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**NOT232-190517**  
Version du 17/05/2019

**Caution : Federal law (USA) restricts this device to sale,  
distribution and use by or on the order of a physician.**

## Instructions for use

# **BIOMESH® CA.B.S.'Air® Light**

**NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT**  
**Expandable parietal prosthesis with inflatable balloon for visceral surgery**  
**STERILE PRODUCT FOR SINGLE USE**  
**EXTRA-PERITONEAL USE**

### **DESCRIPTION**

The **BIOMESH® CA.B.S.'Air® Light** prosthesis consists of:

- A parietal prosthesis made of a plug in polypropylene and two polyethylene terephthalate and ePTFE fixation threads crimped with stainless steel needles.
- A balloon which allows the deployment of the parietal prosthesis (by using the syringe included in the medical device).

### **IMPLANTED MATERIALS :**

Polypropylene, Polyethylene terephthalate, Biocompatible dye (DC Green #6), ePTFE,

### **MATERIALS IN CONTACT WITH THE PATIENT DURING THE SURGERY :**

Stainless steel, Poly vinyl chloride.

Nor human neither animal origin. Non resorbable.

### **INDICATIONS**

The **BIOMESH® CA.B.S.'Air® Light** is intended for the use in all forms of hernia repair requiring reinforcement with non-absorbable support material

### **PERFORMANCES**

The **BIOMESH® CA.B.S.'Air® Light** medical device is particularly designed to be extra-peritoneally implanted for the mini invasive surgical treatment (small incision) of all form of hernia.

Biocompatible non resorbable product. They have the advantage of having a very high resistance to suture, have a great flexibility and allow an optimum and rapid integration and colonisation.

### **CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

### **UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to involve possible undesirable effects:

- |                     |                           |                  |
|---------------------|---------------------------|------------------|
| - Discomfort/Pain   | - Recurrence              | - Adhesions      |
| - Erosions          | - irritation nearby organ | - Infection      |
| - Inflammation      | - Fistula formation       | - Seroma         |
| - mesh deformation  | - Hematoma                | - Mesh migration |
| - Allergic reaction | - Foreign body reaction   |                  |

### **PRECAUTIONS FOR USE**

**BIOMESH® CA.B.S.'Air® Light** devices are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and trained on the use of the product (knowledge of anatomy and visceral surgery).

**IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**  
**THE BALLOON MUST BE IMPERATIVELY WITHDRAWN**



As indicated on the product labelling, the **BIOMESH® CA.B.S.'Air® Light** implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

### **STORAGE OF THE PROSTHESIS**

To be stored in a dry place away from light and at room temperature in its original packaging.

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

## **SURGICAL PROCEDURE**

<b>SURGICAL PROCEDURE</b>					
<b>1-Installation:</b> - The prosthesis equipped with its balloon crosses completely the hernial defect, -its positioning becomes then entirely retro muscular		<b>2-Deployment:</b> - The expansion and the spreading out of the prosthesis are obtained by inflating the balloon, (30ml), volume of a syringe. - Parietal fixing is carried out by two points thanks to the pre-positionned sutures.		<b>3- Prosthesis in place:</b> -The prosthesis is fixed, -This balloon must imperatively be withdrawn	

### **SURGICAL TECHNIQUE FOR UMBILICAL HERNIA REPAIR**

Hemicircumferential periumbilical incision centred on the hernia. Dissection of the peritoneal sac, enabling its abdominal reintegration by clearing sufficient space behind the aponeurosis for the prosthesis to be implanted. The folded prosthesis is slipped through the aponeurotic orifice of the hernia in order to position it in the preperitoneal space.. Inflation of the expansion balloon of the prosthesis while ensuring that the prosthesis remains centred relative to the hernia orifice. Fixation of the prosthesis using prepositioned sutures passing through the aponeurosis. Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned. If possible, aponeurotic approximation in front of the prosthesis without tension. Subcutaneous and cutaneous closure with absorbable suture material.

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