IMPLANT BIOMESH®

NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT EXTRAPERITONEAL

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>		3
de	<u>Gebrauchsanweisung</u>	Seite	4
it	<u>Istruzioni per l'uso</u>	Pagina	5
es	<u>Instrucciones de uso</u>	Pagina	6
pt	<u>Nota de instruções</u>	Pagina	7
ελ	<u>Οδηγίες χρήσης</u>	σελίδα	8
bg	<u>Инструкция за употреба</u>	Страница	9
da	<u>Brugsvejledning</u>	Side	10
et	<u>Kasutusjuhend</u>	lehel	11
fi	<u>Käyttöohjeet</u>	Sivu	12
hu	<u>Felhasználási Útmutató</u>	Oldal	13
lv	<u>Lietošanas instrukcija</u>	Lappuse	14
nl	<u>Gebruiksaanwijzing</u>	Pagina	15
ro	<u>Instrucțiuni de Utilizare</u>	Pagină	16
tu	<u>Kullanım Kilavuzu</u>	Sayfa	17

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NOT029/190625

Version du 25/06/2019

This release is the last update of the instructions of use and replace the previous edition

IMPLANT BIOMESH®
NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT
SINGLE USE STERILE PRODUCT

DESCRIPTION

BIOMESH® meshes are non resorbable parietal reinforcement implants.

IMPLANTED MATERIALS

Mesh : polypropylene (knitted monofilament)

Marked thread : Polyethylene terephthalate tinted with a biocompatible dye (See on the package).

Origin neither human nor animal – Non resorbable.

INDICATIONS

Repair and parietal reinforcement for inguinals and crurals hernias and eventrations

PERFORMANCE

Meshes are particularly adapted either to coelioscopic surgery, or to conventional surgery.

The BIOMESH® is designed for extraperitoneal implantation only.

Meshes, biocompatible and non resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a shape memory, a very easy dimensioning, a very high resistance to suture, a great flexibility and allow an optimum integration and colonisation.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects:

- Discomfort
- Infection
- Recurrence

USAGE PRECAUTIONS

BIOMESH® meshes are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labeling, the BIOMESH® mesh is for single use only. It can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

STORAGE

To be stored in a dry place away from light and at room temperature in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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