LIFT[®]Kit

STERILE UROGENITAL IMPLANTS: SUBURETHRAL SUPPORT TAPE

en	Instructions for use Surgical Procedure	Page	2 26
fr	Notice d'instructions Technique Opératoire	Page	4 26
de	Gebrauchsanweisung Operationstechnik	Seite	6 26
it	Istruzioni per l'uso Tecnica Operatoria	Pagina	8 26
es	Instrucciones de uso Técnica Operatoria	Pagina	10 26
pt	Nota de instruções Técnica Operatoria	Pagina	12 26
ελ	Οδηγίες χρήσης Χειρουργική Διαδικασία	σελίδα	14 26
hu	Használati útmutató Sebészeti eljárás	oldal	16 26
nl	Gebruiksaanwijzing Chirurgische procedure	Page	18 26
pl	Instrukcja obsługi Technika operacyjna	Strona	20 26
sk	Návod k použítí Chirurgický postup	Strana	22 26
cs	Pokyny k použítí Chirurgický postup	Strana	24 26

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NOT219/190618

Version du 18/06/2019

This release is the last update of the instructions for use and replaces the previous edition

LIFT®Kit

Suburethral support tape

DESCRIPTION

The LIFT® kit comprises a tape and a winged guide, 3 ancillary devices (1 pair for the transobturator technique and 1 for the retropubic and suprapubic approach).

They are single use and delivered sterile.

The LIFT® prosthesis consists of a tape of knitted monofilament polypropylene, packaged in a plastic sheath.

It is prolonged by self-tightening loops for placing.

The ancillary devices have different shapes depending on the chosen approach: 1 curved needle for the retropubic and suprapubic approach, a pair of helical needles for the in-out and out-in transobturator approaches.

The end of the ancillary needles has a slot for fixing self-tightening loops.

A single-use, sterile guide is supplied with the implant to secure the passage of the ancillary device during the in-out transobturator approach.

When placing the self-tightening loop in the slot, always make sure it is as tight as possible and that the thread is in the slot designed for this purpose so that the prosthesis does not come apart from the ancillary when passing in the patient's tissues. (See figure 1 below)


Fig 1

MATERIALS

- **Prosthesis:** polypropylene (knitted monofilament)
 - **Thread:** Polyethylene terephthalate + biocompatible dye + dimethyl siloxane
 - **Ancillary device:** Stainless steel needle, polyoxymethylene handle.
 - **Guide:** Liquid crystal polymer + biocompatible dye.
- Not human nor animal origin – Non absorbable.**

INDICATIONS

The LIFT® prosthesis is designed for use as a suburethral support **tape** for the surgical treatment by vaginal approach in women.

PERFORMANCE

These biocompatible, not biodegradable prostheses cause reactive fibrosis that produces additional urethral support.

These prostheses are easy to insert due to their specially designed shape, excellent flexibility and optimum, rapid colonization and integration.

CONTRAINDICATIONS

Do not implant in the following cases:

- Allergy to any one of the components
- For women receiving anticoagulant therapy.
- Urinary tract infection
- Pregnant women
- Skeletally-immature children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

SIDE EFFECTS

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, urethra, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection

Over-correction may lead to temporary or permanent obstruction of the lower urinary tract and to retention of urine.

The known risks of surgery for the treatment of incontinence include hematoma, acute or chronic pain, scarification of the vagina, retraction or release of the mesh, infection, erosion, de novo urgency or voiding difficulties, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

USAGE PRECAUTIONS

The LIFT®Kit device should be implanted by qualified surgeons only (knowledge of gynaecological and urological anatomy and surgery) who are extremely familiar with the device.

The LIFT® Kit implant, guide and ancillary devices are delivered sterile (ethylene oxide sterilization). Check that the packaging is intact before using (do not use if the labels and/or protective packaging are damaged). Do not use if the device is damaged.

The plastic sheath that covers the sling must be removed and must not be implanted.

Do not reuse or resterilize the implant, the guide and ancillary devices.

IMPORTANT: DO NOT REUSE – DO NOT RESTERILIZE

As indicated on the product labelling, the implant, guide and ancillary devices are designed for single use. They should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

STORAGE

To be stored in a dry place away from light, at room temperature, in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be done following recommendations of ISO 12891-1:2015 « Implants for surgery – Retrieval and analysis of surgical implants » Part 1: « Retrieval and Handling ».

Any explanted device must be send back, for analysis, following the current protocol. This protocol is available on demand to COUSIN BIOTECH. It is important to note that any implant that should not have been cleaned and disinfected before expedition must be contained in a sealed package.

The elimination of explanted medical device must be conducted in accordance with standards in the country for the disposal of infectious hazards waste.

The elimination of a non-implanted device is not the subject of specific recommendations.

INFORMATION REQUEST AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible.

In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time .

For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor, or COUSIN BIOTECH at contact@cousin-biotech.com.

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
Transobturator Methode

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Via transobturatriz

Installation of the patient, anaesthesia and infiltration performed as specified on the page named «Preliminary Information 1-2-3»

1. Trans-obturator incisions of approximately 1 cm each, positioned at least 1 cm above the horizontal line passing by the urethral meatus at 1 cm from the ischiopubic branch.
 2. Vaginal incision and urethro-vaginal dissection.
 3. Passage of heliocoid ancillary by way of the trans-obturator (blue handle). (the tunnelling device is inserted into one of the trans-obturator incisions).
 4. • Hooking of LIFT® bandage onto needle.
 5. • Passage of the tape; and heterolateral passage: performance of the 3 previous operations on the heterolateral side.
 6. • The tape is deployed by pulling on the thread. There should remain a space of about one centimetre between the tape and the urethra.
 7. • Cutting of the tape, flush with the skin at the level of the incisions. Hemostasis is checked after tensioning of the tape.
- Closure of the vaginal wall by a crossed suture with absorbable thread and closure of the cutaneous incision.

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