INTRASPINE®

DEVICE FOR THE INTER SPINOUS SPACE WITH LAMINAR SUPPORT

| | | | |
|----|---------------------------------------|--------|----|
| en | Instructions for use | Page | 2 |
| fr | Notice d'instructions | Page | 6 |
| de | Gebrauchsanweisung | Seite | 10 |
| it | Istruzioni per l'uso | Pagina | 14 |
| es | Instrucciones de uso | Pagina | 18 |
| pt | Nota de instruções | Página | 22 |
| ελ | Οδηγίες χρήσης | σελίδα | 26 |

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This release is the last update of the instructions of use and replace the previous edition

Instructions for use

IntraSPINE®

DEVICE FOR THE INTER-SPINOUS SPACE WITH LAMINAR SUPPORT

SINGLE USE STERILE PRODUCT

DESCRIPTION

The IntraSPINE® prosthesis is an interlaminar medical device that comprises a dimethyl siloxane (silicon) wedge (available in several sizes) adapted to the intervertebral space and two ligaments in polyethylene terephthalate.

The wedge is covered with woven polyethylene terephthalate with an area covered by sheeting on the front face in the medullary zone. Tension is applied to the ligament by a titanium loop (TA6V).

MATERIALS

IntraSPINE® prosthesis: Dimethyl siloxane, Polyethylene terephthalate, Titanium.

IntraSPINE® ancillaries and trial prostheses: Stainless steel, Acetal copolymer, Dimethyl siloxane .

Non-human and non-animal origin – Non-resorbable.

INDICATIONS

The IntraSPINE® prosthesis is recommended in lumbar spine and sacral area (from L1 to S1 levels) in the following cases:

- Arthropathic facet-syndrome
- Foraminal stenosis
- Degenerative discopathy
- Intervertebral ligament insufficiency

PERFORMANCE

The IntraSPINE® prosthesis is an elastic support device with a cushioning effect. Its interlaminar implantation enables it to be close to the centre of rotation, thereby producing the cushioning effect and providing tension on the interspinous ligamentary structures due to the elasticity of the supra-spinous wedge.

The principal function of the ligaments surrounding the upper and lower vertebrae in relation to the area receiving the prosthesis is to assist the posterior ligament during flexion. The second ligament stabilizes the interlaminar device on its site.

IntraSPINE® has been validated for a long term implantation in the human body.

CONTRAINDICATIONS

1) Do not use the prosthesis in the following cases:

- Allergy to one of the components
- Growing child
- Infected site
- Pregnant woman

2) Factors likely to compromise successful implantation:

- Severe osteoporosis
- Major deformations of the spinal column
- Local bone tumors
- Systemic or metabolic disorders
- Infectious diseases
- Obesity
- Drug addictions
- Intense physical activity, for example competitive sports or hard labor
- Surgery at more than 2 levels with IntraSPINE Prosthesis

UNDESIRABLE SECONDARY EFFECTS

1) Complications likely related to the device:

- Infection of the prosthesis, rupture of the ligament, loosening, fixed wing breakage, permanent elongation of the ligament, displacement of the wedge, spondylolisthesis to the dynamically stabilized segment, persistent pain attributed to the dynamically stabilized level and required revision surgery with interbody fusion at that level;

Spinous process fractures: intraoperative, post-operative or during device removal (removal due to infection) and can be associated with dislocation of the implant;

Late reoperation due to absence of recovery;

In certain cases the medical instrument had to be removed and pedicle screw fusion had to be performed.

2) Complications unlikely related to the device but related to surgery or to the patient's conditions:

- Infection
- Pseudomeningocele, fistula, persistent CSF leakage, meningitis
- Loss of neurological function,
- Cauda equina syndrome, neuropathy, neurological deficits, arachnoiditis, and/or muscle loss, compression around nerves and/or pain
- Urinary retention or loss of bladder control or other types of urological system compromise
- Fracture, Bone loss or decrease in bone density
- Herniated nucleus pulposus, disc disruption or degeneration at, above, or below the level of surgery: Adjacent segment degeneration
- Non-union or pseudarthrosis, delayed union. Mal-union,
- Cessation of any potential growth of the operated portion of the spine
- Loss of or increase in spinal mobility or function
- Hemorrhage, hematoma, edema, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, wound infection, damage to blood vessels, or other types of cardiovascular system compromise
- Development of respiratory problems,

MRI SAFETY:

The implants are composed of non ferromagnetic materials and present a geometry non susceptible to generate induced current.

Moreover, as they are fixed to tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan.

Their safety, in particular in terms of heating and migration of implant has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded to a conditional compatibility for MRI scan until 1.5 Tesla.

As a precautionary measure, it is recommended to avoid MRI scans within the 48h of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential.

It has to be noted that the devices which present a high contrast with the biological environment can generate « artifacts » that has to be taken into account for the perfect execution and interpretation of imaging exams.

For this purpose, it has to be recommended to the patient who has this implant to warn as far as possible the concerned health professionals (radiologists and radiology operators) about the presence of this implant before these exams.

POSTOPERATIVE PRECAUTIONS

The patient must be informed of the postoperative cares and of precautions.

The use of a brace could be necessary but is not mandatory. This decision is the surgeon one.

PRECAUTIONS FOR USE

Before usage, verify the integrity of the prosthesis, trial prostheses, ancillaries and packaging.

Do not use in the event of deterioration of the prosthesis and/or the labels, and/or the trial prostheses and/or the ancillaries and/or the packaging.

Do not use if the prosthesis has exceeded its use-by date.

Unpacking and manipulation of the implant must be carried out following asepsis standards in order to guarantee the sterility.

When opening the carton box, the surgeon has to make sure that the lot number, the model and the size of the prosthesis are equal to the data on the label of the inner pouch. Otherwise the use of the prosthesis is strictly prohibited.

Prevent the prosthesis, the trial prostheses and ancillaries from coming in to contact with objects that may alter their surface.

Operation conditions: Temperature from +32°C up to 42°C. Implants should be resistant to biologic fluids and the discharge of body tissues, which they contact during use. Therefore, the IntraSPINE® prosthesis has been validated for a long term implantation and, currently, Cousin Biotech has 10 years of feedback concerning these products.

The IntraSPINE® prosthesis, trial prostheses and ancillaries must only be used by a qualified surgeon who is trained in spinal surgery and to the surgical technique of IntraSPINE®. The information contained in this notice is necessary but not sufficient to control the surgical technique.

The IntraSPINE® is made with soft materials, any contact with harder one may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials.

The prescription of the device is decided on by the surgeon, the only person qualified to do so.

Topping off is not contra-indicated but special attention is required. IntraSPINE® is made with soft materials, any contact with harder one may alter its mechanical behaviour and lifetime. It is mandatory to avoid contact with other materials

Use of the system of COUSIN BIOTECH ancillaries and specific trial prostheses is required for implantation of the IntraSPINE® prosthesis.

Implants must be handled and stored with great care, in a dry place, protected from light and at room temperature.

Transportation: Transportation of the devices should be done in the closed vehicles in accordance with the rules for those transport vehicles. There are no special transportation requirements for the device.

IMPORTANT : DO NOT REUSE AND DO NOT RESTERILISE THE INTRASPINE® PROSTHESIS

As indicated on the product labelling, the prosthesis is designed for single use. It should in no event by reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be done following recommendations of ISO 12891-1:2015 « Implants for surgery – Retrieval and analysis of surgical implants » Part 1: « Retrieval and Handling ».

Any explanted device must be send back, for analysis, following the current protocol. This protocol is available on demand to COUSIN BIOTECH. It is important to note that any implant that should not have been cleaned and disinfected before expedition must be contained in a sealed package.

The elimination of explanted medical device must be conducted in accordance with standards in the country for the disposal of infectious hazards waste.


The elimination of a non-implanted device is not the subject of specific recommendations.

Warning about ancillaries and trial prostheses

The IntraSPINE® ancillaries are Class I medical devices, intended for temporary use and are re-usable.

The IntraSPINE® trial prostheses are Class IIa medical devices, intended for temporary use and are re-usable.

Unlike the IntraSPINE® prosthesis, ancillaries and trial prostheses (included with the instrument kit/set) are furnished unsterile.

KITBIOINTRA / RCBINBIOIU

Before first use

Ancillaries and trial prostheses must be:

- Cleaned in the appropriate way
- Sterilised in an autoclave

DECONTAMINATION AND CLEANING OF ANCILLARIES

Decontamination and cleaning are performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

This step has to take place before **the first utilization** and **immediately after the utilization** to avoid adhesion of particles or dry secretion to the instruments

The detachable instruments have to be dismantled

Products advised for the cleaning

- an enzymatic detergent neutral
- Or products adapted for the cleaning (neutral or soft alkaline) of chirurgical instruments in compliance with the regulations in effect.

- **Warning:** Don't use corrosive factor or caustic cleaning product

2 process are possible

A- Chart Automated cleaning in a disinfectant washer with manual Pre-cleaning:

| | Step | Step Description | Step Instruction | Accessories | Duration |
|----------------|------|-----------------------|---|--|---|
| Cleaning Steps | 1 | Contamination Removal | Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush | <ul style="list-style-type: none"> - Tap water - Soft bristle brush (Do not utilize metal cleaning brushes) | Until all visible soils removed |
| | 2 | Pre-soak | Prepare a detergent solution at ambient temperature (15-25°C / 59-77°F). Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position. | <ul style="list-style-type: none"> - Detergent - Water jet pistol | Until product is visually clean |
| | 3 | Ultrasonic cleaning | Clean the instrument in the ultrasonic tray including the solution. Be careful every unclench connector or screws and bolts by the vibrations will be reassure. | <ul style="list-style-type: none"> - Ultrasonic tray | Minimum duration 10 minutes, this duration can vary according product |
| | 4 | Soak | Soak the instrument in tap water during 30 seconds. After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position. | <ul style="list-style-type: none"> - Reverse osmosis water | Minimum 30 seconds of soak |
| | 5 | Automated Washer | Place entire disassembled device into the automated washer | <ul style="list-style-type: none"> - Automated Washer - Disassembled Device | Minimum total cycle time: 39 minutes |

| Instructions of cleaning in washer / disinfectant | | | |
|---|------------|---------------------|-------------------------|
| Cycles parameters | Time | Minimum temperature | Type of detergent/water |
| Pre-cleaning | 2 minutes | <45°C / <113°F | Tap water |
| Cleaning | 5 minutes | 55°C / 131°F | Detergent |
| Draining | 2 minutes | <45°C / <113°F | Tap water |
| Thermic draining | 5 minutes | 90°C / 194°F | Reverse osmosis water |
| Drying | 25 minutes | > 70°C / > 158°F | Not applicable |

| | | | | | |
|-----------------|---|--------------|--|--|-------------------------------|
| Finishing Steps | 6 | Final Rinse | Rinse carefully with distilled water | <ul style="list-style-type: none"> - Distilled water | Minimum 1 minute |
| | 7 | Final Drying | Dry device utilizing medical quality filtered air | <ul style="list-style-type: none"> - Medical quality filtered air | Until product is visually dry |
| | 8 | Tidying | The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization. | | |

B- Chart Just manual cleaning :

| | Step | Step Description | Step Instruction | Accessories | Duration |
|--------------------|------|--------------------------|--|--|--|
| Cleaning Steps | 1 | Contamination Removal | Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush | <ul style="list-style-type: none"> - Tap water - Soft bristle brush (Do not utilize metal cleaning brushes) | Until all visible soils removed |
| | 2 | Drying | Dry the device utilizing a dry non-shedding wipe. Medical quality filtered air may be utilized if available | <ul style="list-style-type: none"> - Non-scheduling wipe - Medical quality filtered compressed air | Until product is visually dry |
| Disinfection Steps | 3 | Disinfection Application | Prepare a low-foaming neutral enzymatic detergent solution, using tap water (15-25°C / 59-77°F). Soak the instrument in the open position (if possible) during minimum 1 minute. During the immersion, actuate mobile parts of the instrument minimum 3 times (if possible), so that the detergent can access all parts of the instrument. | <ul style="list-style-type: none"> - Detergent | Minimum duration 1 minute, this duration can vary according product |
| | 4 | Drying | After the soaking time, take out the instrument and wipe it with a disposable cloth. Then, place the instrument in a new bath of enzymatic detergent solution using warm tap water. | <ul style="list-style-type: none"> - Disposable cloth - Detergent | Until product is visually dry |
| | 5 | Manual Disinfection | Brush all surface of the instrument using a soft brush for 2 minutes. Adjust brushing time if needed. Actuate mobile parts of the instrument 3 times (when applicable) and use a syringe to access all difficult areas. Use a volume of 60 mL of the detergent solution. | <ul style="list-style-type: none"> - Soft bristle brush - Detergent - Syringe | Manual cleaning time duration is complete when the device's surface, joints, and crevices have been manually cleaned |
| | 6 | Rinse and finition | Soak the instrument in reverse osmosis water for 1 minute. Use a syringe and 60mL of reverse osmosis water to access all difficult areas. Repeat soaking 2 additional times using fresh water | <ul style="list-style-type: none"> - Reverse osmosis water - Syringe | Minimum 1 minute |
| | 7 | Final Rinse | Rinse carefully with distilled water | <ul style="list-style-type: none"> - Distilled water | Minimum 1 minute |
| | 8 | Final Drying | Dry device utilizing medical quality filtered air | <ul style="list-style-type: none"> - Medical quality | Until product is |

| | | | | |
|---|---------|--|--------------|--------------|
| | | | filtered air | visually dry |
| 9 | Tidying | The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization. | | |

Note: In case of **suspection** of **Creutzfeld-Jakob Disease (CJD)**, incinerate the ancillaries and the fixation pliers
 In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

STERILIZATION OF ANCILLAIRES

Sterilization is performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

1- The sterilization in autoclave is to make in a **specific container** (fenced basket, Tray, Plastic tray placed in a peelable pouches autoclavable).

NB: It **isn't recommended** to realize the sterilization with peelables pouches **autoclavable** without a container Ancillaries can pierce the films components the pouches.

2-A sterilization **BY AUTOCLAVE** is advised following this protocol:

| Steam sterilization parameters | Values | Values |
|---|--------------------|--------------------|
| Cycle Type | Pre-vacuum | Pre-vacuum |
| Set Point Temperature | 134°C / 273°F | 132°C / 269.6°F |
| Cycle Time | Minimum 3 minutes | Minimum 4 minutes |
| Dry Time | Minimum 20 minutes | Minimum 20 minutes |
| Cool Time (inside and outside the autoclave chambers) | Minimum 40 minutes | Minimum 40 minutes |

3 - Remark: After each cycle of cleaning/sterilization, insure of all the good working of the ancillaries: verify the integrity and the correct operation, of the locking systems (clipping, prehension) without excessive play.

4- If the ancillaries utilization is different, after sterilization stock them in a tray on a peelable pouches.

In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent sterilization process for which it will be fully responsible.

STORAGE PRECAUTIONS

Ancillaries and trial prostheses must be stored in their dedicated container or equivalent packaging to prevent any deterioration.

Keep in a dry place away from light, at room temperature, in its original packaging.

INFORMATION REQUEST AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible.

In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time .

For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor, or COUSIN BIOTECH at spine@cousin-biotech.com.

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