# BIOMESH®

## CA.B.S'Air® Composite

### NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT

---

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	4
de	<u>Gebrauchsanweisung</u>	Seite	6
it	<u>Istruzioni per l'uso</u>	Pagina	
es	<u>Instrucciones de uso</u>	Página	8
pt	<u>Nota de instruções</u>	Página	10
ελ	<u>Οδηγίες χρήσης</u>	Σελίδα	12
da	<u>Brugsvejledning</u>	Side	14
fi	<u>Käyttöohjeet</u>	Sivu	16
nl	<u>Gebruiksaanwijzing</u>	Pagina	18
ro	<u>Instrucțiuni de utilizare</u>	Pagină	20
tu	<u>Kullanım Talimatları</u>	Sayfa	22

#### COUSIN BIOTECH

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**NOT188 / 180828**

Version du 28/08/2018

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

## Instructions for use

## NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT

**BIOMESH®CA.B.S.'AIR®COMPOSITE**

EXPANDABLE PARIETAL PROSTHESIS WITH INFLATABLE BALLOON FOR VISCERAL SURGERY

**SINGLE USE STERILE PRODUCT****DESCRIPTION**

The BIOMESH® CA.B.S.'Air® Composite prosthesis consists of:

- A parietal prosthesis made of an ePTFE and polypropylene plug, two or four polyethylene terephthalate with biocompatible colorant (D&C Green #6) and ePTFE fixation threads crimped with stainless steel needles (see on package).
- A balloon which allows the deployment of the parietal prosthesis (except for CABST1R05T).

**IMPLANTED MATERIALS :**

Polypropylene, Polyethylene terephthalate, ePTFE, Biocompatible colorant (D&C Green #6).

**MATERIALS IN CONTACT WITH THE PATIENT DURING THE SURGERY :**

Stainless steel, Poly vinyl chloride.

No human nor animal origin. Non resorbable.

**INDICATIONS:**

BIOMESH® CA.B.S.'Air® Composite is intended for use in all forms of hernia repair requiring reinforcement with non absorbable support material.

**PERFORMANCES**

The BIOMESH® CA.B.S.'Air® Composite prosthesis has a non-adhesive side, which is designed to avoid or reduce visceral adherences during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect for the mini invasive surgical treatment of inguinal, crural, umbilical and Linea Alba hernias.

Products, biocompatible and non resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a very high resistance to suture, have a great flexibility and allow a optimum and rapid integration and colonisation.

If necessary, the surgeon has the possibility of adding two additional points of fixing, according to his appreciation.

**CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to generate possible undesirable effects:

- Recurrence
- Infection
- Discomfort
- Adhesions

**PRECAUTIONS FOR USE**

BIOMESH® CA.B.S.'Air® Composite prosthesis are delivered sterile (ethylene oxide sterilization).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

**IMPORTANT : DO NOT REUSE - DO NOT RESTERILIZE****THIS BALLOON MUST IMPERATIVELY BE WITHDRAWN**

As specified on the product's labeling, the BIOMESH® CA.B.S.'Air® Composite mesh is for single use only. It can not be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

**STORAGE**

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

**IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

<b>SURGICAL PROCEDURE (except for CABST1R05T which comes without balloon)</b>					
<p><b>1-Installation:</b></p> <ul style="list-style-type: none"> <li>- The prosthesis equipped with its balloon crosses completely the hernial defect,</li> <li>-its positioning becomes then entirely retro parietal (intraperitoneal)</li> </ul>		<p><b>2-Deployment:</b></p> <ul style="list-style-type: none"> <li>- The expansion and the spreading of the prosthesis are obtained by inflating the balloon,</li> <li>- Parietal fixation is carried out by two sutures thanks to the pre-positioned threads</li> </ul>		<p><b>3- Prosthesis in place:</b></p> <ul style="list-style-type: none"> <li>-The prosthesis is fixed,</li> <li>-The balloon must imperatively be withdrawn</li> </ul>	

#### Description of the prosthesis:

This prosthesis is in a shape of a flat plug with 2 fixation threads crimped with needles (4 needles) for the product references CABST1R05T and CABS AIRC7, and with 4 fixation threads crimped with needles (8 needles) for the CABS AIRC09 product reference.

The aim of its innovating expansion balloon concept is to make dissection easier and to secure the unfolding of the prosthesis. This balloon is located inside the plug; it is filled-up with 30cm<sup>3</sup> of air. It acts as a pneumatic Jack ensuring complete crease-free unfolding of the prosthesis.

Its double-sided design combines a smooth ePTFE side which reduces viscera adhesions with a knitted monofilament polypropylene side which allows better colonisation of the peritoneum.

#### Sizes of the prosthesis:

The medical device is a sterile kit.

The product range is composed of:

CABS AIR 5cm diameter: for small hernias and ventral hernias on trocar site

CABS AIR 7cm diameter: for umbilical and linea alba hernias, as well as small ventral hernias < 2cm

CABS AIR 9cm diameter: for umbilical and linea alba hernias, as well as small ventral hernias that measure between 2 and 5cm

#### Installation of the patient :

Patient lying on his back

Both legs kept in the body's axis

Each arm placed on an arm wrestler in a 90° angle from the body's axis

#### 1-Incision :

The cutaneous incision (2 to 4cm according to the size of the defect) is centred on the hernia.

#### 2-Dissection and resection of the peritoneal sac :

The hernial sac is identified and its content is re-introduced back inside the abdominal cavity. The sac is then left opened in order to allow the Intra-peritoneal prosthesis to be put in place.

Next step consists in liberating epiploico-peritoneal adhesions around the edges of the hernia so that the prosthesis can fully expand.

#### 3-Positioning and fixating the prosthesis :

BIOMESH® CA.B.S.'Air® Composite is introduced through the hernial orifice after being rolled up into a tube shape or folded in an umbrella-like way.

The balloon is filled up in order to allow the expansion of the prosthesis and the placement of the prosthesis is checked so that it remains centred above the hernial orifice.

The prosthesis is fixed with the two pre-mounted threads that will transfix the peritoneum and the muscular wall.

Because of its Intraperitoneal placement, two additional points should be used to ensure secure fixation of the prosthesis and to avoid any epiploic tissue being caught in between the prosthesis and the wall.

The balloon is deflated and removed.

#### 4-Closure of the incision :

Musculo-aponeurotic closure in front of the prosthesis with resorbable suture thread. Sub-cutaneous and then cutaneous closure, normally without drainage.

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