BIOMESH[®] A2 A3 A4

NON RESORBABLE EXTRAPERITONEAL PARIETAL REINFORCEMENT IMPLANT

en	<u>Instructions for use</u>	Page 2
fr	<u>Notice d'instructions</u>	Page 3
de	<u>Gebrauchsanweisung</u>	Seite 4
it	<u>Istruzioni per l'uso</u>	Pagine 5
es	<u>Instrucciones de uso</u>	Pagina 6
ελ	<u>Οδηγίες χρήσης</u>	σελίδα 7
cs	<u>Navod k použití</u>	Strana 8
sk	<u>Pokyny k použití</u>	Strana 9

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NOT069/180725

Version du 25/07/2018

This release is the last update of the instructions of use and replace the previous edition

BIOMESH A2 A3 A4
Extraperitoneal reinforcement implant
STERILE PRODUCT FOR SINGLE USE

DESCRIPTION

BIOMESH® meshes are parietal reinforcement implants.

MATERIALS

Polyethylene terephthalate (see on the package)

Polyethylene terephthalate deep-dyed with dimethyl siloxane (see on the package)

Marked thread : Polyethylene terephthalate tinted with a biocompatible dye (See on the package).

Origin neither human nor animal – Non resorbable.

INDICATIONS

Repair and parietal reinforcement for inguinal, crural and ventral hernias.

PERFORMANCES

Meshes are adapted to both laparoscopic and conventional surgery.

The BIOMESH® device is designed for extraperitoneal implantation.

Meshes, biocompatible and non resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a shape memory, being cut-able, have a very high resistance to suture, have a great flexibility and allow an optimum and rapid integration and colonisation.

CONTRAINDICATIONS

Do not use in the following cases :

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy
- Surgery via a vaginal approach

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to involve possible undesirable effects:

- Discomfort
- Infection
- Recurrence

PRECAUTIONS FOR USE

BIOMESH® meshes are delivered sterile.

Before any use, inspect the integrity of packaging and device (of which blister / pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a **qualified** surgeon and trained to the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT : DO NOT REUSE - DO NOT RESTERILIZE

As indicated on the product labelling, the implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

STORAGE OF THE PROSTHESIS

To be stored in a dry place away from light and at room temperature in its original packaging.

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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