# 4D *Ventral*®

## SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT EXTRAPERITONEAL

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	4
de	<u>Gebrauchsanweisung</u>	Seite	6
it	<u>Istruzioni per l'uso</u>	Pagina	8
es	<u>Instrucciones de uso</u>	Pagina	10
pt	<u>Instruções de utilização</u>	Página	12
ελ	<u>Οδηγίες χρήσης</u>	Σελίδα	14
bg	<u>Инструкции за употреба</u>	страница	16
cs	<u>Návod k použití</u>	Strana	18
da	<u>Brugsvejledning</u>	side	20
fi	<u>Käyttöohjeet</u>	Sivu	22
nl	<u>Gebruiksaanwijzing</u>	Pagina	24
pl	<u>Instrukcja obsługi</u>	Strona	26
ro	<u>Instrucțiuni de utilizare</u>	Pagină	28
sk	<u>Návod na použitie</u>	Strana	30

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**NOT245/200204**

Version du 04/02/2020

Caution: Federal law (USA) restricts this device to sale,  
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT  
**EXTRAPERITONEAL**  
SINGLE USE STERILE PRODUCT

**DESCRIPTION**

4D Ventral® meshes are extraperitoneal partially resorbable parietal reinforcement implants.

**IMPLANTED MATERIALS**

Polypropylene (non-resorbable knitted monofilament) – Poly L Lactic Acid (resorbable knit)  
Origin neither human nor animal – Semi resorbable.

**INDICATIONS**

Repair of ventral hernia or other fascial defects that require the addition of an extraperitoneal reinforcing or bridging material to obtain the desired surgical result.

**PERFORMANCE**

The 4D Ventral® implant is designed for extraperitoneal implantation only.  
Meshes are biocompatible and partially resorbable. They have the advantage of being cut-able, of having a very high resistance to suture, of having a great flexibility.

**CONTRAINDICATIONS**

Do not implant intraperitoneally.

Do not use in the following cases:

- Allergy to any of the components
- Growing children
- Infected site
- Anticoagulant therapy
- Pregnancy

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
- Infection
- Inflammation
- Fitsula formation
- Recurrence
- Adhesion formation
- Erosion
- Seroma
- Hematoma
- Mesh migration
- mesh deformation
- Allergic reaction
- Foreign body reaction
- irritation nearby organ

**USAGE PRECAUTIONS**

4D Ventral® implants are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. Mesh fixation points should be at least 1cm from the edge of the mesh with 1cm spacing between fixation points.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

**IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

As specified on the product's labeling, the 4D Ventral® implants are for single use only. They can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

**EXPLANTATION AND ELIMINATION OF DEVICES**

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

## **STORAGE**

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

## **INFORMATION REQUESTS AND COMPLAINTS**

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country, and inform COUSIN BIOTECH immediately.

For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

## **IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.



## **SURGICAL PROCEDURE**

### **Extra-Peritoneal Hernia repair**

Example: Retrorectus mesh repair of a ventral hernia (Rives-Stoppa technique).

This technique utilizes the hernia sac to separate the mesh from the intra-abdominal contents.

- Superior to the umbilicus, dissection is performed above the posterior rectus fascia and under the rectus muscle.
- Below the umbilicus, dissection occurs in the preperitoneal space due to the lack of a posterior rectus sheath.
- The 4D Ventral® implant is placed in the newly formed pre-peritoneal space. The 4DVentral mesh has an oriented elasticity:

- Fixation of the mesh if needed, to the muscle layer above
- The peritoneum must be meticulously checked and no defect in the peritoneum should be left open to prevent any contact between the mesh and intra-abdominal contents.
- Redon's drains under aspiration are usually placed. The anterior sheath is closed. If tension is present, relaxing incisions may be used
- Skin closure.
- Using the Redon's drains, fill the pre-peritoneal space with an antiseptic solution and clamp the drains for 1 hour.
- A belt can be used to wrap the patient's abdomen for a few days to reinforce the abdominal wall.

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