SACROMESH SOFT PROLAPS'[®]

IMPLANT FOR THE TREATMENT OF PROLAPSE BY SACROCOLPOPEXY

en	Instructions for use	Page	2
fr	Notice d'instructions	Page	3
de	Gebrauchsanweisung	Seite	4
it	Istruzioni per l'uso	Pagina	5
es	Instrucciones de uso	Página	6
pt	Nota de instruções	Página	7
ελ	Οδηγίες χρήσης	Σελίδα	8
fi	Käyttöohjeet	Sivu	9
hu	Felhasználási útmutató	Oldal	10
lt	Naudojimo instrukcija	Puslapis	11
nl	Gebruiksaanwijzing	Pagina	12
pl	Instrukcja obsługi	Strona	13
ro	Instrucțiuni de utilizare	Pagină	14
sv	Bruksanvisning	Sida	15

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Made in France

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NOT123 / 200316
Version du 16/03/2020

This release is the last update of the instructions for use and replaces the previous edition

**SACROMESH SOFT PROLAPS®
UROGENITAL IMPLANT
Implant for the treatment of prolapse by sacrocolpopexy**

DESCRIPTION

Sacromesh Soft Prolaps® consist of knitted polypropylene monofilament.
There are intended for sacrocolpopexy of the pelvic prolapse of the women, by abdominal approach.

MATERIALS

100% polypropylene: single thread knitted ,weight 38 g/m² .
Marked thread : Polyethylene terephthalate tinted with a biocompatible dye (See on the package).
Origin neither human, nor animal _ Non resorbable.

INDICATIONS

Sacromesh Soft Prolaps® prosthesis are designed for use as reinforcement and support of tissues when treating genital prolapse surgically by sacrocolpopexy.
They can be used in either conventional or laparoscopic surgery.

PERFORMANCES

The biocompatible and non-resorbable prosthesis have the advantage of being easy to install, show good resistance to sutures and traction, with excellent flexibility, optimal and rapid integration and colonization.

CONTRAINDICATIONS

Do not use the implant in the following circumstances:

- Allergy to any of the components
- In a place with progressive or latent infection or signs of tissue necrosis
- For women receiving anticoagulant therapy
- Urinary or genital tract infection or non-aseptic conditions
- Pregnancy
- Growing children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

SIDE EFFECTS

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, urethra, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection

Over-correction may lead to temporary or permanent obstruction of the lower urinary tract and to retention of urine.

The known risks of surgery for the treatment of incontinence include hematoma, acute or chronic pain, scarification of the vagina, retraction or release of the mesh, infection, erosion, de novo urgency or voiding difficulties, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

SURGICAL PROCEDURE

Surgery under general anaesthetic.

Abdominal incision or installation of trocars.

Exposure of the lumbosacral angle.

Dissection of the prerectal fascia after incision of the anterior wall of Douglas' cul-de-sac. The dissection continues laterally as far as the posterior elevator muscles, using a vaginal valve if necessary.

The using mesh is fixed as distally as possible on the pelvic floor.

Dissection of the anterior face of vagina.

The anterior mesh is fixed by a running stitch, avoiding transfixation. Depending on the strip used and the surgeon's preference, the meshes is passed laterally relative to the uterine isthmus. The uterus is held in place, or two branches are branches are passed on either side of it (for precut mesh or mesh cut longitudinally by the surgeon).

Fastening the anterior and posterior mesh to the lumbosacral angle without excessive tension (presacral ligament)

Peritonization must be as complete as possible to avoid contact between the prosthesis and the intraperitoneal constituents.

PRECAUTIONS FOR USE

Sacromesh Soft Prolaps® devices must **always** be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in gynaecological and urological surgery) who has a good knowledge of the device, its intended use, the instruments and the surgical technique.

Sacromesh Soft Prolaps® implants are delivered sterile (sterilised using ethylene oxide), inspect the packaging to be sure it is intact (do not use if the labels and/or protective packaging are damaged).

Do not use if the device is damaged or out of date.

Do not reuse and do not resterilize the prosthesis.

STORAGE OF THE PROSTHESIS

To be stored in a dry place away from sunlight, at room temperature, in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

RETRIEVAL AND DISPOSAL OF DEVICES

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

INFORMATION REQUESTS AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country, and inform COUSIN BIOTECH immediately.

For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

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