IMPLANT BIOMESH®

NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT EXTRAPERITONEAL

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	4
de	<u>Gebrauchsanweisung</u>	Seite	6
it	<u>Istruzioni per l'uso</u>	Pagina	8
es	<u>Instrucciones de uso</u>	Pagina	10
pt	<u>Instruções de utilização</u>	Pagina	12
ελ	<u>Οδηγίες χρήσης</u>	σελίδα	14
bg	<u>Инструкции за употреба</u>	Страница	16
da	<u>Brugsvejledning</u>	Side	18
et	<u>Kasutusjuhend</u>	lehel	20
fi	<u>Käyttöohjeet</u>	Sivu	22
hu	<u>Felhasználási Útmutató</u>	Oldal	24
lv	<u>Lietošanas instrukcija</u>	Lappuse	26
nl	<u>Gebruiksaanwijzing</u>	Pagina	28
ro	<u>Instructiuni de Utilizare</u>	Pagină	30
tu	<u>Kullanım Kilavuzu</u>	Sayfa	32
sv	<u>Bruksanvisning</u>	Sida	34

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NOT029/200312

Version du 12/03/2020

This release is the last update of the instructions of use and replace the previous edition

IMPLANT BIOMESH®
NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT
SINGLE USE STERILE PRODUCT

DESCRIPTION

BIOMESH® meshes are non resorbable parietal reinforcement implants.

IMPLANTED MATERIALS

Mesh : polypropylene (knitted monofilament)

Marked thread : Polyethylene terephthalate tinted with a biocompatible dye (See on the package).

Origin neither human nor animal – Non resorbable.

INDICATIONS

Repair and parietal reinforcement for inguinal, femoral and ventral hernias.

PERFORMANCE

Meshes are particularly adapted either to laparoscopy surgery, or to conventional surgery.

The BIOMESH® is designed for extraperitoneal implantation only.

Meshes, biocompatible and non resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a shape memory, a very easy dimensioning, a very high resistance to suture, a great flexibility and allow an optimum integration and colonisation.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects:

- Discomfort
- Infection
- Recurrence

USAGE PRECAUTIONS

BIOMESH® meshes are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labeling, the BIOMESH® mesh is for single use only. It can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

STORAGE

To be stored in a dry place away from light and at room temperature in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be performed in accordance with the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis, following the current protocol. This protocol is available on demand from COUSIN BIOTECH. It is important to note that any implant that is not to be cleaned and disinfected before expedition must be contained in a sealed package. Explanted medical devices must be eliminated in accordance with the standards applicable in the country for the disposal of infectious hazardous waste. There are no specific recommendations for the elimination of a non-implanted device. The second paragraph is in case of an information request or complaint.

INFORMATION REQUEST AND COMPLAINT

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, he/she must inform COUSIN BIOTECH of this as soon as possible. In the case of failure of an implant or if it contributed to causing serious adverse effects for the patient, the health centre must follow the legal procedures in its country and inform COUSIN BIOTECH of this as quickly as possible. For any correspondence, please quote the reference and batch number, the reference coordinates and give a comprehensive description of the incident or complaint. Brochures, documentation and the surgical technique are available on request from COUSIN BIOTECH and its distributors.

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