# SOFT LIFT®

## STERILE UROGENITAL IMPLANTS: SUBURETHRAL SUPPORT TAPE

|    |                               |          |    |       |
|----|-------------------------------|----------|----|-------|
| en | Instructions for use          | Page     | 2  |       |
|    | <b>Surgical Procedure</b>     |          | 35 | 36-38 |
| fr | Notice d'instructions         | Page     | 5  |       |
|    | <b>Technique Opératoire</b>   |          | 35 | 36-38 |
| de | Gebrauchsanweisung            | Seite    | 8  |       |
|    | <b>Operationstechnik</b>      |          | 35 | 36-38 |
| it | Istruzioni per l'uso          | Pagina   | 11 |       |
|    | <b>Tēcnica Operatoria</b>     |          | 35 | 36-38 |
| es | Instrucciones de uso          | Pagina   | 14 |       |
|    | <b>Técnica Operatoria</b>     |          | 35 | 36-38 |
| pt | Nota de instruções            | Pagina   | 17 |       |
|    | <b>Técnica operatória</b>     |          | 35 | 36-38 |
| da | Brugsdvejledning              | Side     | 20 |       |
|    | <b>Operativ teknik</b>        |          | 35 | 36-38 |
| lt | Naudojimo instrukcija         | Puslapis | 23 |       |
|    | <b>Operatyvinė technika</b>   |          | 35 | 36-38 |
| nl | Gebruiksaanwijzing            | Pagina   | 26 |       |
|    | <b>Chirurgische procedure</b> |          | 35 | 36-38 |
| ro | Instrucțiuni de utilizare     | Pagină   | 29 |       |
|    | <b>Procedură chirurgicală</b> |          | 35 | 36-38 |
| sv | Bruksanvisning                | Sida     | 32 |       |
|    | <b>Kirurgiskt ingrepp</b>     |          | 35 | 36-38 |

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**NOT131/200316**

Version du 16/03/2020

This release is the last update of the instructions for use and replaces the previous edition

## Suburethral support tape SOFT LIFT® And tunnelling device associated (FANCOFIL, FANCOFEN)

### DESCRIPTION

The SOFT LIFT® device consists of a woven, single-thread polypropylene tape, extended at either end with traction treads for installation and peroperative tension adjustment.

The special COUSIN BIOTECH tunnelling device (FANCOFIL, FANCOFEN) has to be used to implant the SOFT LIFT® device.

A guide with implant is provided to guide the passage of the ancillary through the transobturator approach.

**Origin neither human, nor animal - Non Resorbable.**

### MATERIALS

• **Prosthesis:** Tape: 100% polypropylene (woven, single-thread)  
Traction threads: Polyethylene terephthalate

• **Tunnelling device:** Stainless steel, acetal copolymer

• **Guide:** Liquid crystal polymer + biocompatible dye.

### INDICATIONS

The SOFT LIFT® prosthesis is designed for use as a urethra support tape for the surgical treatment of feminine stress incontinence.

### CONTRA INDICATIONS

Not to be implanted in the following circumstances:

- Allergy to any of the component materials
- For woman receiving anticoagulant therapy
- Urinary tract infection
- Pregnancy
- Growing children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

### SIDE EFFECTS

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, urethra, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection

Over-correction may lead to temporary or permanent obstruction of the lower urinary tract and to retention of urine.

The known risks of surgery for the treatment of incontinence include hematoma, acute or chronic pain, scarification of the vagina, retraction or release of the mesh, infection, erosion, de novo urgency or voiding difficulties, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected

### WARNINGS

The SOFT LIFT® device must **always** be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in gynaecological and urological surgery) who has a good knowledge of the device, its intended use, the instruments and the surgical technique.

The SOFT LIFT® device is delivered sterile (sterilised using ethylene oxide), inspect the packaging to be sure it is intact (do not use if the labels and/or protective packaging are damaged). Do not use if the device is damaged or out-of-date.

Do not reuse and do not resterilize the prosthesis.

Guides can only be used once and should not be implanted.

To extract the tape of the packaging :

Roll the exterior white rubber band to free the green threads. Then, pull gently one of the green threads so that it slides and comes out of the tube.

### STORAGE

To be stored in a dry place away from sunlight, at room temperature, in its original packaging.

### EXPLANTATION AND DISPOSAL OF DEVICES

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

### INFORMATION REQUESTS AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

### DECONTAMINATION AND CLEANING OF ANCILLARIES

Decontamination and cleaning are performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

This step has to take place before **the first utilization** and **immediately after the utilization** to avoid adhesion of particles or dry secretion to the instruments

**The detachable instruments have to be dismantled**

↳ **Products advised for the cleaning**

- an enzymatic detergent neutral

- Or products adapted for the cleaning (neutral or soft alkaline) of surgical instruments in compliance with the regulations in effect.

- **Warning:** Don't use corrosive factor or caustic cleaning product

2 process are possible

**A- Chart Automated cleaning in a disinfectant washer with manual Pre-cleaning:**

|                | Step | Step Description      | Step Instruction                                                                                                                                                                                                                            | Accessories                                                                                                                                | Duration                                                              |
|----------------|------|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------|
| Cleaning Steps | 1    | Contamination Removal | Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush                                                                                                              | <ul style="list-style-type: none"> <li>Tap water</li> <li>Soft bristle brush<br/><b>(Do not utilize metal cleaning brushes)</b></li> </ul> | Until all visible soil is removed                                     |
|                | 2    | Pre-soak              | Prepare a detergent solution at ambient temperature (15-25°C / 59-77°F). Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position.      | <ul style="list-style-type: none"> <li>Detergent</li> <li>Water jet pistol</li> </ul>                                                      | Until product is visually clean                                       |
|                | 3    | Ultrasonic cleaning   | Clean the instrument in the ultrasonic tray including the solution.<br><b>Be careful every unclench connector or screws and bolts by the vibrations will be reassured.</b>                                                                  | <ul style="list-style-type: none"> <li>Ultrasonic tray</li> </ul>                                                                          | Minimum duration 10 minutes, this duration can vary according product |
| Cleaning Steps | 4    | Soak                  | Soak the instrument in tap water during 30 seconds. After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position. | <ul style="list-style-type: none"> <li>Reverse osmosis water</li> </ul>                                                                    | Minimum 30 seconds of soak                                            |
|                | 5    | Automated Washer      | Place entire disassembled device into the automated washer                                                                                                                                                                                  | <ul style="list-style-type: none"> <li>Automated Washer</li> <li>Disassembled Device</li> </ul>                                            | Minimum total cycle time: 39 minutes                                  |

| Instructions of cleaning in washer / disinfectant |            |                     |                         |
|---------------------------------------------------|------------|---------------------|-------------------------|
| Cycles parameters                                 | Time       | Minimum temperature | Type of detergent/water |
| Pre-cleaning                                      | 2 minutes  | <45°C / <113°F      | Tap water               |
| Cleaning                                          | 5 minutes  | 55°C / 131°F        | Detergent               |
| Draining                                          | 2 minutes  | <45°C / <113°F      | Tap water               |
| Thermic draining                                  | 5 minutes  | 90°C / 194°F        | Reverse osmosis water   |
| Drying                                            | 25 minutes | > 70°C / > 158°F    | Not applicable          |

|                 |   |              |                                                                                                                                                                  |                                                                                |                               |
|-----------------|---|--------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------|-------------------------------|
| Finishing Steps | 6 | Final Rinse  | <b>Rinse</b> carefully with <b>distilled water</b>                                                                                                               | <ul style="list-style-type: none"> <li>Distilled water</li> </ul>              | Minimum 1 minute              |
|                 | 7 | Final Drying | <b>Dry</b> device utilizing medical quality filtered air                                                                                                         | <ul style="list-style-type: none"> <li>Medical quality filtered air</li> </ul> | Until product is visually dry |
|                 | 8 | Tidying      | The instruments have to be placed <b>in the right position</b> inside the trays as described <b>on the position map fixed</b> on the tray for the sterilization. |                                                                                |                               |

**B- Chart Just manual cleaning :**

|                    | Step | Step Description         | Step Instruction                                                                                                                                                                                                                                                                                                                           | Accessories                                                                                                                                | Duration                                                                                                             |
|--------------------|------|--------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------|
| Cleaning Steps     | 1    | Contamination Removal    | Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush                                                                                                                                                                                                             | <ul style="list-style-type: none"> <li>Tap water</li> <li>Soft bristle brush<br/><b>(Do not utilize metal cleaning brushes)</b></li> </ul> | Until all visible soil is removed                                                                                    |
|                    | 2    | Drying                   | Dry the device utilizing a dry non-shedding wipe. Medical quality filtered air may be utilized if available                                                                                                                                                                                                                                | <ul style="list-style-type: none"> <li>Non-shedding wipe</li> <li>Medical quality filtered compressed air</li> </ul>                       | Until product is visually dry                                                                                        |
| Disinfection Steps | 3    | Disinfection Application | Prepare a low-foaming neutral enzymatic detergent solution, using tap water (15-25°C / 59-77°F). Soak the instrument in the open position (if possible) during minimum 1 minute. During the immersion, actuate mobile parts of the instrument minimum 3 times (if possible), so that the detergent can access all parts of the instrument. | <ul style="list-style-type: none"> <li>Detergent</li> </ul>                                                                                | Minimum duration 1 minute, this duration can vary according product                                                  |
|                    | 4    | Drying                   | After the soaking time, take out the instrument and wipe it with a disposable cloth. Then, place the instrument in a new bath of enzymatic detergent solution using warm tap water.                                                                                                                                                        | <ul style="list-style-type: none"> <li>Disposable cloth</li> <li>Detergent</li> </ul>                                                      | Until product is visually dry                                                                                        |
|                    | 5    | Manual Disinfection      | Brush all surface of the instrument using a soft brush for 2 minutes. Adjust brushing time if needed. Actuate mobile parts of the instrument 3 times (when applicable) and use a syringe to access all difficult areas. Use a volume of 60 mL of the detergent solution.                                                                   | <ul style="list-style-type: none"> <li>Soft bristle brush</li> <li>Detergent</li> <li>Syringe</li> </ul>                                   | Manual cleaning time duration is complete when the device's surface, joints, and crevices have been manually cleaned |
|                    | 6    | Rinse and                | Soak the instrument in reverse osmosis water for 1                                                                                                                                                                                                                                                                                         | <ul style="list-style-type: none"> <li>Reverse osmosis</li> </ul>                                                                          | Minimum 1 minute                                                                                                     |

|  |   |              |                                                                                                                                                                  |                                |                               |
|--|---|--------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------|-------------------------------|
|  |   | finition     | minute.<br>Use a syringe and 60mL of reverse osmosis water to access all difficult areas.<br>Repeat soaking 2 additional times using fresh water                 | - Syringe                      |                               |
|  | 7 | Final Rinse  | <b>Rinse</b> carefully with <b>distilled water</b>                                                                                                               | - Distilled water              | Minimum 1 minute              |
|  | 8 | Final Drying | <b>Dry</b> devise utilizing medical quality filtered air                                                                                                         | - Medical quality filtered air | Until product is visually dry |
|  | 9 | Tidying      | The instruments have to be placed <b>in the right position</b> inside the trays as described <b>on the position map fixed</b> on the tray for the sterilization. |                                |                               |

**Note:** In case of **suspicion of Creutzfeld-Jakob Disease (CJD)**, incinerate the ancillaries and the fixation pliers

In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

#### **STERILIZATION OF ANCILLAIRES**

Sterilization is performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

1- The sterilization in autoclave is to make in a **specific container** (fenced basket, Tray, Plastic tray placed in a peelable pouches autoclavable).

**NB:** It **isn't recommended** to realize the sterilization with peelables pouches **autoclavable** without a container → Ancillaries can pierce the films components the pouches.

2-A sterilization **BY AUTOCLAVE** is advised following this protocol:

| Steam sterilization parameters                        | Values             | Values             |
|-------------------------------------------------------|--------------------|--------------------|
| Cycle Type                                            | Pre-vacuum         | Pre-vacuum         |
| Set Point Temperature                                 | 134°C / 273°F      | 132°C / 269.6°F    |
| Cycle Time                                            | Minimum 3 minutes  | Minimum 4 minutes  |
| Dry Time                                              | Minimum 20 minutes | Minimum 20 minutes |
| Cool Time (inside and outside the autoclave chambers) | Minimum 40 minutes | Minimum 40 minutes |

3-Remark: After each cycle of cleaning/sterilization, insure of all the good working of the ancillaries: verify the integrality and the correct operation, of the locking systems (clipping, prehension) without excessive play.

4- If the ancillaries utilization is different, after sterilization stock them in a tray on a peelable pouches.

In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent sterilization process for which it will be fully responsible.

#### **PACKAGING, STORAGE AND STORING OF TUNNELLING DEVICE SOFT LIFT®**

The tunnelling device must be stored conditioned in order to avoid any deterioration.

Sterilization as well as the preceding stages such as cleaning and decontamination, are under the responsibility of the centers of care which must use a material and methods validated (Reduction of the risk of transmission of ATNC, circular DGS/R13/2011/449 of the 01/12/2011).

COUSIN BIOTECH recommends a steam sterilization. The instruments are compatible with a steam treatment with 134°C during 18min.

Always make sure that ancillaries are complete and in good working order after each cleaning/ sterilization cycle and before use.

#### **HANDLING OF TUNNELLING DEVICE SOFT LIFT®**

Before use, check that the packaging is not damaged in any way that might suggest compromise of the contents.

The tunnelling device must be stored in a container specifically designed for this use in equivalent packaging.

To prevent that the instruments are not in contact with objects which can deteriorate their surface.

Proceed to a visual examination of each instrument, in order to detect any possible deterioration.

The tunnelling device SOFT LIFT® are only suitable for use with original COUSIN BIOTECH instruments.

Moreover, the instruments must be used by a surgeon correctly trained to the implantation of prostheses SOFT LIFT®.

#### **Use of tunnelling device SOFT LIFT® during the intervention**

Refer to the following operative technique instructions for further information.

SOFT LIFT® tunnelling device are designed for use in implantation of the SOFT LIFT® prosthesis, using the traction threads at the ends of the tape.

During the procedure, these threads are attached to the ends of needles which are used to position the tape according to the pelvic anatomy prior to implantation of the device.

SOFT LIFT® tunnelling device come in different shapes according to the route of implantation: curved needles for the vaginal or supra-pubic route, and spiral needles for the in-out and out-in trans-obturator routes.

#### **Warnings on the tunnelling device**

-Tunnelling device SOFT LIFT® are medical devices of Class I, intended for a temporary use or reusable.

-Clean and decontaminate tunnelling device SOFT LIFT® after use, employing an appropriate method (to reduce the risk of non-conventional transmissible agent transmission).

If there is a possibility that the patient has contracted Creutzfeldt-Jakob Disease (CJD), proceed with an appropriate cleaning and decontamination method.

If the patient has contracted CJD, the instruments must be incinerated.

#### **IMPORTANT**

**For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.**

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