BIOMESH[®] A2 A3 A4

NON RESORBABLE EXTRAPERITONEAL PARIETAL REINFORCEMENT IMPLANT

en	Instructions for use	Page 2
fr	Notice d'instructions	Page 4
de	Gebrauchsanweisung	Seite 6
it	Istruzioni per l'uso	Pagine 8
es	Instrucciones de uso	Pagina 10
ελ	Οδηγίες χρήσης	σελίδα 12
cs	Pokyny k použití	Strana 14
sk	Návod na použitie	Strana 16
nl	Gebruiksaanwijzing	Pagina 18

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CE
1639

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NOT069/200331
Version du 31/03/2020

This release is the last update of the instructions of use and replace the previous edition

BIOMESH A2 A3 A4
Extraperitoneal reinforcement implant
STERILE PRODUCT FOR SINGLE USE

DESCRIPTION

BIOMESH® meshes are parietal reinforcement implants.

MATERIALS

Polyethylene terephthalate (see on the package)

Polyethylene terephthalate deep-dyed with dimethyl siloxane (see on the package)

Marked thread : Polyethylene terephthalate tinted with a biocompatible dye (See on the package).

Origin neither human nor animal – Non resorbable.

INDICATIONS

Repair and parietal reinforcement for inguinal, crural and ventral hernias.

PERFORMANCES

Meshes are adapted to both laparoscopic and conventional surgery.

The BIOMESH® device is designed for extraperitoneal implantation.

Meshes, biocompatible and non resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a shape memory, being cut-able, have a very high resistance to suture, have a great flexibility and allow an optimum and rapid integration and colonisation.

CONTRAINDICATIONS

Do not use in the following cases :

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy
- Surgery via a vaginal approach

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to involve possible undesirable effects:

- Discomfort
- Infection
- Recurrence

PRECAUTIONS FOR USE

BIOMESH® meshes are delivered sterile.

Before any use, inspect the integrity of packaging and device (of which blister / pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a **qualified** surgeon and trained to the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT : DO NOT REUSE - DO NOT RESTERILIZE

As indicated on the product labelling, the implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

STORAGE OF THE PROSTHESIS

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

“EXPLANTATION AND ELIMINATION OF DEVICES”

Explantation and handling should be performed in accordance with the recommendations of ISO 12891-1:2015 “Implants for surgery – Retrieval and analysis of surgical implants” Part 1: “Retrieval and Handling”. Any explanted device must be returned for analysis, following the current protocol. This protocol is available on demand from COUSIN BIOTECH. It is important to note that any implant that is not to be cleaned and disinfected before expedition must be contained in a sealed package. Explanted medical devices must be eliminated in accordance with the standards applicable in the country for the disposal of infectious hazardous waste. There are no specific recommendations for the elimination of a non-implanted device. The second paragraph is in case of an information request or complaint.

“INFORMATION REQUEST AND COMPLAINTS”

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, he/she must inform COUSIN BIOTECH of this as soon as possible. In the case of failure of an implant or if it contributed to causing serious adverse effects for the patient, the health centre must follow the legal procedures in its country, and inform COUSIN BIOTECH of this as quickly as possible. For any correspondence, please quote the reference and batch number, the reference coordinates and give a comprehensive description of the incident or complaint. Brochures, documentation and the surgical technique are available on request from COUSIN BIOTECH and its distributors.

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