PREMIUM®

NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	5
de	<u>Gebrauchsanweisung</u>	Seite	8
it	<u>Istruzioni per l'uso</u>	Pagina	11
es	<u>Instrucciones de uso</u>	Pagina	14
pt	<u>Nota de instruções</u>	Pagina	17
ελ	<u>Οδηγίες χρήσης</u>	Σελίδα	20
fi	<u>Käyttöohjeet</u>	Sivu	23
nl	<u>Gebruiksaanwijzing</u>	Pagina	26
ro	<u>Instrucțiuni de Utilizare</u>	Pagină	29
et	<u>Kasutusjuhend</u>	lehel	32
tu	<u>Kullanım Kilavuzu</u>	Sayfa	35

COUSIN BIOTECH

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Version du 03/04/2020

R_x
only

Caution : Federal law (USA) restricts this device to sale,
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

PREMIUM®
NON RESORBABLE PARIETAL REINFORCEMENT IMPLANT
SINGLE USE STERILE PRODUCT

DESCRIPTION

PREMIUM® meshes are non resorbable parietal reinforcement implants.

IMPLANTED MATERIALS

• Mesh : polypropylene (knitted monofilament). Blue thread (depending on references): monofilament polypropylene, non absorbable, colorant [phthalocyaninato(2-)] copper

Origin neither human nor animal – Non resorbable.

INDICATIONS

Repair of hernia or other fascial defects that require the addition of a reinforcing or bridging material to obtain the desired surgical result.

PERFORMANCE

Meshes are particularly adapted either to laparoscopic surgery, or to conventional surgery.

The PREMIUM® is designed for extraperitoneal implantation only.

Meshes are biocompatible and non resorbable, They have the advantage of having a shape memory, being cut-able, a very high resistance to suture, having a great flexibility and allow an optimum integration and colonisation.

CONTRAINDICATIONS

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

UNDESIRABLE SIDE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects:

- | | | | |
|---------------------|----------------|----------------------|-----------|
| - Discomfort / pain | - Recurrence | - Adhesion formation | - Erosion |
| - Infection | - Inflammation | - Fistsula formation | - Seroma |

USAGE PRECAUTIONS

PREMIUM® meshes are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. Mesh fixation points should be at least 1cm from the edge of the mesh with 1cm spacing between fixation points.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labeling, the PREMIUM® mesh is for single use only. It can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

STORAGE

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

DESCRIPTION OF THE SURGICAL PROCEDURES:**Open Inguinal hernia repair : Lichtenstein**

- The patient is placed in the supine position
- Local or general anaesthesia
- Iliac incision of approximately 5 cm
- Dissection is performed accordingly to the Lichtenstein technique
- The hernia sac is identified and reduced
- The transversalis fascia can be tightened with a non-absorbable suture. With direct hernias it helps to avoid any immediate post-operative recurrence
- PREMIUM® mesh can be cut to fit the patient's anatomy
- The prosthesis is progressively put in place
 - begin by placing PREMIUM® on the pubic spine
 - the upper part of the prosthesis is placed on the conjoint tendon

- carry placing the prosthesis towards the deep inguinal ring so that the first arm is placed behind
- then the second arm is placed so that it encircles the spermatic cord
- Fix the PREMIUM® mesh with non-resorbable sutures
 - apex is sutured to the pubic tubercle
 - the lower border of the mesh is sutured to the free edge of the inguinal ligament
 - continuous suture extends up just medial to the anterior superior iliac spine
 - suture the two tails of the mesh together around the spermatic cord
 - the infero-medial corner of the mesh is attached overlapping the pubic tubercle.
 - the mesh is anchored to the conjoined tendon by interrupted sutures
- the cord then rests on the prosthesis
- suction drain can be placed beneath the external oblique aponeurosis, especially in large inguinal hernias, where an extensive dissection was performed
- the aponeurosis of external oblique is then closed with absorbable sutures
- before closure of the surgical incision, its edges are infiltrated with a long-acting local anesthetic, such as Naropein®
- During peri-operative care of the patient, prophylactic antibiotics is usually given for 48 – 72 hours postoperatively

Laparoscopic procedure (TAPP or TEP)

1/ Example of a TAPP repair

- Local or general anaesthesia
- A pneumoperitoneum is created in the usual fashion
- intra-abdominal findings are reported [intra-abdominal pathology and inguinal hernia defects and sacs].
- The two additional Trocars are inserted under direct vision.
- Dissection is initiated and carried out following classical TAPP technique
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected.
- The PREMIUM® mesh is rolled like a cigarette and inserted through the Trocar into the intra-abdominal cavity and deployed over the inguinal region. The blue thread on the oval mesh is the anatomic landmark for the pubic bone. For the meshes with flaps, the cross shows the upper flap. The PREMIUM® mesh is attached to Cooper's Ligament, around and lateral to the Inferior Epigastric Vessels using tacks, staples or other mesh fixating device, according to the Surgeon's usual technique. The operator should be meticulous to avoid the triangle of doom, triangle of pain, iliac vessels and to place the fixations lateral to the inguinal ring.
- The operator should check the PREMIUM® mesh is well anchored to the surrounding structures.
- The peritoneum is closed meticulously and no defect between the peritoneum and the abdominal wall should be left open. In addition, it should cover the entire Mesh.
- The trocars are removed under direct vision. The fascia of the sub-umbilical trocar site is closed as needed.

2/ Example of a TEP repair

- General anaesthesia
- Create the pro-peritoneal space.
- The 10 mm trocar is inserted.
- Create the pneumo-pro-peritoneum.
- Two 5 mm ports are inserted under direct vision.
- If a peritoneal defect is created at this time it should be immediately closed.
- The anatomical landmarks are identical to when performing a TAPP repair.
- The indirect inguinal hernia sac should be dissected carefully from the Spermatic Cord. Direct hernia sacs are easily dissected. This is done by gentle traction with atraumatic graspers.
- The PREMIUM® mesh is rolled like a cigarette and inserted into the pro-peritoneal space and deployed over the inguinal region. The blue thread on the oval mesh is the anatomic landmark for the pubic bone. For the meshes with flaps, the cross shows the upper flap.
- The PREMIUM® mesh is first attached or secured to Cooper's Ligament, and the superior aspect of the pubic ramus. The PREMIUM® mesh is then fixed on the posterior aspect of the Linea Alba. The anchoring is continued around and lateral to the Inferior Epigastric Vessels. A few fixations are used to fix the PREMIUM® mesh lateral to the internal ring.
- The PREMIUM® mesh should gently and generously cover the Iliac Vessels without major gaps.
- The operator should check the Mesh is well anchored to the surrounding structures. The Mesh should not migrate and remain in place.
- The pro-peritoneal space will be checked for any peritoneal defects. If any, they should be closed. The 5 mm trocars are removed under direct vision. The pro-peritoneal space will collapse.
- If there are any questions about a missed peritoneal defect, a completion laparoscopy should be performed.
- Trocars are removed and the fascial defect closed with the appropriate suture. The skin edges are approximated in the usual manner.

Open Extra-Peritoneal Ventral Hernia repair

Example: Retrorectus mesh repair: Rives-Stopppa technique.

This technique utilizes the hernia sac to separate the mesh from the intra-abdominal contents. Superior to the umbilicus, dissection is performed above the posterior rectus fascia and under the rectus muscle.

Below the umbilicus, dissection occurs in the preperitoneal space due to the lack of a posterior rectus sheath. The PREMIUM® mesh is placed in the newly formed pre-peritoneal space, and fixated to the muscle layer above. The peritoneum must be meticulously checked and no defect in the peritoneum should be left open to prevent any contact between the mesh and intra-abdominal contents. Redon's drains under aspiration are usually placed.

The anterior sheath is closed. If tension is present, relaxing incisions may be used.

Skin closure. Using the Redon's drains, fill the pre-peritoneal space with an antiseptic solution and clamp the drains for 1 hour. A belt can be used to wrap the patient's abdomen for a few days to reinforce the abdominal wall.

Description of the prosthesis:

PREMIUM® meshes are constructed from knitted polypropylene monofilament with a diameter of approximately 0.1 mm.

PREMIUM® meshes are specifically indicated for laparoscopic and open surgery.

PREMIUM® exists in different sizes and shape, it is adapted to surgeon's technique and to patient's anatomy.

Size of the prostheses:

The PREMIUM® product range is composed of meshes of various sizes and shapes :

- Pre-cut meshes with keyhole for LICHTENSTEIN inguinal hernia repair
- Oval, pre-cut and rectangular meshes for TAPP/TEP inguinal hernia repair
- Square and rectangular meshes for extra-peritoneal open ventral hernia repairs

EXPLANTATION AND DISPOSAL OF DEVICES

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUESTS AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

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[İçindekiler](#)

This image shows a single page of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a full page of blank white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page, providing a template for writing or drawing. There are no margins, text, or other markings on the paper.

CE
1639

en fr de it es pt	- Symbols used on labelling - Symboles utilisés sur l'étiquette - Bei Etiketten verwandete Symbole - Simboli utilizzati sull'etichetta - Símbolos utilizados en el etiquetaje - Símbolos usados na etiqueta	ελ fi nl ro et tu	- Σύμβολα που χρησιμοποιούνται στις ετικέτες - Merkinnäissä käytetyt symbolit - Op de etikettering gebruikte symbolen - Simboluri utilizate la etichetare - Märgistusel kasutatud sümbolid - Etikette kullananl semboller	 e-IFU www.cousin-biotech.com/ifu  IFU can be downloaded through the QR code. A hard copy can be sent within 7 days on request by email to ifurequest@cousin-biotech.com or by using the order form on our website.
----------------------------------	--	----------------------------------	--	---

					
en fr de it es pt ελ fi nl ro et tu	Batch number Numéro de lot Chargennummer Numero di lotto Número de lote Número do lote Αριθμός παρτίδας Eränumero Batchnummer Numar lot Partii number Parti numarası	Reference on the brochure Référence du catalogue Verweis auf der Broschüre Referimento del catalogo Referencia al folleto Referência do catálogo Αριθμός καταλόγου Viittaus esitteessä Referentie in de brochure Referinte in brosură Viide brošüüris Katalog numarası	Manufacturer Fabricant Hersteller Produttore Fabricante Fabricante Κατασκευαστής Valmistaja Fabrikant Producator Tootja Üretici	Caution (see instructions for use) Attention, voir notice d'instructions Siehe Gebrauchsanweisung Vedere manuale istruzioni Veanse las instrucciones de uso Atenção, consultar nota de instruções Δείτε τις οδηγίες χρήσης Katso käyttöohjeet Zie gebruiksinstructies Pentru utilizare vezi instructiuni Ettevaatust (vaata kasutusjuhendit) Kullanım talimatlarına başvurun	Do not re-use Ne pas réutiliser Nicht wiederverwenden Non riutilizzare No volver a utilizar Não reutilizar Να μην επαναχρησιμοποιείται Älä käyttää uudelleen Niet hergebruiken A nu se reutiliza Mitte taaskasutada Yeniden kullanmayın

			
en fr de it es pt ελ fi nl ro et tu	Keep in a dry place Conserver au sec Trocken lagern Conservare in un luogo asciutto Manténgase en un lugar seco Manter em lugar seco Φυλάξτε το προϊόν μακριά από την υγρασία Säilytä kuivassa Op een droge plaats bewaren A se pastra la loc uscat Opbevaars på et tørt sted Kuru bir yerde saklayın	Keep away from sunlight A stocker à l'abri de la lumière du soleil Vor Sonnenlicht geschützt lagern Tenere lontano dalla luce solare Mantener alejado de la luz solar Manter ao abrigo da luz solar Προφυλάξτε το προϊόν από το ηλιακό φως Pidä poissa auringonvalolta Vrij van zonlicht houden A se pastra ferit de lumina soarelui Hoida eemal päikesevalgusest İşığa maruz bırakmayın	Do not use if the packaging is damaged Ne pas utiliser si l'emballage est endommagé Nicht verwenden, wenn Packung beschädigt ist Non utilizzare se la confezione è danneggiata No usar si el embalaje está dañado Não usar se a embalagem estiver danificada Να μην χρησιμοποιείται, εάν η συσκευασία έχει υποστεί κάποια βλάβη Älä käyttää, mikäli pakkaus on vaurioitunut Niet gebruiken indien de verpakking beschadigd is A nu se utiliza dacă ambalajul este deteriorat Mitte kasutada, kui pakend on kahjustatud Ambalaj hasarlıysa kullanmayın

			
en fr de it es pt ελ fi nl ro et tu	Use before: year and month Utiliser jusque : année et mois Verbrauchen bis : Jahr und Monat Utilizzare entro e non oltre: anno e mese Usese antes de: año y mes Utilizar até : ano e mês Ημερομηνία λήξης: έτος και μήνας Käytettävä ennen: vuosi ja kuukausi Te gebruiken vóór: jaar en maand A se utiliza înainte de: an luna Kasutada enne: aasta ja kuu Son kullanım tarihi: Yıl ve ay	Do not re-sterilize Ne pas restériliser Nicht erneut sterilisieren Non ristilizzare No volver e esterilizar Não reesterilizar Να μην επαναποστεριώνεται Älä steriloi uudelleen Niet opnieuw steriliseren A nu se reesteriliza Mitte korduvalt steriliseeda Yeniden sterilize etmeyin	Sterilized by ethylene oxide Produit stérile. Méthode de stérilisation : oxyde d'éthylène Mit Ethylenoxid sterilisiert Sterilizzato all'ossido di etilene Esterilizado con óxido de etileno Produto estéril. Método de esterilização: óxido de etileno Στείρο προϊόν. Μέθοδος αποστείρωσης: αιθυλενοξειδίο Steriloitu etyleenioksidilla Gesteriliseerd met ethyleenoxide Sterilizat cu oxid etilena Etüleenoksiid-steriliseerimine Etilen oksit ile sterilize edilmştir