# UNIVERSAL POLYMER CERCLAGE

## Sterile self-locking polymer cerclage system

---

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**NOT278 / 200417**

Version du 17/04/20

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

# UNIVERSAL POLYMER CERCLAGE

Sterile self-locking polymer cerclage system  
**SINGLE USE STERILE PRODUCT**

## **1-DESCRIPTION**

The "Sterile self-Locking Polymer Cerclage system" consists of a braided cerclage polymer cable with metal anchor plate and self-locking mechanism.

## **2- MATERIALS**

- Cable: polyethylene terephthalate multifilament braided cable, non absorbable, implanted
- Anchor: Titanium TA6V, non resorbable, implanted
- Wire guide : polyethylene terephthalate multifilament braided cable, non absorbable, not implanted

Origin neither human nor animal – Biocompatible. Non resorbable.

## **3-INDICATIONS**

- Repair of long bone fractures due to trauma or reconstruction;
- Reattachment of the greater trochanter in total hip arthroplasty, surface replacement arthroplasty, or other procedures involving trochanteric osteotomy

## **4-PERFORMANCE**

The role of cerclage device (with or without bone plate) is to hold the fragments of the bone in position till the bone heals.

In details, it works according three sequences:

- Ease of cerclage device setting around the bone (soft and conformable materials)
- Ease of hand pre-setting during initial fracture reduction to stabilise multi-fragmented bone elements
- Reliable tightening of the bone during healing process

The device is designed to be permanently implanted

## **5-CONTRAINDICATIONS**

- Skeletally immature patients;
- Active or suspected infection, either systemic or localized, in or around the implant site;
- Patient conditions, mental or neurological, that would tend to impact the patient's ability to follow physician's instructions during the post-operative healing phase;
- Compromised vascularity that would inhibit adequate blood supply to the fracture or operative site;
- Insufficient quality or quantity of bone, which would inhibit rigid device fixation;
- Any disease affecting the support and fixation of the prosthesis;

## **6-UNDESIRABLE SIDE EFFECTS**

- Cable failure, early or late loosening, wear, and change in position of one or more components
- Malunion of fracture / fracture displacement due to inadequate reattachment and/or early weight bearing;
- Progressive bone resorption and osteolysis
- Infection
- Neurological injury
- Vascular injuries due to the passer or the cable
- Foreign body reaction
- Allergy to a component of the implant
- Migration of the cerclage material

These undesirable side effects may occur a reoperation of the patient

## **7-POST-OPERATIVE PROCEDURES**

- The patient should be released from the hospital with complete instructions and warnings regarding further exercises and therapies coming from medical staff. The patient should be encouraged to report any unusual changes to his/her physician.
- The patient should be cautioned to govern his/her activity level. Postoperative therapy should be structured to prevent excessive loading.

## **8-PRECAUTIONS**

- The patient must be advised of the short- and long-term limitations of the procedure and the need to protect the implant from full weight bearing until adequate fixation and healing have occurred. Excessive activity and trauma affecting the reattachment and/or fixation have been implicated in failure of this procedure by loosening, fracture, and/or wear of the implants. Loosening of the components can result in increased production of wear particles, as well as damage to the bone, making successful revision surgery more difficult. The patient should be cautioned to limit activities, protect the bone from unreasonable stresses, and to follow the instructions of the physician with respect to follow-up care and treatment. The patient should be warned of the surgical risks and made aware of possible adverse effects. The patient should be warned that the device cannot and will not replicate the flexibility, strength, reliability or durability of normal healthy bone, that the implant can break or become damaged, particularly as a result of strenuous activity or trauma, and that the device has a finite expected service life and may need to be replaced in the future.
- Transient bacteremia can occur in daily life. Dental manipulation, endoscopic examination, and other minor surgical procedures have been associated with transient bacteremia. To prevent infection at the implant site, it may be advisable to use antibiotic prophylaxis before and after such procedures.
- An implant should never be re-used. Even though an implant may appear undamaged, previous handling and in-service stresses may have created imperfections that would reduce the service life of the implant.
- Patients should be questioned regarding sensitivity to metals. If there is a question pertaining to the patient's tolerance for titanium, appropriate testing should be performed
- The following conditions, singularly or concurrently, tend to impose severe loading on the affected bone, thereby placing the patient at higher risk for failure: obesity, heavy labor, active sports participation, high general activity level, likelihood of falls, alcohol or drug addiction, and other disabilities
- Cerclage cable breakage is reported in the literature. Possible causes of breakage include, but are not limited to, over tensioning the cable, fretting, and delayed or non-union of the fracture, osteotomy, or fusion site. In many instances, adverse effects may be clinically related rather than implant related. The surgical and post-operative management of the patient must be carried out with due consideration for all existing conditions of increased risk. Mental attitudes or personality disorders resulting in patient's failure to adhere to his/her physician's orders might delay post-operative recovery and aggravate adverse effects.
- Over-tensioning can create the potential for the cable to cut into the bone. Under-tensioning can cause inadequate fixation and increased fatigue stress.
- Improper placement, positioning and fixation of these devices may result in unusual stress conditions reducing the service life of the implant. Proper handling of any implant is mandatory.
- Use care in utilizing the Cable Passer and other instrumentation to avoid damaging neurovascular structures. Exercise caution when positioning and tightening the Cerclage Cables to avoid impinging or tethering neurovascular structures and causing attendant significant sequelae
- Prior to surgical use, a visual inspection of each implant for possible imperfections should be performed. Damaged or altered implants may produce undesirable stresses and cause defects, which could lead to failure.
- Before any use, inspect the integrity of packaging and device (of which pelables sachets). Do not use in the event of deterioration of the device and/or the packaging.
- Do not use if the device is out of date.
- This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and orthopaedic surgery). The surgeon is to be thoroughly familiar with the surgical procedure, instruments and implant characteristics prior to performing surgery.
- Twisting, kinking, cutting, notching, or scratching the braided surface may reduce the strength, fatigue resistance and/or wear characteristics of the implant system. These, in turn may induce internal stresses that are not obvious to the eye and may lead to fracture of the component.

### **MRI SAFETY:**

The implants are composed of non ferromagnetic materials and present a geometry non susceptible to generate induced currents. Moreover, as they are fixed to bone/or /tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular, in terms of heating and migration of implant has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded to a safety use for MRI scan between 1.5 and 3 Tesla.

NOTE : The device is often implanted with others metal devices. All the implanted devices must be considered for the MRI safety

As a precautionary measure, it is recommended to avoid MRI scans within the 48h of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential.

It has to be noted that the devices which present a high contrast with the biological environment can generate « artifacts » that has to be taken into account for the perfect execution and interpretation of imaging exams.

For this purpose, it has to be recommended to the patient who has this implant to warn as far as possible the concerned health professionals (radiologists and radiology operators) about the presence of this implant before these exams.

## **9-IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

As specified on the product's labeling, the UNIVERSAL POLYMER CERCLAGE device is for single use only. It can not be re-used and/or re-sterilised (potential risks would be and are not limited to loss of the product's sterility, to risk of infection, to loss of the product's efficiency, recurrence)

## **10-DESCRIPTION OF THE SURGICAL PROCEDURES:**

### **STEP 1. CERCLAGE SETTING**

It is important to note the direction of introduction of the Cable Passer around the bone as this affects the orientation of the cable and thus the orientation of the tensioning instrument relative to the user and the incision. Position the appropriate Cable Passer around the bone such that the distal end of the cannula emerges on the operating surgeon's side of the bone.

**Prior to Cerclage fixation, Remove carefully soft tissues around bone (or bone fragments) to secure a direct hard contact between polymer cerclage device and bone surface.**

Pass the Cable Passer around the bone from posterior to anterior. Insert the guide wire into the tip of the Cable Passer until the wire is seen exiting from the shaft of the Cable Passer

**Use care in passing the cable passer to avoid damage to soft tissues, bone surface and neurovascular structures (including muscle insertion line)**

**Remove carefully Cable Passer while Guide Wire remains in place, wrapped around the bone.**

Wrap Polymer Cerclage (anchor plate side) all the way through around the bone.

Cut Guide Wire and remove it entirely.

Pull through completely the open-end wire.

**Avoid wrapping the cables over sharp metal or bone graft surfaces.**

Pass fixed loop into the anchor plate. The auto-fix mechanism of the anchor will keep the connection while moving or rotating the cable.

**Maintain cerclage wires parallel. If twisted wires, unlock the anchoring, untwist wires and restart anchoring.**

### **STEP 2. HAND PRE-SETTING**

Fracture Reduction and Hand pre-setting

For fracture reduction, polymer cerclage must be used complementary to bone clamps. Bone clamps will address proper bone alignment and reduction with multi bone fragment.

If needed, several cables could be used. The order of cable insertion depends upon the fracture pattern, particular use or situation. The cables on each end of the plate may be wrapped around the bone to stabilize the bone while the final reduction is carried out.

The cables can be loosely tensioned (pre-setting) by hand as the self-locking mechanism will hold the tension.

The cable can be softly hand re-tight as required.

**Be aware of sharp age bone structure that could damage cerclage system.**

### **STEP 3. FINAL TIGHTENING**

Insert open-end wire into the Tensioner (first through the front nose then second into the rotational axis)

Position Tensioner in direct contact with buckle eyelet for stability while tensioning.

Note: for optimal tensioning, it is recommended to maintain tensioner well aligned with wire splice (reduce angle).

Rotate manually the Tensioner axis until desired compression load is achieved.

**Self-locking Mechanism will hold tension to allow re-tightening as desired**

To release the Tensioning instrument, De-Rotate (counterclockwise) the axis. That will allow to verify compression load around the bone.

As the self-locking mechanism hold the Tension, you can re-tight the cerclage as needed.

**Note : We recommend 3 Times Tightening process to ensure stable desired compression load**

**Do not cut free cables until complete operation is achieved.**

After removal of the tensioning instrument, use an electric blade to cut off at 2 cm length.

NOTE: The cable may be placed over smooth metal implants, such as bone plates.

NOTE: Complementary to bone plate, clinical literature recommends to adjust minimum of two cerclages devices (at both extremities) before cerclage final tensioning. For long bone fracture, it is described up to four cerclages to improve stability and secure healing process.

NOTE: With bone plate with screws, it is recommended to secure cerclage final tensioning (tightening) after screws fixation. This sequence will allow full operational performance for cerclage.

NOTE : Considering its specific design and for optimal effectiveness, it's recommended to place the polymer cerclage at the exterior surface of the bone plate (outside bone contact).

NOTE: Damaged instruments may lead to improper cable tension or implant position and result in implant failure. Thorough familiarity with the surgical technique is essential to ascertain their proper working condition.

NOTE: It is recommended to tension the self-locking polymer cerclage system with the tensioner developed by Cousin Biotech. The use of another tensioner must be approved by the company Cousin Biotech before use.

## **11-STORAGE PRECAUTION**

Use extreme care in handling and storage of implant components. Cable and anchor plate must be handled with care.

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

## **12-EXPLANTATION AND ELIMINATION OF THE PRODUCT**

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

## **13-INFORMATION REQUEST AND COMPLAINTS**

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country and inform COUSIN BIOTECH immediately.

For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors. If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor.

Patients may declare complaints to the manufacturer, distributor or to the local authorities.

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