PATCHS BIOMESH®

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Date de marquage CE: 02/1997

NOT076 / 200313
Version du 13/03/2020

BIOMESH® PATCH
Neurological patches N3, N3L
Cranial and spinal dura mater substitutes
STERILE PRODUCT FOR SINGLE USE

DESCRIPTION

The BIOMESH® neurological patches come in the form of non-resorbable non-woven meshes with a non-adherent side.

MATERIALS

Polyethylene terephthalate (polyester) soaked with dimethylsiloxane (silicone) and with a dimethylsiloxane (silicone) coating on one side.

No human or animal origin – Non-resorbable.

INDICATIONS

Replacement, protection and surgical reinforcement of the cranial and spinal dura mater.

PERFORMANCES

Products have a smooth, shiny, non-adherent side intended to prevent or limit adhesions to nerve tissue, and a rehabilitable side.

These neurological patches limit possible cerebrospinal fluid leak.

They avoid or limit nerve root and spinal cord compression by any bone graft and by scar tissue.

The neurological patches, which are biocompatible and non-resorbable, cause moderate reactional fibrosis. The neurological patches are very resistant, non-woven with mechanical bonding. They have the advantages of being micro-porous and tight, they are easily sized, have a strong resistance to sutures and are highly flexible.

They facilitate any reoperation by serving as a shield.

CONTRAINDICATIONS:

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Penetrating brain injury (PBI)

- Open Spina bifida

- Growing children

UNDESIRABLE SIDE EFFECTS:

This implant may lead to:

- Cerebrospinal fluid leak
- Cerebrospinal infection
- Extradural or subdural hematoma formation
- Cellular reaction with significant implant encapsulation
- Adhesions

PRECAUTIONS FOR USE

BIOMESH® patches are delivered sterile.

Before any use, ensure that the packaging and device are intact (including peelable pouches).

When attaching the device, it must not be stitched at less than 2 mm from the edge of the mesh.

The patch should not be under or oversized.

Do not use if the labels and/or the device and/or the packaging are damaged.

Do not use if the device is out of date.

COUSIN BIOTECH does not guarantee or recommend any particular brand of attachment device.

This device must be implanted only by a qualified surgeon, trained in the use of the product (knowledge of anatomy and neurological surgery).

After inserting non-absorbable sutures, it is better to use fibrin glue to avoid any cerebrospinal fluid leak.

CAUTION: DO NOT REUSE - DO NOT RESTERILIZE
FACE THE SMOOTH SILICONE SIDE TOWARDS THE BRAIN OR THE DURA MATER

In accordance with the product labelling, BIOMESH® N3 and N3L implants are for single use; they must never be reused or re-sterilized (potential risks include, but are not limited to, loss of sterility of the products, risk of infection, loss of product efficacy, recurrence, etc.)

RETRIEVAL AND DISPOSAL OF DEVICES

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2011 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

STORAGE

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

SURGICAL TECHNIQUE

In the case of replacement and/or reinforcement of the dura mater, the size of the patch should be as close as possible to the size of the loss of substance. The patch should completely cover the defect and extend beyond. The patch should not be stretched to cover the dural defect. Inadequate coverage may predispose the defect to the formation of adhesions and lead to cerebrospinal fluid leakage.

If the patch is too small, the tissue or the material may be subjected to excessive stress, which can lead to leakage or disinsertion around the suture. If the tissue is too large, excessive folds can form and cause undesirable tissue attachments.

The patch can therefore be cut by the surgeon to be perfectly adapted to the defect.

The patch can be sutured with non-absorbable thread equipped with atraumatic needles. When attaching the device, sutures should not be placed less than 2 mm from the edge of the patch. Use fibrin glue to seal the device completely.

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor, or write directly to: contact@cousin-biotech.com.

EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be performed in accordance with the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis, following the current protocol. This protocol is available on demand from COUSIN BIOTECH. It is important to note that any implant that is not to be cleaned and disinfected before expedition must be contained in a sealed package. Explanted medical devices must be eliminated in accordance with the standards applicable in the country for the disposal of infectious hazardous waste. There are no specific recommendations for the elimination of a non-implanted device. The second paragraph is in case of an information request or complaint.

INFORMATION REQUEST AND COMPLAINT

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, he/she must inform COUSIN BIOTECH of this as soon as possible. In the case of failure of an implant or if it contributed to causing serious adverse effects for the patient, the health centre must follow the legal procedures in its country and inform COUSIN BIOTECH of this as quickly as possible. For any correspondence, please quote the reference and batch number, the reference coordinates and give a comprehensive description of the incident or complaint. Brochures, documentation and the surgical technique are available on request from COUSIN BIOTECH and its distributors.

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