CORTICAL FIXATION SYSTEM

STERILE ANCHOR DEVICES

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	4
de	<u>Gebrauchsanweisung</u>	Seite	6
it	<u>Istruzioni per l'uso</u>	Pagina	8
es	<u>Instrucciones de uso</u>	Pagina	10
pt	<u>Nota de instruções</u>	Pagina	12
ελ	<u>Οδηγίες χρήσης</u>	σελίδα	14
bg	<u>Инструкции за употреба</u>	Страница	16
tu	<u>Kullanım Kilavuzu</u>	Sayfa	18

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Version du 08/04/2020

Rx
only

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician

Instructions for use

**CORTICAL FIXATION SYSTEM
STERILE ANCHOR DEVICES
SINGLE USE STERILE PRODUCT**

DESCRIPTION

Cortical Fixation System for tendons and ligaments.

The cortical fixation system comprises a continuous loop for graft suspension and a titanium button for the cortical support.

MATERIALS***- Materials implanted:***

Cortical button: Non-ferromagnetic titanium alloy (Ti TA6V ELI)

Continuous loop: 100% polyethylene terephthalate (PET)

Extended plate for cortical fixation system (optional): Non-ferromagnetic titanium alloy (Ti TA6V ELI)

Neither animal nor human origin. Non-absorbable

- Materials in contact with the patient during the operation:

Traction threads: Siliconized polyethylene terephthalate + biocompatible colorant D&C green#6

INDICATIONS

The Cortical Fixation Systems are used for fixation of tendons and ligaments during orthopaedic reconstruction procedures such as anterior cruciate ligament (ACL) reconstruction.

PERFORMANCE

This system facilitates ligament or tendon repair/reconstruction by graft fixation in the bone.

It has been approved for long-term implantation in the human body.

CONTRAINDICATIONS

- Known hypersensitivity to the implant material. Where material sensitivity is suspected, appropriate tests should be made and sensitivity ruled out prior to implantation.
- Insufficient quantity or quality of bone.
- Blood supply and previous infections which may tend to retard healing.
- Active infection.
- Conditions which tend to limit the patient's ability or willingness to restrict activities or follow directions during the healing period.
- Not to be used for the bone – patella tendon – bone technique.
- Do not implant in children while they are growing.

POSSIBLE ADVERSE EFFECTS

Possible complications linked to the implant:

- 1- Breaking of the loop, breaking of the traction wires, breaking of the button
- 2- Failure when placed under tension
- 3- Migration or upturn of the button
- 4- Reaction to a foreign body

Possible complications linked to insertion of the implant:

- 1- Pain in the incision area or surgical site
- 2- Local and transient inflammatory reactions
- 3- Infection
- 4- Wound dehiscence

WARNING

1- Any decision to remove an implant must take into account the fact that a second surgical operation may involve a risk for the patient. Removal must be followed by suitable postoperative care.

2- The patient must be informed of the limitations inherent in the nature and function of the implant, in particular regarding the weight and constraints placed by the body on the system before the bone has sufficiently healed.

3- Ensure that the bone tunnel has the correct length and diameter.

PRECAUTIONS BEFORE USE

Before use, check the integrity of the implant and the packaging.

Do not use if the implant, labels and/or protective packaging is damaged.

Do not use the product after the expiry date indicated on the label. This could compromise the performance of the system and the patient's safety.

Unpacking and handling of the implant must be carried out in accordance with asepsis standards, to ensure sterility.

When opening the cardboard box, the surgeon must check that the batch number and model match to those indicated on the label on the inside cover. If they do not, use of the implant is strictly forbidden.

The device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

This implant can be used with standard instruments available on the market and specified for orthopaedic surgery for ligament or tendon anchorage.

Surgical instruments are subject to wear during normal use. Instruments that have been used many times or been subject to very high stress are likely to break. Instruments must only be used for their intended purpose. COUSIN BIOTECH recommends that instruments should be regularly inspected for wear and bending.

DO NOT REUSE

DO NOT RESTERILIZE

The implant is for single use, as specified on the label. It must not on any account be reused or re-sterilized (potential risks include, but are not limited to: loss of sterility of the product, risk of infection, loss of effectiveness of the product).

COUSIN BIOTECH does not provide any guarantee or recommendations regarding the use of a specific make of fixation device. The properties of these devices are subject to alteration by the manufacturer, over which COUSIN BIOTECH has no control.

This device must only be implanted by a qualified surgeon trained in the use of the product. Avoid using staples or sharp-edged screws.

SURGICAL TECHNIQUE

Preliminary information

The system is optimized for use of the 4-strand semitendinosus technique with retrodrilling and the all-inside and inside-out techniques.

In case of retrodrilling, COUSIN BIOTECH recommends a diameter of 5mm for the tunnel to insert the implant. In case of a full tunnel, it is preferable to use the additional extended plate.

Use of the device with a tunnel diameter of less than 5mm has not been validated. If the diameter is greater than 5mm, the surgeon can decide whether or not to use the extended plate.

Technique

1. Suture the graft and put it in the loop.
2. Create a bone tunnel to insert the implant with a diameter of 5mm (recommended diameter) and a graft cavity suited to its diameter.
3. Choosing the size of the implant. Total length of the femoral tunnel (between the two cortices) - length of the graft tunnel +5cm. If the result is not equal to the size of the implant, round up to the larger implant size.
4. Pass the threads through the bone tunnel.
5. Pull on one of the sutures until the button comes out of the bone tunnel through the cortical bone.
6. Tip over the button against cortical bone surface by pulling on the graft.
7. Fix the graft in the opposite bone tunnel using the required method.
8. Remove the suture thread(s) from the button.

In case of use of the extended plate

It may be necessary to use an extended plate to increase the contact surface of the button on the cortical bone (full tibial or femoral tunnel).

Repeat steps 1 to 6 above.

7. Remove one of the two sutures
8. Slide the extended plate under the button. Ensure that the remaining suture is on the split side to avoid hindering the interlocking.
9. Fix the graft in the opposite bone tunnel using the required method.
10. Remove the remaining suture thread from the button.

STORAGE PRECAUTIONS

Store in a dry place away from sunlight, at room temperature and in the original packaging.

EXPLANTATION AND DISPOSAL OF DEVICES

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUESTS AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

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