# INTRAMESH® T1

## STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS

en	<u>Instructions for use</u>	page	2
fr	<u>Notices d'instructions</u>	page	3
de	<u>Gebrauchsanweisung</u>	Seite	4
it	<u>Istruzioni per l'uso</u>	Pagina	5
es	<u>Instrucciones de uso</u>	Pagina	6
pt	<u>Nota de instruções</u>	Pagina	7
ελ	<u>Οδηγίες χρήσης</u>	σελίδα	8
bg	<u>Инструкции за употреба</u>	страни	9
da	<u>Brugsvejledning</u>	Side	10
fi	<u>Käyttöohjeet</u>	Sivu	11
hu	<u>Felhasználási útmutató</u>	Pagina	12
nl	<u>Gebruiksaanwijzing</u>	Pagina	13
pl	<u>Instrukcja użycia</u>	Strona	14
ro	<u>Instrucțiuni de Utilizare</u>	Pagină	15
sk	<u>Návod na použitie</u>	Strana	16
cs	<u>Návod k použití</u>	Strana	17
tu	<u>Kullanım Kilavuzu</u>	Sayfa	18
sv	<u>Bruksanvisning</u>	Sida	19

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**NOT009\_200603**

Version du 03/06/2020

This release is the last update of the instructions of use and replace the previous edition

**IMPLANT INTRAMESH® T1**  
**STERILE INTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS**  
**SINGLE USE STERILE PRODUCT**

**DESCRIPTION**

INTRAMESH® T1 implants are parietal reinforcement implants.

**IMPLANTED MATERIALS**

Polypropylene and ePTFE for anti-adherent side (See on the package).

Biocompatible ink with dimethyl siloxane (See on the package).

**Origin neither human nor animal – Non resorbable.**

**INDICATIONS**

Repair and parietal reinforcement for umbilical and ventral hernias.

**PERFORMANCE**

These products can be used in laparoscopic surgery or open surgery. INTRAMESH® T1 products have a no-adhesive side, which is designed to avoid or reduce visceral adhesences during intra peritoneal use or extra peritoneal use in the presence of a peritoneal defect. Marking in black ink, composed of a cross and dots is on the non-adhesive side. The meshes, biocompatible and non resorbable have the advantage of having a shape memory and a great flexibility.

**CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant could generate possible undesirable effects:

- |                   |                      |                     |                            |
|-------------------|----------------------|---------------------|----------------------------|
| - Discomfort/Pain | - Adhesion formation | - Hematoma          | - Foreign Body Reaction    |
| - Infection       | - Fistula formation  | - Mesh migration    | - Irritation near by organ |
| - Recurrence      | - Erosion            | - Mesh deformation  | - Obstruction              |
| - Inflammation    | - Seroma             | - Allergic reaction |                            |

**USAGE PRECAUTIONS**

INTRAMESH® implants are delivered sterile. Before any use, inspect the integrity of packaging and device (of which blister / pelables sachets). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted only by a qualified surgeon and trained to the use of the product (knowledge of anatomy and visceral surgery).

As indicated on the product labelling, the implant, is designed for single use. They should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

**IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

Do not cut out the mesh and suture it at least 1cm of the edge. The trocar to be used must be chosen according to the size of the implant to be used. This is to allow passage without excessive force. Always roll the implant with the ePTFE face inside for the passage in the trocar. For use in keyhole surgery, we recommend the use of non-notched atraumatic forceps.

**STORAGE**

Store in a dry place away from sunlight and at room temperature in its original packaging.

**EXPLANTATION AND DISPOSAL OF DEVICES**

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

**INFORMATION REQUESTS AND COMPLAINTS**

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

**IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor

**IMPLANT INTRAMESH® T1**  
**IMPLANT DE RENFORCEMENT PARIETAL INTRAPERITONEAL STERILE NON RESORBABLE**  
**PRODUIT STERILE A USAGE UNIQUE**

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