# 4D DOME®

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**NOT035 / 200907**

Version du 07/09/2020

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

## **Instructions for use**

# 4DDOME®

## **SEMI-RESORBABLE REINFORCEMENT PARIETAL IMPLANT SINGLE USE STERILE PRODUCT**

### **DESCRIPTION**

The 4DDOME® device consists of a dome and a semi-resorbable reinforcement mesh for anterior and posterior reinforcement of the groin.

### **MATERIAL**

Polypropylene (non resorbable) Poly-L-Lactic Acid (resorbable)

**Origin neither human nor animal origin. Semi-resorbable**

### **INDICATIONS**

Repair and surgical reinforcement of inguinal and femoral hernias.

### **PERFORMANCES**

These implants, biocompatible and semi-resorbable make a reactional fibrosis which replaces the reinforcement after an implantation of 6 months. They have the advantages to have a very resistance to the suture, a quick and optimum integration and colonization and to have a resistance to the abdominal pressure.

**The 4DDOME® device is particularly designed to be implanted extra-peritoneally.**

### **CONTRAINDICATIONS**

Do not implant in the following cases:

- allergy to any of the components
- infected site
- pregnancy
- growing child
- Anticoagulant therapy

### **POSSIBLE UNDESIRABLE EFFECTS**

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
  - Recurrence
  - Adhesions
  - Foreign body reaction
  - Infection
  - Inflammation
  - Fistula formation
  - Mesh deformation
  - Hematoma
  - Mesh migration
  - Allergic reaction
  - Irritation nearby organ
  - Erosions
  - Seroma
  - Lymphocele
- A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes (open or laparoscopic)

### **PRECAUTIONS FOR USE**

The dome should always be installed with an anterior reinforcement.

The dome should not be cut.

4DDOME® devices are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which blister / peel-away pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and trained on the use of the product (knowledge of anatomy and visceral surgery).

### **IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

As indicated on the product labelling, the 4DDOME® implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

### **EXPLANTATION AND ELIMINATION OF DEVICES**

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

### **STORAGE**

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

### **INFORMATION REQUESTS AND COMPLAINTS**

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country, and inform COUSIN BIOTECH immediately.

For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

### **IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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