NAJA™

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CE
1639

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NOT237 / 201001
Version du 01/10/2020

Caution: Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

NAJA™**Ligament System for Spine Stabilization
Single use sterile product****DESCRIPTION**

The NAJA™ device is a spinal implant system dedicated to the correction and stabilization of the spine during fusion of the instrumented levels.

It consists of:

- A braid allowing correction and stabilization of the levels, after tensioning
- A malleable lamella allowing guidance of the braid during surgery and attached to the braid through a shrink sleeve. The lamella and shrink sleeve are not intended for long term implantation and must be removed after placement of the implant.
- A connector designed to be fixed on a rod and holding the tension on the braid

Because of its conception, and in order to obtain symmetric constructions, two versions are available for each size of the implant: one for the left size, one for the right side.

MATERIALS

Braid	Polyethylene Terephthalate (PET)	For a long term implantation
Connector	Titanium alloy Ti-6Al-4V ELI	
Malleable lamella	Stainless Steel	Not for long term implantation
Shrink Sleeve	Polyethylene Terephthalate (PET)	

All used materials are non resorbable and non-human and non-animal origin.

INDICATIONS

The NAJA™ System is a temporary implant for use in orthopedic surgery. The system is intended to provide temporary stabilization as a bone anchor during the development of solid bony fusion and aid in the repair of bone fractures. The indications for use include the following applications:

1. Spinal trauma surgery, used in sublaminar or facet wiring techniques.
2. Spinal reconstructive surgery, incorporated into constructs for the purpose of correction of spinal deformities such as idiopathic and neuromuscular scoliosis in paediatric population from 8 years old.
3. Spinal degenerative surgery, as an adjunct to spinal fusions.

The NAJA™ System may also be used in conjunction with other medical grade implants made of titanium or cobalt chrome alloy whenever “wiring” may help secure the attachment of the other implants.

CONTRAINDICATIONS

Do not use NAJA™ in the following cases:

- General infection or located infection in the intervention area
- Inflammation located in the intervention area
- Pregnancy
- Immunodepression
- Any pathology affecting bone metabolism, such as osteopenia, osteoporosis, metastasis, tumor... can compromise the expected mechanical support for the implant, and therefore, the success of the surgery
- Insufficient tissues or muscles coverage
- Severe fractures such that segments may not be maintained in satisfactory proximate reduction.
- Suspected or confirmed sensibility to at least one of the components
- The overweight or obesity can lead to mechanical overload on the implant, compromising the success of the surgery
- Any mental state or activity level non compatible with the postoperative follow-up should compromise the success of the surgery
- The smoking patients are more likely to present an unsatisfactory bone fusion or non-union
- Any other contraindication to an instrumented surgery of the spine

The described contraindications can be absolute or relative and must be evaluated by the surgeon prior to decide to implant.

POTENTIAL UNDESIRABLE SECONDARY EFFECTS

In addition of the risks linked to an instrumented spine surgery, we can observe, in few cases, the following effects:

- Unusual sensations, discomfort or pain due to the presence of the implant
- Migration, deformation or breakage of one or many components of the implant, can cause neurological, motion complications or pain
- Infection, allergy or inflammatory reaction due to the components, or to the presence of microparticles (metallosis)
- Halting the growth of instrumented levels
- Partial or total lost of correction
- Intraoperative or postoperative damage, fissures or fractures, in case of insufficient bone density
- Pseudoarthrosis
- Pulmonary complications
- Feeding difficulties

- Respiratory problems
- Proximal junctional kyphosis (PJK)
- Wound dehiscence
- Blood loss
- Death

In case of undesirable secondary effect, a second surgery may be necessary, in order to remove or replace one or many implants.

PRECAUTIONS FOR USE

General precautions

NAJA™ devices must be implanted by spine specialists surgeons, trained to the surgical technique. The information contained in this notice is necessary but not sufficient to control the surgical technique.

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

NAJA™ devices are delivered sterile and are single use. They must not be resterilized nor reused (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

NAJA™ devices are made from titanium, they must be associated with titanium or cobalt chrome arthrodesis material. Any association with other materials will not be guaranteed by COUSIN BIOTECH.

NAJA™ devices are used with rods; it is necessary to ensure that the diameter is equal to the size of the implant, indicated on the label and engraved on the connector.

NAJA™ device is available in two versions, one for the left side, the other for the right side. It is necessary to respect this repartition. The version is indicated on the label and engraved on the connector. The letter L indicates left side and the letter R indicates right side.

Implants must be handled and stored with great care, in a dry place, protected from light and at room temperature.

Preoperative precautions

The implantation decision must be taken by the surgeon, after information and agreement of the patient, considering, indications, contraindications and potential secondary effects.

Before the surgery, it is necessary to check the availability and the compatibility of devices which will be associated (size, materials...).

Preoperative precautions

Before opening, it is recommended to check the packaging integrity. Any damage on packaging may jeopardize the implant sterility. In case of damage on the packaging, the implant must not be used, nor resterilized.

Any implant whose expiration date has passed must not be used.

After opening, it is recommended to check the implant integrity.

Any implant whose have a deteriorated surface must not be implanted.

The implant setting up must be done with the dedicated ancillaries, following the surgical technique recommended by COUSIN BIOTECH.

The malleable lamella and the shrink sleeve are not intended to a long term implantation and must be removed after implant set up.

Postoperative precautions

The patient must be informed of the postoperative cares and of precautions.

The use of a brace could be necessary but is not mandatory. This decision is the surgeon one.

The device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

EXPLANTATION AND DISPOSAL OF DEVICES

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUESTS AND COMPLAINTS CLAIMS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

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