# IMPLANT BIOMESH®

## STERILE EXTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANT

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	4
de	<u>Gebrauchsanweisung</u>	Seite	6
it	<u>Istruzioni per l'uso</u>	Pagina	8
es	<u>Instrucciones de uso</u>	Pagina	10
pt	<u>Instruções de utilização</u>	Pagina	12
ελ	<u>Οδηγίες χρήσης</u>	σελίδα	14
bg	<u>Инструкции за употреба</u>	Страница	16
da	<u>Brugsvejledning</u>	Side	18
et	<u>Kasutusjuhend</u>	lehel	20
fi	<u>Käyttöohjeet</u>	Sivu	22
hu	<u>Felhasználási Útmutató</u>	Oldal	24
lv	<u>Lietošanas instrukcija</u>	Lappuse	26
nl	<u>Gebruiksaanwijzing</u>	Pagina	28
ro	<u>Instrucțiuni de Utilizare</u>	Pagină	30
tu	<u>Kullanım Kilavuzu</u>	Sayfa	32
sv	<u>Bruksanvisning</u>	Sida	34

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Date de marquage CE: Oct 1996

**NOT029/210402**  
Version du 02/04/2021

This release is the last update of the instructions of use and replace the previous edition

**IMPLANT BIOMESH®**  
**STERILE EXTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS**  
**SINGLE USE STERILE PRODUCT**

**DESCRIPTION**

BIOMESH® meshes are non-absorbable parietal reinforcement implants.

**IMPLANTED MATERIALS**

Mesh : polypropylene (knitted monofilament) Blue thread (depending on references): monofilament polypropylene, non absorbable, colorant [phthalocyaninato(2-)] copper

**Origin neither human nor animal – Non-absorbable.**

**INDICATIONS**

Repair and parietal reinforcement for inguinal and femoral hernias and eventrations.

**PERFORMANCE**

Meshes are particularly adapted either to laparoscopy surgery, or to conventional open surgery.

**The BIOMESH® is designed for extraperitoneal implantation only.**

The mesh is biocompatible and non-absorbable, causing a fibrosis reaction which takes over the reinforcement after six months of implantation. It has the advantage of having an excellent shape memory, easy to uses dimensions, a very high resistance to suturing, and great flexibility.

The blue thread on preshaped meshes is the anatomic landmark for the pubic bone.

**CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant could generate possible undesirable side effects:

- |                   |                      |                     |                            |
|-------------------|----------------------|---------------------|----------------------------|
| - Discomfort/Pain | - Adhesion formation | - Hematoma          | - Foreign Body Reaction    |
| - Infection       | - Fistula formation  | - Mesh migration    | - Irritation near by organ |
| - Recurrence      | - Erosion            | - Mesh deformation  |                            |
| - Inflammation    | - Seroma             | - Allergic reaction |                            |

**USAGE PRECAUTIONS**

BIOMESH® meshes are delivered sterile (ethylene oxide sterilization).

Before any use, inspect the integrity of packaging and device (of which peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date.

COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

When the surgeon considers it necessary to do so, it is recommended to cut the implant outside the surgical site to avoid any risk of releasing particles into the patient.

This device must be implanted only by a qualified surgeon trained on the use of the product (knowledge of anatomy and visceral surgery).

**IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**

As specified on the product's labeling, the BIOMESH® mesh is for single use only. It can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

**STORAGE**

Store in a dry place away from sunlight and at room temperature in its original packaging.

**IMPORTANT**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

**EXPLANTATION AND ELIMINATION OF DEVICES**

Explantation and handling should be performed in accordance with the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis, following the current protocol. This protocol is available on demand from COUSIN BIOTECH. It is important to note that any implant that is not to be cleaned and disinfected before expedition must be contained in a sealed package. Explanted medical devices must be eliminated in accordance with the standards applicable in the country for the disposal of infectious hazardous waste. There are no specific recommendations for the elimination of a non-implanted device. The second paragraph is in case of an information request or complaint.

### **INFORMATION REQUEST AND COMPLAINT**

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high-quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, he/she must inform COUSIN BIOTECH of this as soon as possible. In the case of failure of an implant or if it contributed to causing serious adverse effects for the patient, the health centre must follow the legal procedures in its country and inform COUSIN BIOTECH of this as quickly as possible. For any correspondence, please quote the reference and batch number, the reference coordinates and give a comprehensive description of the incident or complaint. Brochures, documentation and the surgical technique are available on request from COUSIN BIOTECH and its distributors.

**R<sub>x</sub>**  
only **Caution : Federal law (USA) restricts this device to sale,distribution and use by or on the order of a physician.**

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