4D DOME[®]

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NOT035 / 220301
Version du 01/03/2022

R_x
only

Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

Instructions for use

4DDOME®

SEMI-RESORBABLE REINFORCEMENT PARIETAL IMPLANT SINGLE USE STERILE PRODUCT

DESCRIPTION

The 4DDOME® device consists of a dome and a semi-resorbable reinforcement mesh for anterior and posterior reinforcement of the groin.

MATERIAL

Polypropylene (non resorbable) Poly-L-Lactic Acid (resorbable)

Origin neither human nor animal origin. Semi-resorbable

INDICATIONS

Repair and surgical reinforcement of inguinal and femoral hernias.

PERFORMANCES

These implants, biocompatible and semi-resorbable make a reactional fibrosis which replaces the reinforcement after an implantation of 6 months. They have the advantages to have a very resistance to the suture, a quick and optimum integration and colonization and to have a resistance to the abdominal pressure.

The 4DDOME® device is particularly designed to be implanted extra-peritoneally.

CONTRAINDICATIONS

Do not implant in the following cases:

- allergy to any of the components
- infected site
- pregnancy
- growing child
- Anticoagulant therapy

POSSIBLE UNDESIRABLE EFFECTS

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- Discomfort/Pain
 - Recurrence
 - Adhesions
 - Foreign body reaction
 - Infection
 - Inflammation
 - Fistula formation
 - Mesh deformation
 - Hematoma
 - Mesh migration
 - Allergic reaction
 - Irritation nearby organ
 - Erosions
 - Seroma
 - Lymphocele
- A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes (open or laparoscopic)

PRECAUTIONS FOR USE

The dome should always be installed with an anterior reinforcement.

The dome should not be cut.

4DDOME® devices are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which blister / peel-away pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and trained on the use of the product (knowledge of anatomy and visceral surgery).

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As indicated on the product labelling, the 4DDOME® implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

EXPLANTATION AND ELIMINATION OF DEVICES

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery – Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste.

There are no specific recommendations for the disposal of a non-implanted device.

STORAGE

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

INFORMATION REQUESTS AND COMPLAINTS

Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country, and inform COUSIN BIOTECH immediately.

For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

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[Sommaire](#)

Patient information leaflet 4DDOME®

Item	Patient information leaflet
1	<p>a. The name of the device : 4DDOME®</p> <p>b. The model of the device : 4DDOM24GSR, 4DDOM30GSR, 4DDOM38GSR, 4DDOME24SR, 4DDOME30SR, 4DDOME38SR</p>
2	<p>a. The intended purpose of the device : Repair and surgical reinforcement of inguinal and femoral hernias.</p> <p>b. The kind of patient on whom the device is intended to be used : Adult male or female suffering from a symptomatic ventral hernia necessitating a surgical mesh repair with an extraperitoneal implantation method.</p>
3	<p>Special operating instructions for the use of the device : The device must be implanted only by a qualified surgeon and trained on the use of the product (knowledge of anatomy and visceral surgery).</p>
4	<p>a. The intended performance of the device :</p> <p>These implants, biocompatible and semi-resorbable make a reactional fibrosis which replaces the reinforcement after an implantation of 6 months. They have the advantages to have a very resistance to the suture, a quick and optimum integration and colonization and to have a resistance to the abdominal pressure. The 4D DOME® device is particularly designed to be implanted extra-peritoneally.</p> <p>b. Undesirable side effects : Like any implantable medical device, this implant could generate possible undesirable effects :</p> <ul style="list-style-type: none"> - Discomfort/pain - Infection - Recurrence - Inflammation - Adhesions - Fistula formation - Erosion - Seroma - Hematoma - Mesh migration - Mesh deformation - Allergic reaction - Foreign body reaction - Irritation near by organ - Lymphocele - A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes (open or laparoscopic)

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5	<p>When 4DDOME® Should Not be Used : Do not implant in the following cases: Allergy to any of the components, infected site, pregnancy, growing children, anticoagulant therapy</p>
6	<p>MRI Information : 4DDOME ® are non-conducting and non-magnetic devices. They present no risk in Magnetic Resonance Imaging (MRI) environment.</p>
7	<p>Precautions : The dome should always be installed with an anterior reinforcement. The dome should not be cut. 4D DOME® devices are delivered sterile (ethylene oxide sterilisation). Before any use, inspect the integrity of packaging and device (of which blister / peel-away pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned. This device must be implanted only by a qualified surgeon and trained on the use of the product (knowledge of anatomy and visceral surgery).</p>
8	<p>Materials included in the device: This surgical mesh is composed of Polypropylene (non resorbable) Poly-L-Lactic Acid (resorbable) Origin neither human nor animal origin. Semi-resorbable</p>
9	<p>Reporting a serious incident : If you encounter any serious incident with your 4DDOME® implant you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link serviceclients@cousin-surgery.com or materiovigilance@cousin-biotech.com and to the Therapeutic Goods Administration using the following website: http://www.tga.gov</p>
10	<p>Legal Manufacturer : Cousin Biotech, 8 rue de l'Abbé Bonpain, 59117 Wervicq-Sud, France</p> <p>Manufacturing Facility : Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France</p>

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille