BIORING[®]

REPLACEMENT KIT WITH PORT

en	<u>Instructions for use</u>	Page	2
fr	<u>Notice d'instructions</u>	Page	4
de	<u>Gebrauchsanweisung</u>	Seite	6
it	<u>Istruzioni per l'uso</u>	Pagina	8
es	<u>Instrucciones de uso</u>	Pagina	10
et	<u>Kasutusjuhend</u>	Lehekülg	12
lv	<u>Lietošanas instrukcija</u>	Lappuse	14
lt	<u>Naudojimo instrukcija</u>	Puslapis	16
en	<u>Patient information leaflet (Australian requirements)</u>	Page	18

**COUSIN BIOTECH**

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CE
1639

Date de marquage CE: 02/2009
NOT147/220301
Version du 01/03/2022

ADJUSTABLE GASTRIC BANDING BIORING®

REPLACEMENT KIT WITH PORT

STERILE SINGLE USE PRODUCT

DESCRIPTION

Kit including a catheter rincer, a catheter, a connector, a plug, a huber needle, a port and a connection system. The various parts of the kit were designed for fitting perfectly to the BIORING® adjustable gastric banding.

MATERIALS

The injection site is made of polyoxymethylene and baryum sulphate. The septum is made of silicone. The catheter is made of silicone (Dimethyl Siloxane) and baryum sulphate labelled with a biocompatible ink. The connection system is made of silicone (Dimethyl Siloxane) and baryum sulphate. The plug is made of stainless steel and silicone (Dimethyl Siloxane). The connector is made of stainless steel.

No human nor animal origin. Not resorbable

INDICATIONS

Kit to be used when a BIORING® gastric banding has been implanted and in the event of preoperative and postoperative complications, repeated surgery on patients following complications (infection of the port, leakage, catheter damaged).

MRI SAFETY

The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). This replacement kit is mainly made with dimethyl siloxane and other component which are non-ferromagnetic in these quantities according to literature. However, there is one component is stainless steel. Consequently, this replacement kit is considered MR conditional.

PERFORMANCES

The port is an X-Ray- opaque injection site and it allows per-cutaneous adjustment of the diameter of the gastric banding. The connection system secures the connection between the port and the catheter by an on-clipping system.

The injection site can be:

- buried and sutured to the aponeurosis of the Rectus abdominis muscle under the left costal edge
- buried and sutured sub-cutaneously under the left costal edge
- buried and sutured in pre-sternal position on the Xiphoid apophysis
- buried and sutured placed in a different position up to the liking of the surgeon

Plug: In the event of an infection at the junction between the port and the catheter, the port is replaced with the plug until a new port is fitted. The plug is connected to the free end of the catheter of the BIORING® adjustable gastric banding.

Connector: When the catheter on the BIORING® adjustable gastric banding is damaged (cut, leakage), the connector and the catheter from the kit can be used to replace the damaged system.

Always refer to the instructions for the BIORING® adjustable gastric banding marketed by COUSIN BIOTECH.

The product **must be** implanted only by a qualified surgeon and well formed for the use of it.

CONTRA-INDICATIONS

Do not use the implant in the following circumstances:

- Allergy to any of the components
- Growing children
- Do not implant to fit with a gastric band that is not a BIORING® adjustable gastric banding
- Inflammatory disease or congenital or acquired anomalies of the gastrointestinal tract.
- Severe organic disease.
- Potential bleeding disorder of the upper gastrointestinal tract or patients who regularly or systematically take aspirin or non-steroidal anti-inflammatory drugs.
- Diagnosis or preexisting symptoms or family history of connective tissue autoimmune disease.
- Pregnancy.
- Patients with a drug or alcohol addiction.
- Patients who refuse to accept the dietary restrictions required by the procedure.
- Serious psychiatric disorders.

POSSIBLE COMPLICATIONS

Complications that may arise from the use of an adjustable gastric band include risks associated with the anaesthesia, the surgical technique, and the degree of intolerance of the patient to any foreign body implanted in the body.

EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be done following recommendations of ISO 12891-1:2015 « Implants for surgery – Retrieval and analysis of surgical implants » Partie 1: « Retrieval and Handling ». Any explanted device must be sent back, for analysis, following the current protocol. This protocol is available on demand to COUSIN BIOTECH. It is important to note that any implant that should not have been cleaned and disinfected before expedition must be contained in a sealed package. The elimination of explanted medical device must be conducted in accordance with standards in the country for the disposal of infectious hazardous waste. The elimination of a non-implanted device is not the subject of specific recommendations. The second paragraph is in case of information request or complaints.

INFORMATION REQUEST AND CLAIM

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time. For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim. Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

Perioperative complications (due to the surgical technique):

- Lesions of the spleen
- Perforation of the stomach or oesophagus
- other bleeding

In the immediate postoperative period:

- Nausea, vomiting (especially in the first few days after the intervention, and when the patient eats more than recommended).
- Bloating due to gas
- Dysphagia, dehydration
- Thromboembolic complications
- Infectious complications
- Postoperative bleeding

Delayed postoperative complications:

- Dilation of the upper pouch
- Slippage of the band
- Gastroesophageal reflux
- Pyrosis
- Gastritis of the upper pouch
- Migration of the band
- Displacement of the injection port
- Constrictive obstruction (due for example to oedema, food, an inadequate calibration, slippage of the band, or to torsion or dilation of the upper pouch)
- Deflation of the band (due to a leak in the band, injection port, or connection tube)
- Rapid weight loss may lead to symptoms of malnutrition, anemia, and associated complications
- Infection of the injection port
- Constipation

PRECAUTION FOR USE AND FOR STORAGE

- The components of the implantable port and of the connecting system have to be used only with the BIORING® adjustable gastric banding.
- The band can be inflated immediately using physiological saline (NaCl 0.9%). ONLY use a Huber needle to fill-in or deflate the band. The volume, as well as the interval before filling, is at the discretion of the practitioner/surgeon (maximum volume, see the next table).
- Inflating of the band: maximum volume of physiological saline according to sizes of the bands

Band sizes (diameter : mm)	Maximum volumen of physiological saline (cm³)
Ø 21/24	10
Ø 27/29	10
Ø 36/38	15

An example of band-filling algorithm is available on request.

- To keep in a dry place without any sunlight, and in an ambient atmosphere, in its packaging of origin
- Please check the perfect integrity of the packaging (do not use in case of damaged packaging protection and/or damage of the label)
- Do not use in case of damaged product or if the product has ran out date
- **DO NOT RE-USE** **DO NOT RESTERILIZE**

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Item	Patient information leaflet
1	<p>a. The name of the device : Replacement kit with port</p> <p>b. The model of the device : BCBRIREPOU</p>
2	<p>a. The intended purpose of the device : Kit to be used when a BIORING® gastric banding has been implanted and in the event of preoperative and postoperative complications, repeated surgery on patients following complications (infection of the port, leakage, catheter damaged).</p> <p>b. The kind of patient on whom the device is intended to be used : Persons presenting morbid obesity and weight loss resistance, i.e. persons whose Body Mass Index [BMI, formula = weight in kilograms/(height in metres)²] is greater than or equal to 40 kg/m², or 35 to 40 kg/m² when combined with at least one comorbidity. The surgeon can prescribe the adjustable gastric band at his/her own discretion.</p>
3	<p>Special operating instructions for the use of the device : No special operating instructions for the use of the device</p>
4	<p>a. The intended performance of the device :</p> <p>The port is an X-Ray- opaque injection site and it allows per-cutaneous adjustment of the diameter of the gastric banding The connection system secures the connection between the port and the catheter by an on-clipping system The injection site can be:</p> <ul style="list-style-type: none"> - buried and sutured to the aponevrosis of the Rectus abdominis muscle under the left costal edge - buried and sutured sub-cutaneously under the left costal edge - buried and sutured in pre-sternal position on the Xiphoid apophysis - buried and sutured placed in a different position up to the liking of the surgeon <p>Plug: In the event of an infection at the junction between the port and the catheter, the port is replaced with the plug until a new port is fitted. The plug is connected to the free end of the catheter of the BIORING® adjustable gastric banding.</p>

Item	Patient information leaflet
	<p>Connector: When the catheter on the BIORING® adjustable gastric banding is damaged (cut, leakage), the connector and the catheter from the kit can be used to replace the damaged system.</p> <p>Always refer to the instructions for the BIORING® adjustable gastric banding marketed by COUSIN BIOTECH.</p> <p>The product must be implanted only by a qualified surgeon and well formed for the use of it.</p>
4	<p>b. Undesirable side effects :</p> <p>Complications that may arise from the use of an adjustable gastric band include risks associated with the anaesthesia, the surgical technique, and the degree of intolerance of the patient to any foreign body implanted in the body.</p>
5	<p>When the BIORING ® Replacement kit with port Should Not be Used :</p> <p>Do not use the implant in the following circumstances :</p> <ul style="list-style-type: none"> - Allergy to any of the components - Growing children - Do not implant to fit with a gastric band that is not a BIORING adjustable gastric banding - Inflammatory disease or congenital or acquired anomalies of the gastrointestinal tract. - Severe organic disease. - Potential bleeding disorder of the upper gastrointestinal tract or patients who regularly or systematically take aspirin or non-steroidal anti-inflammatory drugs. - Diagnosis or preexisting symptoms or family history of connective tissue autoimmune disease. - Pregnancy. - Patients with a drug or alcohol addiction. - Patients who refuse to accept the dietary restrictions required by the procedure. - Serious psychiatric disorders.
6	<p>MRI Information :</p> <p>The implant has not been tested for safety and compatibility in the case of postoperative investigation using Magnetic Resonance Imaging (MRI). This replacement kit is mainly made with dimethyl siloxane and other component which are non-ferromagnetic in these quantities according to literature. However, there is one component is stainless steel. Consequently, this replacement kit is considered MR conditional.</p>

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
7	<p>Precautions :</p> <p>The components of the implantable port and of the connecting system have to be used only with the BIORING® adjustable gastric banding. The band can be inflated immediately using physiological saline (NaCl 0.9%). ONLY use a Huber needle to fill-in or deflate the band. The volume, as well as the interval before filling, is at the discretion of the practitioner/surgeon (maximum volume, see the next table). Inflating of the band: maximum volume of physiological saline according to sizes of the bands</p> <table border="0" data-bbox="430 492 1308 716"> <thead> <tr> <th style="text-align: center;">Band sizes (diameter : mm)</th> <th style="text-align: center;">Maximum volumen of physiological saline (cm3)</th> </tr> </thead> <tbody> <tr> <td style="text-align: center;">Ø 21/24</td> <td style="text-align: center;">10</td> </tr> <tr> <td style="text-align: center;">Ø 27/29</td> <td style="text-align: center;">10</td> </tr> <tr> <td style="text-align: center;">Ø 36/38</td> <td style="text-align: center;">15</td> </tr> </tbody> </table> <p>To keep in a dry place without any sunlight, and in an ambient atmosphere, in its packaging of origin Please check the perfect integrity of the packaging (do not use in case of damaged packaging protection and/or damage of the label Do not use in case of damaged product or if the product has ran out date DO NOT RE-USE DO NOT RESTERILIZE</p>	Band sizes (diameter : mm)	Maximum volumen of physiological saline (cm3)	Ø 21/24	10	Ø 27/29	10	Ø 36/38	15
Band sizes (diameter : mm)	Maximum volumen of physiological saline (cm3)								
Ø 21/24	10								
Ø 27/29	10								
Ø 36/38	15								
8	<p>Materials included in the device : The injection site is made of polyoxymethylene and baryum sulphate. The septum is made of silicone. The catheter is made of silicone (Dimethyl Siloxane) and baryum sulphate labelled with a biocompatible ink. The connection system is made of silicone (Dimethyl Siloxane) and baryum sulphate. The plug is made of stainless steel and silicone (Dimethyl Siloxane). The connector is made of stainless steel. No human nor animal origin. Not resorbable</p>								
9	<p>Reporting a serious incident :</p> <p>If you encounter any serious incident with your implant you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link serviceclients@cousin-surgery.com or materiovigilance@cousin-biotech.com and to the Therapeutic Goods Administration using the following website: http://www.tga.gov</p>								
10	<p>Legal Manufacturer : Cousin Biotech, 8 rue de l'Abbé Bonpain, 59117 Wervicq-Sud, France</p> <p>Manufacturing Facility : Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France</p>								

Notes

Notes

Notes

en - Symbols used on labelling
fr - Symboles utilisés sur l'étiquette
de - Erklärung der Symbole auf der Verpackung
it - Simboli usati sull'etichetta
es - Símbolos utilizados en la etiqueta
et - Siltidel kasutatud sümbolid
lv - Etiketēje vaižduojami simboliai
lt - Marķējumā izmantotie simboli

www.cousin-biotech.com/ifu

IFU can be downloaded through the QR code.
A hard copy can be sent within 7 days
on request by email to
ifurequest@cousin-biotech.com
or by using the order form on our website.

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille