STERILE UROGENITAL IMPLANTS: SUBURETHRAL SUPPORT TAPE

en	Instructions for use	Page	2
	Surgical Procedure		38
fr	Notice d'instructions	Page	6
	Technique Opérateur		38
de	Gebrauchsanweisung	Seite	10
	Operationstechnik		38
it	Istruzioni per l'uso	Pagina	14
	Tecnica Operatoria		38
es	Instrucciones de uso	Página	18
	Técnica Operatoria		38
pt	Instruções de utilização	Página	22
	Técnica operatória		38
da	Brugsdvejledning	Side	26
	Operativ teknik		38
nl	Gebruiksaanwijzing	Pagina	30
	Chirurgische procedure		38
ro	Instrucțiuni de utilizare	Pagină	34
	Procedură chirurgicală		38

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Date de marquage CE : Aout 2002
NOT007/221019

Version du 19/10/2022

This release is the last update of the instructions for use and replaces the previous edition

SUBURETHRAL SUPPORT TAPE LIFT® **And tunnelling device associated (FANCOFIL, FANCOFEN)**

DESCRIPTION

The LIFT® device consists of a woven, single-thread polypropylene tape, extended at either end with traction treads for installation and peroperative tension adjustment.

The special COUSIN BIOTECH tunnelling device (FANCOFIL, FANCOFEN) has to be used to implant the LIFT® device.

Origin neither human, nor animal - Non Resorbable.

MATERIALS

- **Prosthesis:** Tape: 100% polypropylene (woven, single-thread)
Traction threads: Polyethylene terephthalate
- **Tunnelling device:** Stainless steel, acetal copolymer

INDICATIONS

The LIFT® prosthesis is designed for use as a urethra support tape for the surgical treatment of feminine stress incontinence.

CONTRA INDICATIONS

Not to be implanted in the following circumstances:

- Allergy to any of the component materials
- For woman receiving anticoagulant therapy
- Urinary tract infection
- Pregnancy
- Growing children

Any patient who could possibly be pregnant should be warned that pregnancy may abolish the benefit of the procedure.

SIDE EFFECTS

As with all medical implants, this product may lead to side effects, such as:

Perforation or damage to blood vessels, bladder, urethra, rectum, intraperitoneal organs or nerves.

As with all implants, local irritation and/or a reaction to the foreign body may occur.

At tissue level, reactions to the implant may include vaginal extrusion, erosion through the urethra or surrounding tissues, displacement of the device, fistula or inflammation. If any such reaction occurs, it may be necessary to completely remove the sling.

As with all foreign bodies, the polypropylene sling may aggravate a pre-existing infection

Over-correction may lead to temporary or permanent obstruction of the lower urinary tract and to retention of urine.

The known risks of surgery for the treatment of incontinence include hematoma, acute or chronic pain, scarification of the vagina, retraction or release of the mesh, infection, erosion, de novo urgency or voiding difficulties, dyspareunia, migration of the device, the total failure of the intervention leading to relapse. One or several repeated surgeries could also be necessary to treat those complications, as some complications are not always completely corrected.

WARNINGS

The LIFT® device must **always** be implanted by a qualified surgeon (familiar with the relevant anatomy and experienced in gynaecological and urological surgery) who has a good knowledge of the device, its intended use, the instruments and the surgical technique.

The LIFT® device is delivered sterile (sterilised using ethylene oxide), inspect the packaging to be sure it is intact (do not use if the labels and/or protective packaging are damaged). Do not use if the device is damaged or out-of-date.

Do not reuse and do not resterilize the prosthesis.

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

As specified on the product's labeling, the implants are for single use only. They can not be re-used and/or re-sterilised (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence).

Storage

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

EXPLANTATION AND DISPOSAL OF DEVICES

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery - Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

INFORMATION REQUESTS AND COMPLAINTS

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

DECONTAMINATION AND CLEANING OF ANCILLARIES

Decontamination and cleaning are performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

This step has to take place before **the first utilization** and **immediately after the utilization** to avoid adhesion of particles or dry secretion to the instruments

The detachable instruments have to be dismantled

Products advised for the cleaning

- an enzymatic detergent neutral
- Or products adapted for the cleaning (neutral or soft alkaline) of chirurgical instruments in compliance with the regulations in effect.

- **Warning:** Don't use corrosive factor or caustic cleaning product.

2 process are possible

Chart Automated cleaning in a disinfectant washer with manual Pre-cleaning:

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	<ul style="list-style-type: none"> - Tap water - Soft bristle brush (Do not utilize metal cleaning brushes) 	Until all visible soils are removed
	2	Pre-soak	Prepare a detergent solution at ambient temperature (15-25°C / 59-77°F). Afterwards lumina, threads, joints and gaps of the instruments have to be flushed with a water jet pistol for a minimum time of 10 seconds for each position.	<ul style="list-style-type: none"> - Detergent - Water jet pistol 	Until product is visually clean
	3	Ultrasonic cleaning	Clean the instrument in the ultrasonic tray including the solution. Be careful every unclench connector or screws and bolts by the vibrations will be reassured.	<ul style="list-style-type: none"> - Ultrasonic tray 	Minimum duration 10 minutes, this duration can vary according to product
Cleaning Steps	4	Soak	Soak the instrument in tap water during 30 seconds. After, ultrasonic treatment the lumina, threads, joints and gaps of the instrument have to be flushed again with a water jet pistol for a minimum time of 10 seconds for each position.	<ul style="list-style-type: none"> - Reserve osmosis water 	Minimum 30 seconds of soak
	5	Automated Washer	Place entire disassembled device into the automated washer	<ul style="list-style-type: none"> - Automated Washer - Disassembled Device 	Minimum total cycle time: 39 minutes

Instructions of cleaning in washer / disinfectant			
Cycles parameters	Time	Minimum temperature	Type of detergent/water
Pre-cleaning	2 minutes	<45°C / <113°F	Tap water
Cleaning	5 minutes	55°C / 131°F	Detergent
Draining	2 minutes	<45°C / <113°F	Tap water
Thermic draining	5 minutes	90°C / 194°F	Reserve osmosis water
Drying	25 minutes	> 70°C / > 158°F	Not applicable

Finishing Steps	6	Final Rinse	Rinse carefully with distilled water	<ul style="list-style-type: none"> - Distilled water 	Minimum 1 minute
	7	Final Drying	Dry device utilizing medical quality filtered air	<ul style="list-style-type: none"> - Medical quality filtered air 	Until product is visually dry
	8	Tidying	The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization.		

A- Chart Just manual cleaning :

	Step	Step Description	Step Instruction	Accessories	Duration
Cleaning Steps	1	Contamination Removal	Rinse product room temperature running tap water removing any visible organic material with assistance of a soft bristle brush	<ul style="list-style-type: none"> - Tap water - Soft bristle brush (Do not utilize metal cleaning brushes) 	Until all visible soils removed
	2	Drying	Dry the device utilizing a dry non-shedding wipe. Medical quality filtered air may be utilized if available	<ul style="list-style-type: none"> - Non-schedding wipe - Medical quality filtered compressed air 	Until product is visually dry
Disinfection Steps	3	Disinfection Application	Prepare a low-foaming neutral enzymatic detergent solution, using tap water (15-25°C / 59-77°F). Soak the instrument in the open position (if possible) during minimum 1 minute. During the immersion, actuate mobile parts of the instrument minimum 3 times (if possible), so that the detergent can access all parts of the instrument.	<ul style="list-style-type: none"> - Detergent 	Minimum duration 1 minute, this duration can vary according product
	4	Drying	After the soaking time, take out the instrument and wipe it with a disposable cloth. Then, place the instrument in a new bath of enzymatic detergent solution using warm tap water.	<ul style="list-style-type: none"> - Disposable cloth - Detergent 	Until product is visually dry
	5	Manual Disinfection	Brush all surface of the instrument using a soft brush for 2 minutes. Adjust brushing time if needed. Actuate mobile parts of the instrument 3 times (when applicable) and use a syringe to access all difficult areas. Use a volume of 60 mL of the detergent solution.	<ul style="list-style-type: none"> - Soft bristle brush - Detergent - Syringe 	Manual cleaning time duration is complete when the device's surface, joints, and crevices have been manually cleaned
	6	Rinse and finition	Soak the instrument in reverse osmosis water for 1 minute. Use a syringe and 60mL of reverse osmosis water to access all difficult areas. Repeat soaking 2 additional times using fresh water	<ul style="list-style-type: none"> - Reserve osmosis water - Syringe 	Minimum 1 minute
	7	Final Rinse	Rinse carefully with distilled water	<ul style="list-style-type: none"> - Distilled water 	Minimum 1 minute
	8	Final Drying	Dry devise utilizing medical quality filtered air	<ul style="list-style-type: none"> - Medical quality filtered air 	Until product is visually dry
	9	Tidying	The instruments have to be placed in the right position inside the trays as described on the position map fixed on the tray for the sterilization.		

Note: In case of **suspicion** of **Creutzfeld-Jakob Disease (CJD)**, incinerate the ancillaries and the fixation pliers

In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent cleaning process for which it will be fully responsible.

STERILIZATION OF ANCILLAIRES

Sterilization is performed under the responsibility of the healthcare centers.

The following prescribed methods and materials must be used to reduce the risk of transmission of UTA, (French health ministry - DGS/R13/2011/449 dated 01/12/2011).

1- The sterilization in autoclave is to make in a **specific container** (fenced basket, Tray, Plastic tray placed in a peelable pouches autoclavable).

NB: It **isn't recommended** to realize the sterilization with peelables pouches **autoclavable** without a container. Ancillaries can pierce the films components the pouches.

2-A sterilization **BY AUTOCLAVE** is advised following this protocol:

Steam sterilization parameters	Values	Values
Cycle Type	Pre-vacuum	Pre-vacuum
Set Point Temperature	134°C / 273°F	132°C / 269.6°F
Cycle Time	Minimum 3 minutes	Minimum 4 minutes
Dry Time	Minimum 20 minutes	Minimum 20 minutes
Cool Time (inside and outside the autoclave chambers)	Minimum 40 minutes	Minimum 40 minutes

3-Remark: After each cycle of cleaning/sterilization, insure of all the good working of the ancillaries: verify the integrality and the correct operation, of the locking systems (clipping, prehension) without excessive play.

4- If the ancillaries utilization is different, after sterilization stock them in a tray on a peelable pouches.

In case of no respect of the instructions mentioned above, the healthcare center will have to apply a validated equivalent sterilization process for which it will be fully responsible.

PACKAGING, STORAGE AND STORING OF TUNNELLING DEVICE LIFT®

The tunnelling device must be stored conditioned in order to avoid any deterioration.

Sterilization as well as the preceding stages such as cleaning and decontamination, are under the responsibility of the centers of care which must use a material and methods validated (Reduction of the risk of transmission of ATNC, circular DGS/R13/2011/449 of the 01/12/2011).

COUSIN BIOTECH recommends a steam sterilization. The instruments are compatible with a steam treatment with 134°C during 18min.

Always make sure that ancillaries are complete and in good working order after each cleaning/ sterilization cycle and before use.

HANDLING OF TUNNELLING DEVICE LIFT®

Before use, check that the packaging is not damaged in any way that might suggest compromise of the contents.

The tunnelling device must be stored in a container specifically designed for this use in equivalent packaging.

To prevent that the instruments are not in contact with objects which can deteriorate their surface.

Proceed to a visual examination of each instrument, in order to detect any possible deterioration.

The tunnelling device LIFT® are only suitable for use with original COUSIN BIOTECH instruments.

Moreover, the instruments must be used by a surgeon correctly trained to the implantation of prostheses LIFT®.

Use of tunnelling device LIFT® during the intervention

Refer to the following operative technique instructions for further information.

LIFT® tunnelling device are designed for use in implantation of the LIFT® prosthesis, using the traction threads at the ends of the tape. During the procedure, these threads are attached to the ends of needles which are used to position the tape according to the pelvic anatomy prior to implantation of the device.

LIFT® tunnelling device come in different shapes according to the route of implantation: curved needles for the vaginal or supra-pubic route, and spiral needles for the in-out and out-in trans-obturator routes.

Warnings on the tunnelling device

-Tunnelling device LIFT® are medical devices of Class I, intended for a temporary use or reusable.

-Clean and decontaminate tunnelling device LIFT® after use, employing an appropriate method (to reduce the risk of non-conventional transmissible agent transmission).

If there is a possibility that the patient has contracted Creutzfeldt-Jakob Disease (CJD), proceed with an appropriate cleaning and decontamination method.

If the patient has contracted CJD, the instruments must be incinerated.

IMPORTANT

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

IMPLANT CARD

Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.

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