# ADHESIX<sup>®</sup> BIORING<sup>®</sup>

## REPLACEMENT KIT WITH ADHESIVE PORT

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**NOT178/221128**

Version du 28/11/2022

## ADHESIX® BIORING® ADJUSTABLE GASTRIC BANDING REPLACEMENT KIT WITH ADHESIVE PORT

### DESCRIPTION

- The kit includes the following: rinse catheter, catheter, connector, plug, Huber needle, port and connection system.
- The various components have been designed specifically for use with the **ADHESIX® BIORING®** adjustable gastric banding.

### MATERIALS

The injection port is composed of polyoxymethylene with barium sulphate. The septum is made from dimethylsiloxane. A textile component is fixed to the port by means of a polyethylene terephthalate thread. The woven fabric/mesh comprises polypropylene coated with an adhesive resorbable hydrogel (polyvinylpyrrolidone - Polyethylene glycol).

The septum is made from silicone. The catheter consists of silicone (dimethylsiloxane) and barium sulphate labelled with biocompatible ink. The connection system comprises silicone (dimethylsiloxane) and barium sulphate. The plug is stainless steel and silicone (dimethylsiloxane). The connector is made from stainless steel.

No materials of human or animal origin were used. Semi-resorbable.

### INDICATIONS

The kit is intended for use in the event of pre- and post-surgical complications following implantation of **ADHESIX® BIORING®** gastric banding, and in repeat surgery following complications (port infection, leakage, damaged catheter).

### MRI SAFETY

The implant has not been tested for safety and compatibility with postoperative Magnetic Resonance Imaging (MRI). This replacement kit mostly comprises dimethyl siloxane and other components deemed non-ferromagnetic in these quantities, according to the literature. It also comprises stainless steel. Thus, the use of this replacement kit in an MRI context depends on the conditions of individual investigations.

### PERFORMANCE

- The port is an X-Ray-opaque injection site and allows per-cutaneous adjustment of the gastric banding diameter.
- The connection system clips on to secure the connection between the port and the catheter

The injection port is self-adhesive. Sutures are not required. Simply position it on the muscle as follows:

- insert into the aponeurosis of the rectus abdominis muscle under the left costal edge,
- insert subcutaneously under the left costal edge,
- insert into the presternal region, on the xiphoid apophysis,
- insert in any other position at the surgeon's discretion.

*Plug:* In the event of an infection at the junction between the port and the catheter, the port is replaced with the plug until a new port is fitted. The plug is connected to the free end of the catheter of the **ADHESIX® BIORING®** adjustable gastric banding.

*Connector:* When the catheter on the **ADHESIX® BIORING®** adjustable gastric banding is damaged (cut, leakage), the kit connector and catheter can be used to replace the damaged system.

Always refer to the instructions for the **ADHESIX® BIORING®** adjustable gastric banding marketed by COUSIN BIOTECH.

The product must be implanted **only** by a qualified surgeon experienced in using the **ADHESIX® BIORING®** System.

### CONTRA-INDICATIONS

Do not use the implant in the following circumstances:

- Allergy to any of the components,
- Growing children,
- Do not implant to fit with a gastric band that is not **ADHESIX® BIORING®** adjustable gastric banding.

Inflammatory disease or congenital or acquired anomalies of the gastrointestinal tract,

Severe organic disease,

Potential bleeding disorder of the upper gastrointestinal tract or patients who regularly or systematically take aspirin or non-steroidal anti-inflammatory drugs,

Diagnosis or pre-existing symptoms or family history of connective tissue autoimmune disease,

Pregnancy,

Patients with a drug or alcohol addiction,

Patients who refuse to accept the dietary restrictions imposed by the procedure,

Serious psychiatric disorders.

### POSSIBLE COMPLICATIONS

Complications that may arise from the use of an adjustable gastric band include risks associated with the anaesthesia, surgical technique and the patient's level of intolerance to any foreign body implanted in the body.

**Peri-operative complications** (due to the surgical technique):

- Injury to the spleen
- Perforation of the stomach or oesophagus
- Other types of bleeding

**In the immediate postoperative period:**

- |  |                                   |
|--|-----------------------------------|
| - Nausea, vomiting (especially in the first few days post-procedure and when the patient eats more than recommended) | - Thrombo-embolic complications   |
| - Bloating due to gas  | - Infection-related complications |
| - Dysphagia, dehydration   | - Postoperative bleeding          |

### **Delayed postoperative complications:**

- Dilated upper pouch,
- Band slippage,
- Gastroesophageal reflux,
- Pyrosis,
- Gastritis of the upper pouch,
- Band migration,
- Displacement of the injection port,
- Constrictive obstruction (due for example to oedema, food, inadequate calibration, band slippage or to twisted or dilated upper pouch),
- Band deflation (due to leakage in the band, injection port or connection tube),
- Rapid weight loss may lead to symptoms of malnutrition, anaemia and related complications,
- Infection of the injection port,
- Constipation.

### **EXPLANTATION AND DISPOSAL OF DEVICES**

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery - Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

### **INFORMATION REQUESTS AND COMPLAINTS**

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

### **PRECAUTIONS FOR USE AND STORAGE**

- The components of the implantable port and connecting system must only be used with **ADHESIX® BIORING®** adjustable gastric banding.
- The band can be inflated immediately using physiological saline (0.9% NaCl). ONLY use a **Huber needle** to inflate or deflate the band. The volume and pre-inflation interval are at the discretion of the practitioner/surgeon (see the table below for maximum volume).

Band inflation: maximum volume of physiological saline according to band size:

<b>Band sizes (diameter: mm)</b>	<b>Maximum volume of physiological saline (cm<sup>3</sup>)</b>
Ø 21/24	10
Ø 27/29	10
Ø 36/38	15

- A typical band-inflation algorithm is available on request.
- Store in a dry place, away from light and at room temperature in the original packaging.
- Check the integrity of the packaging (do not use if the labels and/or protective packaging are damaged)
- Do not use if the device is damaged or has passed its expiry date.

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[Inhoudsopgave](#)

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**Patient information leaflet ADHESIX® BIORING®**  
**Replacement kit with adhesive port**

Item	Patient information leaflet
1	<p>a. <b>The name of the device</b> : ADHESIX® BIORING® Replacement kit with adhesive port</p> <p>b. <b>The model of the device</b> : BCBRIREADU</p>
2	<p>a. <b>The intended purpose of the device</b> : Kit to be used when a ASHEDIX® BIORING® gastric banding has been implanted and in the event of preoperative and postoperative complications, repeated surgery on patients following complications (infection of the port, leakage, catheter damaged).</p> <p>b. <b>The kind of patient on whom the device is intended to be used</b> : Persons presenting morbid obesity and weight loss resistance, i.e. persons whose Body Mass Index [BMI, formula = weight in kilograms/(height in metres)<sup>2</sup>] is greater than or equal to 40 kg/m<sup>2</sup>, or 35 to 40 kg/m<sup>2</sup> when combined with at least one comorbidity. The surgeon can prescribe the adjustable gastric band at his/her own discretion.</p>
3	<p><b>Special operating instructions for the use of the device</b> : No special operating instructions for the use of the device</p>
4	<p>a. <b>The intended performance of the device</b> :</p> <ul style="list-style-type: none"> <li>- The port is an X-Ray-opaque injection site and allows per-cutaneous adjustment of the gastric banding diameter.</li> <li>- The connection system clips on to secure the connection between the port and the catheter. The injection port is self-adhesive. Sutures are not required. Simply position it on the muscle as follows: <ul style="list-style-type: none"> <li>- insert into the aponeurosis of the rectus abdominis muscle under the left costal edge,</li> <li>- insert subcutaneously under the left costal edge,</li> <li>- insert into the presternal region, on the xiphoid apophysis,</li> <li>- insert in any other position at the surgeon's discretion.</li> </ul> </li> </ul> <p>Plug: In the event of an infection at the junction between the port and the catheter, the port is replaced with the plug until a new port is fitted. The plug is connected to the free end of the catheter of the ADHESIX® BIORING® adjustable gastric banding.</p> <p>Connector: When the catheter on the ADHESIX® BIORING® adjustable gastric banding is damaged (cut, leakage), the kit connector and catheter can be used to replace the damaged system. Always refer to the instructions for the ADHESIX® BIORING® adjustable gastric banding marketed by COUSIN BIOTECH.</p> <p>The product must be implanted only by a qualified surgeon experienced in using the ADHESIX® BIORING® System.</p>

This leaflet was revised in October 2022

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Item	Patient information leaflet
4	<p><b>b. Undesirable side effects :</b></p> <p>Complications that may arise from the use of an adjustable gastric band include risks associated with the anaesthesia, the surgical technique, and the degree of intolerance of the patient to any foreign body implanted in the body.</p> <p><b>Peri-operative complications (due to the surgical technique) :</b> Injury to the spleen, Perforation of the stomach or oesophagus, Other types of bleeding</p> <p><b>In the immediate postoperative period :</b> Nausea, vomiting (especially in the first few days post-procedure and when the patient eats more than recommended), thrombo-embolic complication, bloating due to gas, dysphagia, dehydration, infection-related complications, postoperative bleeding</p> <p><b>Delayed postoperative complications :</b> Dilated upper pouch, Band slippage, Gastroesophageal reflux, Pyrosis, Gastritis of the upper pouch, Band migration, Displacement of the injection port, Constrictive obstruction (due for example to oedema, food, inadequate calibration, band slippage or to twisted or dilated upper pouch), Band deflation (due to leakage in the band, injection port or connection tube), Rapid weight loss may lead to symptoms of malnutrition, anaemia and related complications, Infection of the injection port, Constipation.</p>
5	<p><b>When the ADHESIX® BIORING® Replacement kit with port Should Not be Used :</b></p> <p>Do not use the implant in the following circumstances :</p> <ul style="list-style-type: none"> <li>- Allergy to any of the components</li> <li>- Growing children</li> <li>- Do not implant to fit with a gastric band that is not a ADHESIX® BIORING® adjustable gastric banding</li> <li>- Inflammatory disease or congenital or acquired anomalies of the gastrointestinal tract.</li> <li>- Severe organic disease.</li> <li>- Potential bleeding disorder of the upper gastrointestinal tract or patients who regularly or systematically take aspirin or non-steroidal anti-inflammatory drugs.</li> <li>- Diagnosis or preexisting symptoms or family history of connective tissue autoimmune disease.</li> <li>- Pregnancy.</li> <li>- Patients with a drug or alcohol addiction.</li> <li>- Patients who refuse to accept the dietary restrictions required by the procedure.</li> </ul> <p>Serious psychiatric disorders.</p>

This leaflet was revised in October 2022

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

Item	Patient information leaflet								
6	<p><b>MRI Information :</b></p> <p>The implant has not been tested for safety and compatibility with postoperative Magnetic Resonance Imaging (MRI). This replacement kit mostly comprises dimethyl siloxane and other components deemed non-ferromagnetic in these quantities, according to the literature. It also comprises stainless steel. Thus, the use of this replacement kit in an MRI context depends on the conditions of individual investigations.</p>								
7	<p><b>Precautions :</b></p> <p>The components of the implantable port and of the connecting system have to be used only with ADHESIX® BIORING® adjustable gastric banding.</p> <p>The band can be inflated immediately using physiological saline (0.9% NaCl). ONLY use a Huber needle to inflate or deflate the band. The volume and pre-inflation interval are at the discretion of the practitioner/surgeon (see the table below for maximum volume).</p> <p>Band inflation: maximum volume of physiological saline according to band size:</p> <table border="1" data-bbox="379 922 1310 1111"> <thead> <tr> <th>Band sizes (diameter : mm)</th> <th>Maximum volume of physiological saline (cm<sup>3</sup>)</th> </tr> </thead> <tbody> <tr> <td>Ø 21/24</td> <td>10</td> </tr> <tr> <td>Ø 27/29</td> <td>10</td> </tr> <tr> <td>Ø 36/38</td> <td>15</td> </tr> </tbody> </table> <ul style="list-style-type: none"> <li>- A typical band-inflation algorithm is available on request.</li> <li>- Store in a dry place, away from light and at room temperature in the original packaging.</li> <li>- Check the integrity of the packaging (do not use if the labels and/or protective packaging are damaged)</li> <li>- Do not use if the device is damaged or has passed its expiry date.</li> <li>- Do not re-use</li> <li>- Do no resterilize</li> </ul>	Band sizes (diameter : mm)	Maximum volume of physiological saline (cm <sup>3</sup> )	Ø 21/24	10	Ø 27/29	10	Ø 36/38	15
Band sizes (diameter : mm)	Maximum volume of physiological saline (cm <sup>3</sup> )								
Ø 21/24	10								
Ø 27/29	10								
Ø 36/38	15								
8	<p><b>Materials included in the device :</b> The injection port is composed of polyoxymethylene with barium sulphate. The septum is made from dimethylsiloxane. A textile component is fixed to the port by means of a polyethylene terephthalate thread. The woven fabric/mesh comprises polypropylene coated with an adhesive resorbable hydrogel (polyvinylpyrrolidone - Polyethylene glycol). The septum is made from silicone. The catheter consists of silicone (dimethylsiloxane) and barium sulphate labelled with biocompatible ink. The connection system comprises silicone (dimethylsiloxane) and barium sulphate. The plug is stainless steel and silicone (dimethylsiloxane). The connector is made from stainless steel. No materials of human or animal origin were used. Semi-resorbable</p>								

Item	Patient information leaflet
9	<p><b>Reporting a serious incident :</b>            If you encounter any serious incident with your implant ADHESIX<sup>®</sup> BIORING<sup>®</sup> Replacement kit with port you should report it immediately to your doctor. You should also report the incident to Cousin Biotech via the link <a href="mailto:serviceclients@cousin-surgery.com">serviceclients@cousin-surgery.com</a> or <a href="mailto:materiovigilance@cousin-biotech.com">materiovigilance@cousin-biotech.com</a> and to the Therapeutic Goods Administration using the following website: <a href="http://www.tga.gov">http://www.tga.gov</a></p>
10	<p><b>Legal Manufacturer :</b> Cousin Biotech, 8 rue de l'Abbé Bonpain, 59117 Wervicq-Sud, France  <b>Manufacturing Facility :</b> Cousin Biotech, Allée des roses, 59117 Wervicq-Sud, France</p>

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Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z - N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille
Cousin Biotech SAS au capital de 320 000 € - NAF 1396Z- N° TVA FR 34 398 460 261 - Siren 398 460 261 RCS Lille