# BIOMESH<sup>®</sup> CA.B.S'Air<sup>®</sup> Light

STERILE EXTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT  
IMPLANT

---

## EXTRAPERITONEAL

en	<a href="#">Instructions for use</a>	Page	2
fr	<a href="#">Notice d'instructions</a>	Page	4
de	<a href="#">Gebrauchsanweisung</a>	Seite	6
it	<a href="#">Istruzioni per l'uso</a>	Pagina	8
es	<a href="#">Instrucciones de uso</a>	Pagina	10

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**CE**  
**1639**

**NOT232-230223**  
Version du 23/02/2023

**Caution : Federal law (USA) restricts this device to sale,  
distribution and use by or on the order of a physician.**

# BIOMESH® CA.B.S.'Air® Light

## STERILE EXTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS SINGLE USE STERILE PRODUCT EXTRA-PERITONEAL USE

### **DESCRIPTION**

The **BIOMESH® CA.B.S.'Air® Light** prosthesis consists of:

- A parietal prosthesis made of a plug in polypropylene and two polyethylene terephthalate and ePTFE fixation threads crimped with stainless steel needles.
- A balloon which allows the deployment of the parietal prosthesis (by using the syringe included in the medical device).

### **IMPLANTED MATERIALS :**

Polypropylene, Polyethylene terephthalate, Biocompatible dye (DC Green #6), ePTFE,

### **MATERIALS IN CONTACT WITH THE PATIENT DURING THE SURGERY:**

Stainless steel, Poly vinyl chloride.

**Origin neither human nor animal - Non-absorbable.**

### **INDICATIONS**

The **BIOMESH® CA.B.S.'Air® Light** is intended for the use in all forms of hernia repair requiring reinforcement with non-absorbable support material

### **PERFORMANCES**

The **BIOMESH® CA.B.S.'Air® Light** medical device is particularly designed to be extra-peritoneally implanted for the mini invasive surgical treatment (small incision) of all form of hernia.

The mesh is biocompatible and non-absorbable. It has the advantage of having a very high resistance to suturing, and a great flexibility.

### **CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

### **UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant could generate possible undesirable effects:

- |                     |                           |                  |
|---------------------|---------------------------|------------------|
| - Discomfort/Pain   | - Recurrence              | - Adhesions      |
| - Erosions          | - irritation nearby organ | - Infection      |
| - Inflammation      | - Fistula formation       | - Seroma         |
| - mesh deformation  | - Hematoma                | - Mesh migration |
| - Allergic reaction | - Foreign body reaction   | - Lymphocele     |

### **PRECAUTIONS FOR USE**

**BIOMESH® CA.B.S.'Air® Light** devices are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and trained on the use of the product (knowledge of anatomy and visceral surgery).

**IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE  
IT IS IMPERATIVE THAT THE BALLOON IS WITHDRAWN**



As indicated on the product labelling, the **BIOMESH® CA.B.S.'Air® Light** implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

### **STORAGE OF THE PROSTHESIS**

Store in a dry place away from sunlight and at room temperature in its original packaging.

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

## **SURGICAL PROCEDURE**

<b>SURGICAL PROCEDURE</b>					
<b>1-Installation:</b> - The prosthesis equipped with its balloon crosses completely the hernial defect, -its positioning becomes then entirely retro muscular		<b>2-Deployment:</b> - The expansion and the spreading out of the prosthesis are obtained by inflating the balloon, (30ml), volume of a syringe. - Parietal fixing is carried out by two points thanks to the pre-positioned sutures.		<b>3- Prosthesis in place:</b> -The prosthesis is fixed, -This balloon must imperatively be withdrawn	

### **SURGICAL TECHNIQUE FOR UMBILICAL HERNIA REPAIR**

Hemicircumferential periumbilical incision centred on the hernia. Dissection of the peritoneal sac, enabling its abdominal reintegration by clearing sufficient space behind the aponeurosis for the prosthesis to be implanted. The folded prosthesis is slipped through the aponeurotic orifice of the hernia in order to position it in the preperitoneal space.. Inflation of the expansion balloon of the prosthesis while ensuring that the prosthesis remains centred relative to the hernia orifice. Fixation of the prosthesis using prepositioned sutures passing through the aponeurosis.

#### **Be careful not to suture the balloon to the implant when fixing the prosthesis.**

Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned. If possible, aponeurotic approximation in front of the prosthesis without tension. Subcutaneous and cutaneous closure with absorbable suture material.

### **EXPLANTATION AND ELIMINATION OF DEVICES**

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery - Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

### **INFORMATION REQUESTS AND COMPLAINTS**

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

### **IMPLANT CARD**

**Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.**

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