ResCUBE™

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NOT248 / 230503

Version du 03/05/2023

Caution: Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

ResCUBE™

Sterile Ligament System for Spine Stabilisation

Single use sterile product

DESCRIPTION

The ResCUBE™ System is a temporary spinal implant intended to provide temporary stabilization as a bone anchor during the development of solid bony fusion and aid in the repair of bone fractures.

It consists of:

- A ligament allowing correction and stabilization of the levels, after tensioning
- A malleable leader allowing guidance of the ligament during surgery and attached to the ligament by a PET sheath. The lamella and shrink sleeve are not intended for long term implantation and must be removed after placement of the implant.
- A low-profile connector holding the tension on the ligament.

MATERIALS

Ligament	Polyethylene Terephthalate (PET)	For a long term implantation
Connector	Titanium alloy Ti-6Al-4V ELI	
Malleable lamella	Stainless Steel	Not for long term implantation
Shrink sleeve	Polyethylene Terephthalate (PET)	

All used materials are non resorbable and non-human and non-animal origin.

INDICATIONS

The ResCUBE™ System is a temporary implant for use in orthopedic surgery. The system is intended to provide temporary stabilization as a bone anchor during the development of solid bony fusion and aid in the repair of bone fractures. The indications for use include the following applications:

1. Spinal trauma surgery, that can be treated with sublaminar, interspinous or facet wiring techniques.
2. Spinal reconstructive surgery, incorporated into constructs for the purpose of correction of spinal deformities such as scoliosis, kyphosis and spondylolisthesis.
3. Spinal degenerative surgery, as an adjunct to spinal fusions.

CONTRAINDICATIONS

Do not use ResCUBE™ implant in the following cases:

- General infection or located infection in the intervention area
- Inflammation located in the intervention area
- Pregnancy
- Immunodepression
- Any pathology affecting bone metabolism, such as osteopenia, osteoporosis, metastasis, tumor... can compromise the expected mechanical support for the implant, and therefore, the success of the surgery
- Insufficient tissues or muscles coverage
- Suspected or confirmed sensibility to at least one of the components
- The overweight or obesity can lead to mechanical overload on the implant, compromising the success of the surgery
- Any mental state or activity level non compatible with the postoperative follow-up should compromise the success of the surgery
- Smoking patients are more likely to present an unsatisfactory bone fusion or non-union
- Any other contraindication to an instrumented surgery of the spine

The described contraindications can be absolute or relative and must be evaluated by the surgeon prior to decide to implant.

POTENTIAL UNDESIRABLE SECONDARY EFFECTS

In addition to the risks linked to an instrumented spine surgery, the following effects can be observed, in rare cases:

- Unusual sensations, discomfort or pain due to the presence of the implant
- Migration, deformation or breakage of one or many components of the implant, can cause neurological, motion complications or pain
- Infection, allergy or inflammatory reaction due to the components, or to the presence of microparticles (metallosis)
- Halting the growth of instrumented levels
- Partial or total lost of correction
- Intraoperative or postoperative damage, fissures or fractures, in case of insufficient bone density
- Pseudoarthrosis
- Pulmonary complications
- Feeding difficulties
- Respiratory problems
- Proximal junctional kyphosis (PJK)
- Wound dehiscence
- Blood loss

- Death

In case of undesirable secondary effect, a second surgery may be necessary, in order to remove or replace one or several implants.

MRI SAFETY

The implants are composed of non-ferromagnetic materials and present a geometry non susceptible to generate induced currents. Moreover, as they are fixed to bone or tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular in terms of heating and migration of implant has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded to a safety use for MRI scan until 1.5 or 3 Tesla.

As a precautionary measure, it is recommended to avoid MRI scans within the 48h of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential.

It has to be noted that the devices which present a high contrast with the biological environment can generate « artifacts » that has to be taken into account for the perfect execution and interpretation of imaging exams.

For this purpose, it has to be recommended to the patient who has this implant to warn as far as possible the concerned health professionals (radiologists and radiology operators) about the presence of this implant before these exams.

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

In accordance with the labelling of this product, the implants are for single use; they must not in any case be reused or re-sterilised (potential risks include, but are not limited to, loss of sterility of the products, risk of infection, loss of product efficacy, recurrence, etc.)

PRECAUTIONS FOR USE

General precautions

ResCUBE™ devices must be implanted by spine specialists surgeons, trained to the surgical technique. The information contained in this notice is necessary but not sufficient to control the surgical technique.

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

ResCUBE™ devices are delivered sterile and are single use. They must not be resterilized nor reused (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

Implants must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature.

Preoperative precautions

The implantation decision must be taken by the surgeon, after information and agreement of the patient, considering, indications, contraindications and potential secondary effects.

Before the surgery, it is necessary to check the availability and the compatibility of devices which will be associated (size, materials...).

Peroperative precautions

Before opening, it is recommended to check the packaging integrity. Any damage on packaging may jeopardize the implant sterility. In case of damage on the packaging, the implant must not be used, nor resterilized.

Any implant whose expiration date has passed must not be used.

After opening, it is recommended to check the implant integrity.

Any implant whose have a deteriorated surface must not be implanted.

The implant setting up must be done with the dedicated ancillaries, following the surgical technique recommended by COUSIN BIOTECH.

The malleable metal lamella and the shrink sleeve are not intended to a long term implantation and must be removed after implant set up.

Postoperative precautions

The patient must be informed of the postoperative cares and of precautions.

The device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be done following recommendations of ISO 12891-1:2015 « Implants for surgery - Retrieval and analysis of surgical implants » Part 1: « Retrieval and Handling ».

Any explanted device must be send back, for analysis, following the current protocol. This protocol is available on demand to COUSIN BIOTECH. It is important to note that any implant that should not have been cleaned and disinfected before expedition must be contained in a sealed package.

The elimination of explanted medical device must be conducted in accordance with standards in the country for the disposal of infectious hazards waste.

The elimination of a non-implanted device is not the subject of specific recommendations.

INFORMATION REQUEST AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible.

In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time.

For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

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