NAJA™

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NOT262 / 230503

Version du 03/05/2023

Caution : Federal law (USA) restricts this device to sale,
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

NAJA™

Ligament System for Spine Stabilization

Single use sterile product

DESCRIPTION

The NAJA™ device is a spinal implant system dedicated to the correction and stabilization of the spine during fusion of the instrumented levels.

It consists of:

- A braid allowing correction and stabilization of the levels, after tensioning
- A malleable lamella allowing guidance of the braid during surgery and attached to the braid through a shrink sleeve. The lamella and shrink sleeve are not intended for long term implantation and must be removed after placement of the implant.
- A connector designed to be fixed on a rod and holding the tension on the braid

Because of its conception, and in order to obtain symmetric constructions, two versions are available for each size of the implant: one for the left size, one for the right side.

MATERIALS

Braid	Polyethylene Terephthalate (PET)	For a long term implantation
Connector	Titanium alloy Ti-6Al-4V ELI	
Malleable lamella	Stainless Steel	Not for long term implantation
Shrink Sleeve	Polyethylene Terephthalate (PET)	

All used materials are non resorbable and non-human and non-animal origin.

INDICATIONS

The NAJA™ System is a temporary implant for use in orthopedic surgery. The system is intended to provide temporary stabilization as a bone anchor during the development of solid bony fusion and aid in the repair of bone fractures. The indications for use include the following applications:

Spinal reconstructive surgery, incorporated into constructs for the purpose of correction of spinal deformities such as idiopathic and neuromuscular scoliosis in paediatric population from 8 years old to adulthood.

The NAJA™ System may also be used in conjunction with other medical grade implants made of titanium or cobalt chrome alloy whenever "wiring" may help secure the attachment of the other implants.

CONTRAINDICATIONS

Do not use NAJA™ in the following cases:

- General infection or located infection in the intervention area
- Inflammation located in the intervention area
- Adult population (because no clinical evidence available)
- Pregnancy
- Immunodepression
- Any pathology affecting bone metabolism, such as osteopenia, osteoporosis, metastasis, tumor... can compromise the expected mechanical support for the implant, and therefore, the success of the surgery
- Insufficient tissues or muscles coverage
- Severe fractures such that segments may not be maintained in satisfactory proximate reduction.
- Suspected or confirmed sensibility to at least one of the components
- The overweight or obesity can lead to mechanical overload on the implant, compromising the success of the surgery
- Any mental state or activity level non compatible with the postoperative follow-up should compromise the success of the surgery
- The smoking patients are more likely to present an unsatisfactory bone fusion or non-union
- Any other contraindication to an instrumented surgery of the spine

The described contraindications can be absolute or relative and must be evaluated by the surgeon prior to decide to implant.

POTENTIAL UNDESIRABLE SECONDARY EFFECTS

In addition of the risks linked to an instrumented spine surgery, we can observe, in few cases, the following effects:

- Unusual sensations, discomfort or pain due to the presence of the implant
- Migration, deformation or breakage of one or many components of the implant, can cause neurological, motion complications or pain
- Infection, allergy or inflammatory reaction due to the components, or to the presence of micoparticles (metallosis)
- Halting the growth of instrumented levels
- Partial or total lost of correction
- Intraoperative or postoperative damage, fissures or fractures, in case of insufficient bone density

- Pseudoarthrosis
- Pulmonary complications
- Feeding difficulties
- Respiratory problems
- Proximal junctional kyphosis (PJK)
- Wound dehiscence
- Blood loss
- Death

In case of undesirable secondary effect, a second surgery may be necessary, in order to remove or replace one or many implants.

MRI SAFETY

The implants are composed of non-ferromagnetic materials and present a geometry non susceptible to generate induced currents. Moreover, as they are fixed to bone or tissues, they are unlikely to be mobilised. A priori they can be considered compatible with an MRI scan. Their safety, in particular in terms of heating and migration of implant has been evaluated through bibliographic data by comparison with data available on devices with similar composition, shape and use. This evaluation concluded to a safety use for MRI scan until 1.5 or 3 Tesla.

As a precautionary measure, it is recommended to avoid MRI scans within the 48h of the implant placement, and to inform the person in charge of the scan of the recent implant placement, if such examination is essential.

It has to be noted that the devices which present a high contrast with the biological environment can generate « artifacts » that has to be taken into account for the perfect execution and interpretation of imaging exams.

For this purpose, it has to be recommended to the patient who has this implant to warn as far as possible the concerned health professionals (radiologists and radiology operators) about the presence of this implant before these exams.

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

In accordance with the labelling of this product, the implants are for single use; they must not in any case be reused or re-sterilised (potential risks include, but are not limited to, loss of sterility of the products, risk of infection, loss of product efficacy, recurrence, etc.)

PRECAUTIONS FOR USE

General precautions

NAJA™ devices must be implanted by spine specialists surgeons, trained to the surgical technique. The information contained in this notice is necessary but not sufficient to control the surgical technique.

IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE

NAJA™ devices are delivered sterile and are single use. They must not be resterilized nor reused (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, relapse).

NAJA™ devices are made from titanium, they must be associated with titanium or cobalt chrome arthrodesis material. Any association with other materials will not be guaranteed by COUSIN BIOTECH.

NAJA™ devices are used with rods; it is necessary to ensure that the diameter is equal to the size of the implant, indicated on the label and engraved on the connector.

NAJA™ device is available in two versions, one for the left side, the other for the right side. It is necessary to respect this repartition. The version is indicated on the label and engraved on the connector. The letter L indicates left side and the letter R indicates right side.

Implants must be handled and stored with great care, in a dry place, protected from light and at room temperature.

Preoperative precautions

The implantation decision must be taken by the surgeon, after information and agreement of the patient, considering, indications, contraindications and potential secondary effects.

Before the surgery, it is necessary to check the availability and the compatibility of devices which will be associated (size, materials...).

Peroperative precautions

Before opening, it is recommended to check the packaging integrity. Any damage on packaging may jeopardize the implant sterility. In case of damage on the packaging, the implant must not be used, nor resterilized.

Any implant whose expiration date has passed must not be used.

After opening, it is recommended to check the implant integrity.

Any implant whose have a deteriorated surface must not be implanted.

The implant setting up must be done with the dedicated ancillaries, following the surgical technique recommended by COUSIN BIOTECH.

The malleable lamella and the shrink sleeve are not intended to a long term implantation and must be removed after implant set up.

Postoperative precautions

The patient must be informed of the postoperative cares and of precautions.

The use of a brace could be necessary but is not mandatory. This decision is the surgeon one.

The device has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the device in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

EXPLANTATION AND ELIMINATION OF DEVICES

Explantation and handling should be done following recommendations of ISO 12891-1:2015 « Implants for surgery – Retrieval and analysis of surgical implants » Part 1: « Retrieval and Handling ».

Any explanted device must be sent back, for analysis, following the current protocol. This protocol is available on demand to COUSIN BIOTECH. It is important to note that any implant that should not have been cleaned and disinfected before expedition must be contained in a sealed package.

The elimination of explanted medical device must be conducted in accordance with standards in the country for the disposal of infectious hazards waste.

The elimination of a non-implanted device is not the subject of specific recommendations.

INFORMATION REQUEST AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible.

In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time .

For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

[Table of contents](#)

Instructions for use of the ancillary

en

Single use NAJA™ Tensioner NAJA™ ancillaries

DESCRIPTION

This ancillary is intended for the placement of the NAJA™ ligament system for spinal stabilization.

This ancillary is a Class Is medical device, intended for single use and is sterile supplied.

MATERIALS

The different elements of the ancillary are made from:

- IXEF®

- Stainless steel (for metallic components)

All the materials used are non-absorbable and are neither of human nor animal origin.

PRECAUTIONS FOR USE

The ancillary must be stored in its original packaging, which ensures safety.

Carry out a visual examination of each instrument before each use in order to detect any possible deterioration

The instruments must not be used except with elements corresponding to the original COUSIN BIOTECH system.

Use of specific ancillary COUSIN BIOTECH system is required for implantation of the implant.

In addition, the instruments must not be used except by a surgeon with appropriate training in the placement of NAJA™ implants.

Sterile device (Sterilized using Gamma irradiation)

Store the ancillary in its original protective packaging, and do not remove them from the packaging until immediately before use.

Prior to use, check the product expiration date and verify the integrity of the sterile packaging. Do not use, if the package is damaged.

Single-use device

According to the labelling of this product, the ancillary is a single use device.

Products intended for single-use must not be re-used. Re-use or reprocessing (e.g. cleaning and resterilization) may compromise the structural integrity of the device and/or lead to device failure which may result in patient injury, illness or death. Furthermore, reuse or reprocessing of single-use devices may create a risk of contamination e.g. due to the transmission of infectious material from one patient to another. This could result in injury or death of the patient or user.

Contaminated ancillary must not be reprocessed. Any ancillary that has been contaminated by blood, tissue, and/or bodily fluids/matter should never be used again and should be handled according to hospital protocol. Even though it may appear undamaged, the ancillary may have small defects and internal stress patterns that may cause material fatigue.

Storage precautions of the ancillary

The device must be handled and stored with great care, in a dry place, protected from sunlight and at room temperature in its original packaging.

INFORMATION REQUEST AND CLAIMS

Following its quality policy, COUSIN BIOTECH is committed to make every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber ...) had a claim or cause of dissatisfaction with a product in terms of quality, safety or performance, he must inform COUSIN BIOTECH as soon as possible.

In case of failure of an implant or if it contributed to cause serious adverse effects for the patient, the health center must follow the legal procedures in his country, and inform COUSIN BIOTECH in the shortest time.

For any correspondence, thank you to specify the reference, batch number, the coordinates of a reference and a comprehensive description of the incident or claim.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

If further information is needed or required, please contact your COUSIN BIOTECH representative or distributor, or COUSIN BIOTECH at spine@cousin-biotech.com.

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