# BIOMESH<sup>®</sup> CA.B.S'Air<sup>®</sup> SR

## SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT

---

en	<a href="#">Instructions for use</a>	Page	2
fr	<a href="#">Notice d'instructions</a>	Page	5
de	<a href="#">Gebrauchsanweisung</a>	Seite	8
it	<a href="#">Istruzioni per l'uso</a>	Pagina	11
es	<a href="#">Instrucciones de uso</a>	Pagina	14
ελ	<a href="#">Οδηγίες χρήσης</a>	Side	17
nl	<a href="#">Gebruiksaanwijzing</a>	Pagina	20

### COUSIN BIOTECH

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**NOT138-230602**

Version du 02/06/2023

**Caution : Federal law (USA) restricts this device to sale, distribution and use by or on the order of a physician.**

This release is the last update of the instructions of use and replace the previous edition

**Instructions for use**

**BIOMESH® CA.B.S.'Air® SR**  
**SEMI RESORBABLE PARIETAL REINFORCEMENT IMPLANT**  
**Expandable parietal prosthesis with inflatable balloon for visceral surgery**  
**STERILE PRODUCT FOR SINGLE USE**

**DESCRIPTION**

The **BIOMESH® CA.B.S.'Air® SR** prosthesis consists of:

- A parietal prosthesis made of a plug in polypropylene and Poly L Lactic acid and two polyethylene terephthalate and ePTFE fixation threads crimped with stainless steel needles. This prosthesis is supplied with or without polypropylene and Poly L Lactic acid mesh (see indication on the packaging).
- A balloon which allows the deployment of the parietal prosthesis (by using the syringe included in the medical device).

**IMPLANTED MATERIALS :**

Polypropylene, Polyethylene terephthalate, Biocompatible dye (DC Green #6), ePTFE, Poly L Lactic acid .

**MATERIALS IN CONTACT WITH THE PATIENT DURING THE SURGERY :**

Stainless steel, Poly vinyl chloride.

No human neither animal origin. Semi resorbable.

**INDICATIONS**

Repair and parietal reinforcement of inguinal, umbilical and femoral hernias.

**PERFORMANCES**

The **BIOMESH® CA.B.S.'Air® SR** medical device is particularly designed to be extra-peritoneally implanted for the mini invasive surgical treatment.

Products, biocompatible and semi resorbable, cause reactional fibrosis taking over the reinforcement after six months of implantation. They have the advantage of having a very high resistance to suture, have a great flexibility and allow an optimum and rapid integration and colonisation.

**CONTRAINDICATIONS**

Do not use in the following cases:

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant is susceptible to generate possible undesirable side effects that may lead to reoperation:

- |   |                           |                  |
|---|---------------------------|------------------|
| - Discomfort/Pain   | - Recurrence              | - Adhesions      |
| - Erosions  | - irritation nearby organ | - Infection      |
| - Inflammation  | - Fistula formation       | - Seroma         |
| - mesh deformation  | - Hematoma                | - Mesh migration |
| - Allergic reaction   | - Foreign body reaction   | - Lymphocele     |
| - A small additional risk of male infertility has been identified after bilateral groin hernia repair using meshes (open or laparoscopic) |                           |                  |

**PRECAUTIONS FOR USE**

**BIOMESH® CA.B.S.'Air® SR** devices are delivered sterile (ethylene oxide sterilisation).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging.

Do not use if the device is out of date. COUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and formed with the use of the product (knowledge of anatomy and visceral surgery).

**IMPORTANT : DO NOT REUSE - DO NOT RESTERILIZE**  
**THIS BALLOON MUST BE IMPERATIVELY WITHDRAWN**

As indicated on the product labelling, the **BIOMESH® CA.B.S.'Air® SR** implant is designed for single use. It should in no event be reused and/or resterilized (potential risks include but are not limited to: loss of product sterility, risk of infection, loss of effectiveness, recurrence).

**EXPLANTATION AND ELIMINATION OF DEVICES**

Devices should be retrieved and handled according to the recommendations of ISO 12891-1:2015 "Implants for surgery - Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling".

Any explanted device must be sent back for analysis, following the current protocol. This protocol is available on request from COUSIN BIOTECH. It is important to note that any implant that must not be cleaned or disinfected prior to despatch must be contained in sealed packaging.

The removed medical device must be disposed of in accordance with the country's standards for the disposal of infectious waste. There are no specific recommendations for the disposal of a non-implanted device.

**STORAGE OF THE PROSTHESIS**

To be stored in a dry place away from sunlight and at room temperature in its original packaging.

**INFORMATION REQUESTS AND COMPLAINTS**


Following its quality policy, COUSIN BIOTECH is committed to making every effort to produce and supply a high quality medical device. However, if a health professional (client, user, prescriber, etc.) has a complaint or cause for dissatisfaction with a product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In case of failure of an implant or if it has contributed to causing a serious adverse reaction in the patient, the health centre must follow the legal procedures in that country, and inform COUSIN BIOTECH immediately.

For any correspondence, please specify the reference, batch number, contact person's details and a comprehensive description of the incident or complaint.

Brochures, documentation and surgical technique are available on request from COUSIN BIOTECH and its distributors.

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor

**SURGICAL PROCEDURE**

<b>SURGICAL PROCEDURE</b>					
The prosthesis should always be installed with an anterior reinforcement for inguinal and femoral hernia.					
<p><b>1-Installation:</b></p> <ul style="list-style-type: none"> <li>- The prosthesis equipped with its balloon crosses completely the hernial defect,</li> <li>-its positioning becomes then entirely retro muscular</li> </ul>		<p><b>2-Deployment:</b></p> <ul style="list-style-type: none"> <li>- The expansion and the spreading out of the prosthesis are obtained by inflating the balloon, (30ml), volume of a syringe.</li> <li>- Parietal fixing is carried out by two points thanks to the pre-positionned switched threads.</li> </ul>		<p><b>3- Prosthesis in place:</b></p> <ul style="list-style-type: none"> <li>-The prosthesis is fixed,</li> <li>-This balloon must imperatively be withdrawn</li> </ul>	

**SURGICAL TECHNIQUE FOR INGUINAL HERNIA REPAIR**

Horizontal inguinal incision of about 5 cm, beginning immediately above the pubic spine. This incision extends down to the level of the external abdominal oblique muscle. Placement of retractors. Short incision of the external abdominal oblique muscle at the top of the external orifice of the inguinal canal. Exploration of the inguinal canal enables 1) direct and indirect inguinal hernias to be differentiated, and 2) the dissection and removal of prehernia lipomas.

**DIRECT INGUINAL HERNIA:**

Simple exposure of the patulous transversalis fascia. Opening of the transversalis fascia. Retraction of the preperitoneal fatty tissue in order to create the space for a correct centred placement of the prosthesis. The prosthesis is introduced folded into the preperitoneal space that was created. Inflation of the balloon in order to expand the prosthesis while ensuring its proper centring. Fixation of the prosthesis with two prepositioned sutures that are passed through the posterior wall of the inguinal canal. Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned. The sutures are tied without tension.

**Be careful not to suture the balloon to the implant when fixing the prosthesis.**

Tension-free closure of the transversalis fascia in front of the prosthesis by a continuous, absorbable suture.

**EXTERNAL OBLIQUE HERNIA:**

Opening of the spermatic cord. Detection of the peritoneal sac of the hernia. The parts of the spermatic cord are located and isolated. Deep dissection of the peritoneal sac in order to facilitate its reintegration into the abdominal cavity and to expose sufficient preperitoneal space for the correct placement of the prosthesis. After retraction of the peritoneal sac of the hernia, the prosthesis is placed folded in the preperitoneal space. Inflation of the balloon enabling expansion of the prosthesis while ensuring that it remains well-centred relevant to the deep orifice of the inguinal canal. Fixation of the prosthesis with two prepositioned sutures that are passed through the posterior wall of the inguinal canal without tension.

**Be careful not to suture the balloon to the implant when fixing the prosthesis.**

Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned.

The wall must be reinforced by placing a second prosthesis around the spermatic cord using the Lichtenstein technique. Closure of the incision of the external abdominal oblique muscle with absorbable suture material. Subcutaneous approximation and cutaneous closure with a continuous intradermal, absorbable suture.

**SURGICAL TECHNIQUE FOR UMBILICAL HERNIA REPAIR**

Hemicircumferential periumbilical incision centred on the hernia. Dissection of the peritoneal sac, enabling its abdominal reintegration by clearing sufficient space behind the aponeurosis for the prosthesis to be implanted. The folded prosthesis is slipped through the aponeurotic orifice of the hernia in order to position it in the preperitoneal space. Inflation of the expansion balloon of the prosthesis while ensuring that the prosthesis remains centred relative to the hernia orifice. Fixation of the prosthesis using prepositioned sutures passing through the aponeurosis.

**Be careful not to suture the balloon to the implant when fixing the prosthesis.**

Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned. If possible, aponeurotic approximation in front of the prosthesis without tension. Subcutaneous and cutaneous closure with absorbable suture material.

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