# BIOMESH<sup>®</sup> CA.B.S'Air<sup>®</sup>

## STERILE EXTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS

---

### EXTRAPERITONEAL

en	<a href="#">Instructions for use</a>	Page	2
fr	<a href="#">Notice d'instructions</a>	Page	5
de	<a href="#">Gebrauchsanweisung</a>	Seite	8
it	<a href="#">Istruzioni per l'uso</a>	Pagina	11
es	<a href="#">Instrucciones de uso</a>	Pagina	14
ελ	<a href="#">Οδηγίες χρήσης</a>	Σελίδα	17

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**NOT013 / 230626**

Version du 26/06/2023

Caution : Federal law (USA) restricts this device to sale,  
distribution and use by or on the order of a physician.

This release is the last update of the instructions of use and replace the previous edition

**BIOMESH® CA.B.S.'Air®**  
**STERILE EXTRAPERITONEAL NON-RESORBABLE PARIETAL REINFORCEMENT IMPLANTS**  
**STERILE PRODUCT FOR SINGLE USE**  
**EXTRA-PERITONEAL**

**DESCRIPTION**

The BIOMESH® CA.B.S.'Air® prosthesis consists of (according to the references on the box):

- A parietal prosthesis made of a plug and a polypropylene mesh and two polyethylene terephthalate and ePTFE fixation threads crimped on to stainless steel needles.
- A balloon which allows the deployment of the parietal prosthesis

**IMPLANTED MATERIALS**

Polypropylene, Polyethylene terephthalate, Biocompatible dye, ePTFE.

Materials in contact with the patient during surgery

Stainless steel, Poly vinyl chloride.

**Origin neither human nor animal - Non-absorbable.**

**INDICATIONS**

The BIOMESH® CA.B.S.'Air® prosthesis is intended for the use in all forms of hernia repair requiring reinforcement with non-absorbable support material.

**PERFORMANCES**

The BIOMESH® CA.B.S.'Air® prosthesis is particularly designed to be implanted **extra-peritoneally** for the mini invasive surgical treatment (small incision) of all form of hernia.

The BIOMESH® CA.B.S.'Air® is biocompatible non resorbable products. The mesh has a very high resistance to sutures with a great flexibility.

**CONTRAINDICATIONS**

Do not use in the following cases :

- Allergy to any of the components
- Infected site
- Pregnancy
- Growing children
- Anticoagulant therapy

**UNDESIRABLE SIDE EFFECTS**

Like any implantable medical device, this implant could generate possible undesirable effects:

- |                   |                      |                     |                            |
|-------------------|----------------------|---------------------|----------------------------|
| - Discomfort/Pain | - Adhesion formation | - Hematoma          | - Foreign Body Reaction    |
| - Infection       | - Fistula formation  | - Mesh migration    | - Irritation near by organ |
| - Recurrence      | - Erosion            | - Mesh deformation  | - Lymphocele               |
| - Inflammation    | - Seroma             | - Allergic reaction |                            |

**PRECAUTIONS FOR USE**

The BIOMESH® CA.B.S.'Air® prosthesis are delivered sterile (ethylene oxide sterilization).

Before any use, inspect the integrity of packaging and device (of which blister / peelable pouches). Do not use in the event of deterioration of the device and/or the packaging. Do not use if the device is out of date. COUSIN

BIOTECHCOUSIN BIOTECHCOUSIN BIOTECHCOUSIN BIOTECH does not offer any guarantee or recommendation as far as the use of a particular type of means of fixation is concerned.

This device must be implanted **only** by a qualified surgeon and trained with the use of the product (knowledge of anatomy and visceral surgery).

**IMPORTANT: DO NOT REUSE - DO NOT RESTERILIZE**  
**IT IS IMPERATIVE THAT THE BALLOON IS WITHDRAWN**

As specified on the product's labelling, the BIOMESH® CA.B.S.'Air® mesh is for single use only. It can not be re-used and/or re-sterilized (potential risks would be and are not limited to: loss of the product's sterility, risk of infection, loss of the product's efficiency, recurrence)

**STORAGE OF THE PROSTHESIS**


Store in a dry place away from sunlight and at room temperature in its original packaging.

**Important**

For more information about the use of this product, please contact your COUSIN BIOTECH representative or distributor.

**SURGICAL PROCEDURE**

BIOMESH® CA.B.S.'Air® is not intended to be used in the laparoscopic and endoscopic surgery. The prosthesis should always be installed with an anterior reinforcement for inguinal and crural hernia.

SURGICAL PROCEDURE					
<p><b>1-Installation:</b></p> <ul style="list-style-type: none"> <li>- The prosthesis equipped with its balloon crosses completely the hernial defect,</li> <li>-its positioning becomes then entirely retro muscular</li> </ul>		<p><b>2-Deployment:</b></p> <ul style="list-style-type: none"> <li>- The expansion and the spreading out of the prosthesis are obtained by inflating the balloon, (30ml), volume of a syringe.</li> <li>- Parietal fixing is carried out by two points thanks to the pre-positionned switched threads.</li> </ul>		<p><b>3- Prosthesis in place:</b></p> <ul style="list-style-type: none"> <li>-The prosthesis is fixed,</li> <li>-This balloon must imperatively be withdrawn</li> </ul>	

**Surgical technique for inguinal hernia repair**

Horizontal inguinal incision of about 5 cm, beginning immediately above the pubic spine. This incision extends down to the level of the external abdominal oblique muscle. Placement of retractors. Short incision of the external abdominal oblique muscle at the top of the external orifice of the inguinal canal. Exploration of the inguinal canal enables 1) direct and indirect inguinal hernias to be differentiated, and 2) the dissection and removal of prehernia lipomas.

**Direct inguinal hernia:**

Simple exposure of the patulous transversalis fascia. Opening of the transversalis fascia. Retraction of the preperitoneal fatty tissue in order to create the space for a correct centred placement of the prosthesis. The prosthesis is introduced folded into the preperitoneal space that was created. Inflation of the balloon in order to expand the prosthesis while ensuring its proper centring. Fixation of the prosthesis with two prepositionned sutures that are passed through the posterior wall of the inguinal canal. Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned. The sutures are tied without tension.

Be careful not to suture the balloon to the implant when fixing the prosthesis.

Tension-free closure of the transversalis fascia in front of the prosthesis by a continuous, absorbable suture.

**External abdominal oblique hernia:**

Opening of the spermatic cord. Detection of the peritoneal sac of the hernia. The parts of the spermatic cord are located and isolated. Deep dissection of the peritoneal sac in order to facilitate its reintegration into the abdominal cavity and to expose sufficient preperitoneal space for the correct placement of the prosthesis. After retraction of the peritoneal sac of the hernia, the prosthesis is placed folded in the preperitoneal space. Inflation of the balloon enabling expansion of the prosthesis while ensuring that it remains well-centred relevant to the deep orifice of the inguinal canal. Fixation of the prosthesis with two prepositionned sutures that are passed through the posterior wall of the inguinal canal without tension.

**Be careful not to suture the balloon to the implant when fixing the prosthesis.**

Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned. The wall must be reinforced by placing a second prosthesis around the spermatic cord using the Lichtenstein technique. Closure of the incision of the external abdominal oblique muscle with absorbable suture material. Subcutaneous approximation and cutaneous closure with a continuous intradermal, absorbable suture.

### **Surgical technique for umbilical hernia repair**

Hemicircumferential periumbilical incision centred on the hernia. Dissection of the peritoneal sac, enabling its abdominal reintegration by clearing sufficient space behind the aponeurosis for the prosthesis to be implanted. The folded prosthesis is slipped through the aponeurotic orifice of the hernia in order to position it in the preperitoneal space. Inflation of the expansion balloon of the prosthesis while ensuring that the prosthesis remains centred relative to the hernia orifice. Fixation of the prosthesis using prepositioned sutures passing through the aponeurosis.

**Be careful not to suture the balloon to the implant when fixing the prosthesis.**

Exsufflation of the balloon, which is removed while ensuring that the prosthesis remains perfectly positioned. If possible, aponeurotic approximation in front of the prosthesis without tension. Subcutaneous and cutaneous closure with absorbable suture material.

### **EXPLANTATION AND ELIMINATION OF DEVICES**

Explantation and handling should be carried out in accordance with ISO 12891-1:2015 "Implants for surgery - Retrieval and analysis of surgical implants" Part 1: "Retrieval and Handling". Any explanted device must be returned for analysis as per the current protocol. This protocol is available from COUSIN BIOTECH on request. It is important to note that any implant that has not been cleaned and disinfected prior to dispatch must be placed in a sealed package. Explanted medical devices must be disposed of in accordance with local regulations governing the disposal of infectious hazardous waste. The disposal of non-implanted devices is not subject to specific recommendations. The second paragraph covers information requests or complaints.

### **INFORMATION REQUESTS AND COMPLAINTS**

In accordance with its quality policy, COUSIN BIOTECH is committed to producing and supplying high-quality medical devices. However, if health professionals (client, user, prescriber, etc.) wish to lodge a complaint or report dissatisfaction with a given product in terms of quality, safety or performance, they must inform COUSIN BIOTECH as soon as possible. In the event of implant failure or if the implant is implicated in serious adverse effects experienced by the patient, the health centre must follow the legal procedures in the country in question, and inform COUSIN BIOTECH without delay. The reference, batch number, reference contact details and a comprehensive description of the incident or complaint must be stated in all correspondence. Brochures, documentation and surgical techniques are available on request from COUSIN BIOTECH and its distributors.

### **IMPLANT CARD**

**Be sure to scan your implant card immediately upon receipt so that you can keep track of it in case of loss.**

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